

# **V7 Ultrasound System**

## **DICOM Conformance Statement**

**Revision 2.0**  
**System Version 1.01.00**

Created: May 3, 2022  
Last Updated: Aug 16, 2022



## 0 COVER PAGE

Company Name: SAMSUNG MEDISON CO., LTD

Produce Name: V7

Version: 1.01.00

Internal Document Number:

Date: Aug 16, 2022

| Division    | Department         | Name          | Date       | Signature |
|-------------|--------------------|---------------|------------|-----------|
| Authored by | System development | Eun-kwang Kim | 2022/08/16 |           |
| Reviewed by |                    |               |            |           |
| Approved by | System development | Ho-Ki Hong    | 2022/08/16 |           |

# 1 CONFORMANCE STATEMENT OVERVIEW

V7 implements the necessary DICOM services to download worklists from information systems, save acquired US images and Structured Reports to a network storage device, CD or DVD, print to a networked hardcopy device and inform the information system about the work actually done.

Table 1-1 provides an overview of the network services supported by V7.

**Table 1-1  
NETWORK SERVICES**

| SOP Classes                          | User of Service (SCU) | Provider of Service (SCP) |
|--------------------------------------|-----------------------|---------------------------|
| <b>Transfer</b>                      |                       |                           |
| Ultrasound Image Storage             | Yes                   | No                        |
| Ultrasound Multi-frame Image Storage | Yes                   | No                        |
| Comprehensive SR                     | Yes                   | No                        |
| <b>Workflow Management</b>           |                       |                           |
| Modality Worklist                    | Yes                   | No                        |
| Storage Commitment Push Model        | Yes                   | No                        |
| Modality Performed Procedure Step    | Yes                   | No                        |
| <b>Print Management</b>              |                       |                           |
| Basic Grayscale Print Management     | Yes                   | No                        |
| Basic Color Print Management         | Yes                   | No                        |
| <b>Query/Retrieve</b>                |                       |                           |
| Study Root Information Model FIND    | Yes                   | No                        |
| Study Root Information Model MOVE    | Yes                   | No                        |

Provide Storage SCP only Q/R service running.

Table 1-2 provides an overview of the Media Storage Application Profiles supported by V7.

**Table 1-2  
MEDIA SERVICES**

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| <b>Compact Disk - Recordable</b>  |                          |                  |
| STD-US-SC-MF-CDR                  | Yes                      | No               |
| <b>DVD</b>                        |                          |                  |
| STD-US-SC-MF-DVD                  | Yes                      | No               |

## 2 TABLE OF CONTENTS

|         |   |    |
|---------|---|----|
| 0       | COVER PAGE.....   | 2  |
| 1       | CONFORMANCE STATEMENT OVERVIEW .....                      | 3  |
| 2       | TABLE OF CONTENTS.....                                    | 4  |
| 3       | INTRODUCTION .....  | 10 |
| 3.1     | REVISION HISTORY.....                                     | 10 |
| 3.2     | AUDIENCE.....   | 10 |
| 3.3     | REMARKS .....   | 10 |
| 3.4     | DEFINITIONS, TERMS AND ABBREVIATIONS .....                | 11 |
| 3.5     | REFERENCES.....   | 11 |
| 4       | NETWORKING .....  | 12 |
| 4.1     | IMPLEMENTATION MODEL.....                                 | 12 |
| 4.1.1   | Application Data Flow.....                                | 12 |
| 4.1.2   | Functional Definition of AE's.....                        | 13 |
| 4.1.2.1 | Functional Definition of Storage Application Entity.....  | 13 |
| 4.1.2.2 | Functional Definition of Workflow Application Entity..... | 13 |
| 4.1.2.3 | Functional Definition of Hardcopy Application Entity..... | 14 |
| 4.1.2.4 | Functional Definition of the Q/R Application Entity ..... | 14 |
| 4.1.2.5 | Functional Definition of the Q/R Application Entity ..... | 14 |
| 4.1.3   | Sequencing of Real-World Activities .....                 | 15 |
| 4.2     | AE SPECIFICATIONS.....                                    | 18 |
| 4.2.1   | Storage Application Entity Specification.....             | 18 |
| 4.2.1.1 | SOP Classes.....  | 18 |
| 4.2.1.2 | Association Policies.....                                 | 18 |
| 4.2.1.3 | Association Initiation Policy.....                        | 19 |
| 4.2.1.4 | Association Acceptance Policy .....                       | 26 |
| 4.2.2   | Workflow Application Entity Specification.....            | 28 |

|         |  |    |
|---------|--|----|
| 4.2.2.1 | <i>SOP Classes</i> .....   | 28 |
| 4.2.2.2 | <i>Association Establishment Policy</i> .....                          | 29 |
| 4.2.2.3 | <i>Association Initiation Policy</i> .....                             | 30 |
| 4.2.2.4 | <i>Association Acceptance Policy</i> .....                             | 41 |
| 4.2.3   | <i>Hardcopy Application Entity Specification</i> .....                 | 41 |
| 4.2.3.1 | <i>SOP Classes</i> .....   | 41 |
| 4.2.3.2 | <i>Association Policies</i> .....                                      | 41 |
| 4.2.3.3 | <i>Association Initiation Policy</i> .....                             | 42 |
| 4.2.3.4 | <i>Association Acceptance Policy</i> .....                             | 50 |
| 4.2.4   | <i>Q/R Application Entity Specification</i> .....                      | 51 |
| 4.2.4.1 | <i>SOP Classes</i> .....   | 51 |
| 4.2.4.2 | <i>Association Establishment Policy</i> .....                          | 51 |
| 4.2.4.3 | <i>Association Initiation Policy</i> .....                             | 52 |
| 4.2.4.4 | <i>Association Acceptance Policy</i> .....                             | 58 |
| 4.2.5   | <i>STORAGE-SCP Application Entity Specification</i> .....              | 58 |
| 4.2.5.1 | <i>SOP Classes</i> .....   | 58 |
| 4.2.5.1 | <i>Association Establishment Policy</i> .....                          | 58 |
| 4.2.5.2 | <i>Association Initiation Policy</i> .....                             | 59 |
| 4.2.5.3 | <i>Association Acceptance Policy</i> .....                             | 59 |
| 4.3     | <i>NETWORK INTERFACE</i> .....   | 62 |
| 4.3.1   | <i>Physical Network Interface</i> .....                                | 62 |
| 4.3.2   | <i>Additional Protocols</i> .....                                      | 62 |
| 4.3.3   | <i>IPv4 and IPv6 Support</i> .....                                     | 62 |
| 4.4     | <i>CONFIGURATION</i> .....   | 62 |
| 4.4.1   | <i>AE Title/Presentation Address Mapping</i> .....                     | 62 |
| 4.4.1.1 | <i>Local AE Titles</i> .....   | 62 |
| 4.4.1.2 | <i>Remote AE Title/Presentation Address Mapping</i> .....              | 63 |
| 4.4.2   | <i>Parameters</i> .....  | 64 |
| 5       | <i>MEDIA INTERCHANGE</i> .....   | 66 |
| 5.1     | <i>IMPLEMENTATION MODEL</i> .....                                      | 66 |
| 5.1.1   | <i>Application Data Flow</i> .....                                     | 66 |
| 5.1.2   | <i>Functional Definition of AEs</i> .....                              | 66 |
| 5.1.2.1 | <i>Functional Definition of Offline-Media Application Entity</i> ..... | 66 |
| 5.1.3   | <i>Sequencing of Real-World Activities</i> .....                       | 66 |
| 5.1.4   | <i>File Meta Information Options</i> .....                             | 66 |
| 5.2     | <i>AE SPECIFICATIONS</i> .....   | 67 |

|         |   |    |
|---------|---|----|
| 5.2.1   | <i>Offline-Media Application Entity Specification</i> .....   | 67 |
| 5.2.1.1 | <i>File Meta Information for the Application Entity</i> ..... | 67 |
| 5.2.1.2 | <i>Real-World Activities</i> .....                            | 67 |
| 6       | SUPPORT OF CHARACTER SETS.....                                | 69 |
| 7       | SECURITY.....   | 70 |
| 7.1     | GENERAL .....   | 70 |
| 7.2     | SECURE TRANSPORT CONNECTION PROFILES .....                    | 70 |
| 8       | ANNEXES.....  | 72 |
| 8.1     | IOD CONTENTS .....  | 72 |
| 8.1.1   | Created SOP Instances.....                                    | 72 |
| 8.1.1.1 | US or US Multiframe Image IOD.....                            | 72 |
| 8.1.1.2 | Comprehensive Structured Report IOD .....                     | 73 |
| 8.1.1.3 | Common Modules.....   | 74 |
| 8.1.1.4 | US or US Multiframe Image Module.....                         | 78 |
| 8.1.1.5 | Comprehensive Structured Report Modules.....                  | 83 |
| 8.1.2   | Used Fields in received IOD by application .....              | 85 |
| 8.1.3   | Attribute mapping.....  | 86 |
| 8.1.4   | Coerced/Modified Fields.....                                  | 87 |
| 8.2     | DATA DICTIONARY OF PRIVATE ATTRIBUTES.....                    | 87 |
| 8.3     | CODED TERMINOLOGY AND TEMPLATES .....                         | 87 |
| 8.4     | STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES.....    | 88 |
| 8.4.1   | US OR US MULTIFRAME IMAGE STORAGE SOP CLASS .....             | 88 |
| 8.5     | PRIVATE TRANSFER SYNTAXES .....                               | 88 |
| 9       | STRUCTURED REPORT TEMPLATES .....                             | 89 |
| 9.1     | OB-GYN REPORT TEMPLATES.....                                  | 90 |
| 9.1.1   | OB-GYN Ultrasound Procedure Report (TID 5000).....            | 90 |
| 9.1.1.1 | Observation ConText (TID 1001) .....                          | 91 |
| 9.1.1.2 | Patient Characteristics (TID 5001) .....                      | 91 |
| 9.1.1.3 | OB-GYN Summary Section (TID 5002).....                        | 92 |
| 9.1.1.4 | OB-GYN Fetal Biometry Ratio Section (TID 5004).....           | 93 |
| 9.1.1.5 | OB-GYN Fetal Biometry Section (TID 5005).....                 | 93 |
| 9.1.1.6 | OB-GYN Fetal Long Bones Section (TID 5006).....               | 95 |
| 9.1.1.7 | OB-GYN Fetal Cranium Section (TID 5007).....                  | 96 |
| 9.1.1.8 | OB-GYN Early Gestation Section (TID 5011) .....               | 97 |

|          |  |     |
|----------|--|-----|
| 9.1.1.9  | <i>OB-GYN Fetal Biophysical Profile Section (TID 5009)</i>                       | 98  |
| 9.1.1.10 | <i>OB-GYN Amniotic Sac Section (TID 5010)</i>                                    | 99  |
| 9.1.1.11 | <i>OB-GYN Pelvis and Uterus Section (TID 5015)</i>                               | 100 |
| 9.1.1.12 | <i>OB-GYN Ovaries Section (TID 5012)</i>   | 101 |
| 9.1.1.13 | <i>OB-GYN Follicles Section (TID 5013)</i>                                       | 102 |
| 9.1.1.14 | <i>OB-GYN Cyst Section Section (TID SM99003)</i>                                 | 103 |
| 9.1.1.15 | <i>OB-GYN Fetal Vascular Measurement Group (TID 5025)</i>                        | 104 |
| 9.1.1.16 | <i>OB-GYN Pelvic Vascular Measurement Group (TID 5026)</i>                       | 104 |
| 9.1.1.17 | <i>OB-GYN Mass and Flow Section (TID SM99005)</i>                                | 105 |
| 9.1.1.18 | <i>OB-GYN User Creation Group Section (TID SM99010)</i>                          | 106 |
| 9.1.2    | <i>DCMR Context Groups used in V7</i>  | 108 |
| 9.1.2.1  | <i>Standard Extended Context Groups in OB-GYN SR</i>                             | 108 |
| 9.1.2.2  | <i>Gestational Age Equations and Tables (Context Group 12013)</i>                | 115 |
| 9.1.2.3  | <i>OB Fetal Body Weight Equations and Tables (Context ID 12014)</i>              | 116 |
| 9.1.2.4  | <i>Fetal Growth Equations and Tables (Context ID 12015)</i>                      | 117 |
| 9.1.2.5  | <i>Estimated Fetal Weight Percentile Equations and Tables (Context ID 12016)</i> | 118 |
| 9.2      | <i>VASCULAR ULTRASOUND REPORT TEMPLATES</i>                                      | 119 |
| 9.2.1    | <i>Vascular Ultrasound Report (TID 5100)</i>                                     | 119 |
| 9.2.1.1  | <i>Observation ConText (TID 1001)</i>  | 124 |
| 9.2.1.2  | <i>Vascular Patient Characteristics (TID 5101)</i>                               | 124 |
| 9.2.1.3  | <i>Vascular Procedure Summary Section (TID 5102)</i>                             | 125 |
| 9.2.1.4  | <i>Vascular Ultrasound Section (TID 5103)</i>                                    | 125 |
| 9.2.2    | <i>Vascular Measurement and Calculation used in Vascular SR</i>                  | 126 |
| 9.3      | <i>ADULT ECHOCARDIOGRAPHY PROCEDURE REPORT TEMPLATES</i>                         | 141 |
| 9.3.1    | <i>Adult Echocardiography Procedure Report (TID 5200)</i>                        | 141 |
| 9.3.1.1  | <i>Adult Echo User Creation Group Section (TID SM99210)</i>                      | 143 |
| 9.3.2    | <i>Cardiac Measurement and Calculation used in Adult Echocardiography SR</i>     | 144 |
| 9.4      | <i>UROLOGY STRUCTURED REPORT TEMPLATE</i>  | 174 |
| 9.4.1    | <i>Urology Ultrasound Report Templates (TID SM99400)</i>                         | 174 |
| 9.4.1.1  | <i>Observation Context (TID 1001)</i>  | 175 |
| 9.4.1.2  | <i>Patient Characteristics (TID 5001)</i>  | 175 |
| 9.4.1.3  | <i>LHW Volume Group (TID 5016)</i>   | 175 |
| 9.4.1.4  | <i>Residual Volume Group (TID SM99401)</i>                                       | 176 |
| 9.4.1.5  | <i>PSA Measurement Group (TID SM99402)</i>                                       | 176 |
| 9.4.1.6  | <i>Vascular Measurement Group (TID 5104)</i>                                     | 177 |
| 9.4.1.7  | <i>Urology User Creation Group Section (TID SM99410)</i>                         | 177 |

|          |   |     |
|----------|---|-----|
| 9.4.2    | <i>Urology Measurement and Calculation used in Urology SR</i> .....         | 177 |
| 9.5      | <b>SMALL PARTS STRUCTURED REPORT TEMPLATE</b> .....                         | 179 |
| 9.5.1    | <i>SmallParts Ultrasound Report Templates (TID SM99500)</i> .....           | 179 |
| 9.5.1.1  | <i>Observation Context (TID 1001)</i> .....                                 | 179 |
| 9.5.1.2  | <i>Patient Characteristics (TID 5001)</i> .....                             | 180 |
| 9.5.1.3  | <i>Small Parts Measruement Group (TID SM99501)</i> .....                    | 180 |
| 9.5.1.4  | <i>Mass Measruement Group (TID SM99102)</i> .....                           | 181 |
| 9.5.1.5  | <i>LHW Volume Group (TID 5016)</i> .....                                    | 182 |
| 9.5.1.6  | <i>Vascular Measurement Group (TID 5104)</i> .....                          | 182 |
| 9.5.1.7  | <i>Small Parts User Creation Group Section (TID SM99510)</i> .....          | 182 |
| 9.5.2    | <i>SmallParts Measurement and Calculation used in Small Parts SR</i> .....  | 183 |
| 9.5.2.1  | <i>Thyroid Measurement and Calculation used in Small Parts SR</i> .....     | 183 |
| 9.5.2.2  | <i>Testis Measurement and Calculation used in Small Parts SR</i> .....      | 185 |
| 9.5.2.3  | <i>Superficial Measurement and Calculation used in Small Parts SR</i> ..... | 186 |
| 9.6      | <b>BREAST IMAGING STRUCTURED REPORT TEMPLATE</b> .....                      | 189 |
| 9.6.1    | <i>Breast Imaging Report (TID 4200)</i> .....                               | 189 |
| 9.6.1.1  | <i>Language of Content Item and Descendants (TID 1204)</i> .....            | 189 |
| 9.6.1.2  | <i>Breast Imaging Report Narrative (TID 4202)</i> .....                     | 189 |
| 9.6.1.3  | <i>Observer Context (TID 1002)</i> .....                                    | 190 |
| 9.6.1.4  | <i>Breast Imaging Report Supplementary Data (TID 4208)</i> .....            | 190 |
| 9.6.1.5  | <i>Breast Imaging Report Finding Section (TID 4206)</i> .....               | 191 |
| 9.6.1.6  | <i>Breast Imaging Procedure Reported (TID 4201)</i> .....                   | 193 |
| 9.6.1.7  | <i>Breast Imaging Assessment (TID 4203)</i> .....                           | 193 |
| 9.6.1.8  | <i>Linear Measurement Template (TID 1400)</i> .....                         | 193 |
| 9.6.1.9  | <i>Area Measurement Template (TID 1401)</i> .....                           | 194 |
| 9.6.1.10 | <i>Volume Measurement Template (TID 1402)</i> .....                         | 194 |
| 9.6.1.11 | <i>Breast User Creation Group Section (TID SM99610)</i> .....               | 194 |
| 9.6.2    | <i>Standard Extended Context Groups in Breast SR</i> .....                  | 195 |
| 9.7      | <b>THYROID IMAGING STRUCTURED REPORT TEMPLATE</b> .....                     | 199 |
| 9.7.1    | <i>Thyroid Imaging Report (TID SM99700)</i> .....                           | 199 |
| 9.7.1.1  | <i>Observer Context (TID 1002)</i> .....                                    | 199 |
| 9.7.1.2  | <i>Thyroid Imaging Finding Section (TID SM99701)</i> .....                  | 200 |
| 9.7.2    | <i>Standard Extended Context Groups in Thyroid SR</i> .....                 | 202 |
| 9.8      | <b>FETAL ECHO STRUCTURED REPORT TEMPLATE</b> .....                          | 205 |
| 9.8.1    | <i>FetalEcho Ultrasound Report (TID 5220)</i> .....                         | 205 |
| 9.8.1.1  | <i>Observation Context (TID 1001)</i> .....                                 | 205 |

|         |  |     |
|---------|--|-----|
| 9.8.1.2 | <i>Cardiovascular Patient Characteristics (TID 3602)</i> .....                     | 206 |
| 9.8.1.3 | <i>Cardiac Ultrasound Fetal Measurement Section (TID 5228)</i> .....               | 206 |
| 9.8.1.4 | <i>Pediatric, Fetal and Congenital Cardiac Ultrasound Section (TID 5222)</i> ..... | 209 |
| 9.8.2   | <i>Standard Extended Context Groups in FetalEcho SR</i> .....                      | 209 |
| 9.9     | <b>PEDIATRIC STRUCTURED REPORT TEMPLATE</b> .....                                  | 217 |
| 9.9.1   | <i>Pediatric Ultrasound Report (TID SM99800)</i> .....                             | 217 |
| 9.9.1.1 | <i>Pediatric User Creation Group Section (TID SM99810)</i> .....                   | 217 |
| 9.9.2   | <i>Standard Extended Context Groups in Pediatric SR</i> .....                      | 218 |
| 9.10    | <b>MUSCULOSKELETAL STRUCTURED REPORT TEMPLATE</b> .....                            | 218 |
| 9.10.1  | <i>Musculoskeletal Ultrasound Report (TID SM99900)</i> .....                       | 218 |

## 3 INTRODUCTION

### 3.1 REVISION HISTORY

| Document Version | System Version | Date of Issue | Author                    | Description                |
|------------------|----------------|---------------|---------------------------|----------------------------|
| 1.00             | 1.00           | Jul 18 2021   | SAMSUNG MEDISON CO., LTD. | Final Text for System 1.0  |
| 2.00             | 1.01           | Aug 16 2022   | SAMSUNG MEDISON CO., LTD. | Final Text for System 1.01 |
|                  |                |               |                           |                            |
|                  |                |               |                           |                            |

### 3.2 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

### 3.3 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with SAMSUNG MEDISON CO., LTD and other vendor's Medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. However, it is not guaranteed to ensure by itself the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between SAMSUNG MEDISON CO., LTD and non – SAMSUNG MEDISON CO., LTD equipment.
- Test procedures should be defined to validate the desired level of connectivity.
- The DICOM Standard will evolve to meet the users' future requirements. SAMSUNG MEDISON CO., LTD is activity involved in developing the standard further and therefore reserves the right to make changes to its products

or to discontinue their delivery.

### **3.4 DEFINITIONS, TERMS AND ABBREVIATIONS**

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM Standard.

Abbreviations and terms are as follows:

|      |   |
|------|---|
| AE   | DICOM Application Entity                    |
| AET  | Application Entity Title                    |
| ASCE | Association Control Service Element         |
| CD-R | Compact Disk Recordable                     |
| FSC  | File-Set Creator                            |
| FSU  | File-Set Updater                            |
| FSR  | File-Set Reader                             |
| IOD  | (DICOM) Information Object Definition       |
| ISO  | International Standard Organization         |
| MPPS | Modality Performed Procedure Step           |
| MSPS | Modality Scheduled Procedure Step           |
| Q/R  | Query and Retrieve                          |
| R    | Required Key Attribute                      |
| O    | Optional Key Attribute                      |
| PDU  | DICOM Protocol Data Unit                    |
| SCU  | DICOM Service Class User (DICOM client)     |
| SCP  | DICOM Service Class Provider (DICOM server) |
| SOP  | DICOM Service-Object Pair                   |
| U    | Unique Key Attribute                        |


### **3.5 REFERENCES**

[DICOM] Digital Imaging and Communications in Medicine (DICOM), NEMA PS 3.1–3.20, 2011

## 4 NETWORKING

### 4.1 IMPLEMENTATION MODEL

#### 4.1.1 Application Data Flow



**Figure 4.1-1**  
**APPLICATION DATA FLOW DIAGRAM**

- The Storage Application Entity sends images, Structured Reports and requests Storage Commitment to a remote AE. It is associated with the local real-world activities "Send Images & Reports" and "Request Commitments". Methods to send SOP Instances(images and Structured Reports) depend on user

configuration, “Send on end exam”, “Send after acquisition” or “Send manually”. “Send manually” mode is performed upon user request for each study or for specific images selected. “Send on end exam” mode starts to send SOP Instances at End Exam for each study. “Send after acquisition” mode starts when the first SOP Instance is acquired for each study and SOP Instances are transferred immediately after acquisition.

If the remote AE is configured as an archive device, the Storage AE will request Storage Commitment and if a commitment is successfully obtained, it will record this information in the local database and displayed it in the Exam List.

- The Workflow Application Entity receives Worklist information from and sends MPPS information to a remote AE. It is associated with the local real-world activities “Update Worklist” and “Acquire Images”. When the “Update Worklist” local real-world activity is performed, the Workflow Application Entity queries a remote AE for worklist items and provides the set of worklist items matching the query request. “Update Worklist” is performed as a result of an operator request or can be performed automatically at specific time intervals. When the “Acquire Images” local real-world activity is performed, the Workflow Application Entity creates and updates Modality Performed Procedure Step instances managed by a remote AE. Acquisition of images will result in automated creation of an MPPS Instance. Completion of the MPPS is performed at End Exam for each study.
- The Hardcopy Application Entity prints images on a remote AE (Printer). It is associated with the local real-world activity “Film Images”. Methods to film Images depend on user configuration and are equal to the Sending images’ of the Storage Application Entity.

## **4.1.2 Functional Definition of AE’s**

### **4.1.2.1 Functional Definition of Storage Application Entity**

The existence of a send-job with associated network destination will activate the Storage AE. An association request is sent to the destination AEs and upon successful negotiation of a Presentation Context, the image or Structured Report transfer is started. If the association cannot be opened, the related send-job is set to an error state and can be restarted by the user via DICOM Spooler interface or automatically. An automatic retry (retry interval, retry count) can be configured using the Setup/DICOM Menu.

### **4.1.2.2 Functional Definition of Workflow Application Entity**

Worklist Update attempts to download a Worklist from a remote node. If the Workflow AE establishes an

association to a remote AE, it will transfer all matching worklist items via the open Association. By default, Worklist Update use "US" for Modality, current date for Scheduled Procedure Step Start Date and blank (displayed "Any") for Scheduled Station AE-Title as query parameters. The results will be displayed in a separate list, which will be cleared with the next Worklist Update.

The Workflow AE performs the creation of an MPPS Instance automatically whenever the first SOP Instance is acquired for each study. The MPPS "Complete" or "Discontinued" states can only be set by "End Exam" for each study.

#### **4.1.2.3 Functional Definition of Hardcopy Application Entity**

The existence of a print-job will activate the Hardcopy AE. An association is established with the printers and the printer's status determined. If the printer is operating normally, the film sheets described within the print-job will be printed. If the printer is not operating normally, the print-job will set to an error state and can be restarted by the user via DICOM Spooler manager interface or automatically. An automatic retry (retry interval, retry count) can be configured using the Setup/DICOM Menu.


#### **4.1.2.4 Functional Definition of the Q/R Application Entity**

The Query function is activated through the user interface when the user selects a Q/R AE to query (from a pre-configured list), then initiates a query. Queries are performed per the study or series. Retrieval is activated through the user interface when the user selects a series for retrieval. A connection to the Q/R AE is established to initiate and monitor the retrieval and the STORAGE-SCP AE receives the retrieved instances.


#### **4.1.2.5 Functional Definition of the Q/R Application Entity**

The STORAGE-SCP AE waits for another application to connect from the presentation address configured for its AE Title. When another application connects, the STORAGE-SCP AE expects it to be a DICOM application. The STORAGE-SCP AE will accept associations with presentation contexts for SOP Classes of the Storage Device. Any images received in such Presentation Contexts will be stored in the system. Provide Storage SCP only Q/R service running.


### 4.1.3 Sequencing of Real-World Activities



**Figure 4.1-2**  
**SEQUENCING CONSTRAINTS – SEND AFTER ACQUISITION**



**Figure 4.1-3**  
**SEQUENCING CONSTRAINTS – SEND ON END EXAM**



**Figure 4.1-4**  
**SEQUENCING CONSTRAINTS – QUERY AND RETRIEVE**

Under normal scheduled workflow conditions, the sequencing constraints are illustrated in Figure 4.1-2, Figure 4.1-3 and Figure 4.1-4

Other workflow situations (e.g. unscheduled procedure steps) will have other sequencing constraints. Printing could equally take place after the images acquired have been stored. Printing could be omitted completely if no printer is connected or hardcopies are not required.

## 4.2 AE SPECIFICATIONS

### 4.2.1 Storage Application Entity Specification

#### 4.2.1.1 SOP Classes

V7 provides Standard Conformance to the following SOP Classes:

**Table 4.2-1**  
**SOP CLASSES FOR AE STORAGE**

| SOP Classes                             | SOP Class UID                 | SCU | SCP |
|---|-------------------------------|-----|-----|
| Ultrasound Image Storage                | 1.2.840.10008.5.1.4.1.1.6.1   | Yes | No  |
| Ultrasound Multi-frame Image Storage    | 1.2.840.10008.5.1.4.1.1.3.1   | Yes | No  |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.33 | Yes | No  |
| Storage Commitment Push Model           | 1.2.840.10008.1.20.1          | Yes | No  |
| Verification                            | 1.2.840.10008.1.1             | Yes | Yes |

#### 4.2.1.2 Association Policies

##### 4.2.1.2.1 General

The DICOM Standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-2**  
**DICOM APPLICATION CONTEXT FOR AE STORAGE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 4.2.1.2.2 Number of Associations

V7 can initiate one or more Associations at a time for each destination to which a transfer request is being processed in the active job queue list.

**Table 4.2-3**  
**NUMBER OF ASSOCIATIONS INITIATED FOR AE STORAGE**

|   |           |
|---|-----------|
| Maximum number of simultaneous Associations | Unlimited |
|---|-----------|

V7 accepts Associations to receive N-EVENT\_REPORT notifications for the Storage Commitment Push Model SOP

Class.

**Table 4.2-4**  
**NUMBER OF ASSOCIATIONS ACCEPTED FOR AE STORAGE**

|   |           |
|---|-----------|
| Maximum number of simultaneous Associations | Unlimited |
|---|-----------|

**4.2.1.2.3 Asynchronous Nature**

V7 does not support asynchronous communications (multiple outstanding transactions over a single Association).

**Table 4.2-5**  
**ASYNCHRONOUS NATURE AS A SCU FOR AE STORAGE**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

**4.2.1.2.4 Implementation Identifying Information**

The implementation information for this Application Entity is:

**Table 4.2-6**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE STORAGE**

|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

**4.2.1.3 Association Initiation Policy**

**4.2.1.3.1 Activity – Send Images and Structured Reports and Requests Commitment**


**4.2.1.3.1.1 Description and Sequencing of Activities**

A user can select exams or images and request them to be sent to some destination. Each request is forwarded to the job queue and processed individually. When the “Send on end exam” or “Send after acquisition” option is active, Stored images and reports will be forwarded to the network job queue for a pre-configured auto-send target destination automatically. For “Send on end exam” and “Send manually” configuration, the system opens an association, sends all images in the study, and closes the association. If “Send after acquisition” is selected, the system handles the association with the Storage SCP Server using the following method.

- a. Open an Association when the image is acquired.

- b. If an error occurs while sending an SOP Instance to the server because there is no longer an open association (server timed-out), attempt to re-establish the association.
- c. When one image is sended, close the open association after SOP Instances remained in that study are sent.

If the remote AE is configured as an archive device, the Storage AE will, after all images and reports have been sent, transmit Storage Commitment request (N-ACTION) over a separate Association. The Storage AE can only receive an N-EVENT-REPORT request in a subsequent association initiated by the SCP.



**Figure 4.2-1**  
**SEQUENCING OF ACTIVITY - SEND IMAGES AND SR DOCUMENTS**

A possible sequence of interactions between the Storage AE and an Image Manager (e.g. a storage or archive device supporting the Storage and Storage Commitment SOP Classes as an SCP) is illustrated in the figure above.

NOTE: The N-EVENT-REPORT must be sent over a separate association initiated by the Image Manager. (See Section 4.2.1.4)

#### 4.2.1.3.1.2 Proposed Presentation Contexts

V7 is capable of proposing the Presentation Contexts shown in the following table.

**Table 4.2-7  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND IMAGES**

| Presentation Context Table              |                                   |                           |                        |      |           |
|---|-----------------------------------|---------------------------|------------------------|------|-----------|
| Abstract Syntax                         |                                   | Transfer Syntax           |                        | Role | Ext. Neg. |
| Name                                    | UID                               | Name List                 | UID List               |      |           |
| Ultrasound Image Storage                | 1.2.840.10008.5.1                 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU  | None      |
|   | .4.1.1.6.1                        | JPEG Lossy Baseline       | 1.2.840.10008.1.2.4.50 |      |           |
|   |                                   | JPEG Lossless             | 1.2.840.10008.1.2.4.70 |      |           |
| Ultrasound Multi-frame Image Storage    | 1.2.840.10008.5.1                 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU  | None      |
|   | .4.1.1.3.1                        | JPEG Lossy Baseline       | 1.2.840.10008.1.2.4.50 |      |           |
|   |                                   | JPEG Lossless             | 1.2.840.10008.1.2.4.70 |      |           |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.1<br>.4.1.1.88.33 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU  | None      |
| Storage Commitment Push Model           | 1.2.840.10008.1.2                 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU  | None      |
|   | 0.1                               | Explicit VR Little Endian | 1.2.840.10008.1.2.1    |      |           |
| Verification                            | 1.2.840.10008.1.1                 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU  | None      |
|   |                                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    |      |           |

Presentation Contexts for Ultrasound Image Storage and Ultrasound Multi-frame Image Storage will be proposed for the “STORAGE” device configured in Setup/Connectivity/DICOM.

A Presentation Context for Comprehensive Structured Report Storage will be proposed for the “STORAGE SR” device configured in Setup/ Connectivity/DICOM.

A Presentation Context for Storage Commitment Push Model will be proposed for the “SC” device configured in Setup/ Connectivity/DICOM.

A Presentation Context for Verification will be proposed when a user press the “Verify” button for a configured device.

#### 4.2.1.3.1.3 SOP Specific Conformance Image & Comprehensive Structured Report Storage SOP Classes

All Image & Structured Report Storage SOP Classes supported by the Storage AE exhibit the same behavior, except where stated, and are described together in this section.

**Table 4.2-8  
STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning                   | Error Code             | Behavior   |
|----------------|-----------------------------------|------------------------|--|
| Success        | Success                           | 0000                   | The SCP has successfully stored the SOP Instance. If all SOP Instances succeed, the job is marked as complete. |
| Refused        | Out of Resources                  | A700-A7FF              | The association is aborted using A-ABORT and the send job is marked as failed. The status is logged.           |
| Error          | Data Set does not match SOP Class | A900-A9FF              | Same as "Refused" above.   |
| Error          | Cannot Understand                 | C000-CFFF              | Same as "Refused" above.   |
| Warning        | Coercion of Data Elements         | B000                   | Image transmission is considered successful.   |
| Warning        | Data Set does not match SOP Class | B007                   | Same as "Warning" above.   |
| Warning        | Elements Discards                 | B006                   | Same as "Warning" above.   |
| *              | *                                 | Any other status code. | Same as "Refused" above.   |

The Behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-9  
STORAGE COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. |
| Association aborted by the SCP or network layers | The Send job is marked as failed.  |

A failed send job can be restarted by user interaction. The system can be configured to automatically resend failed jobs if a transient status code is received. The delay between resending failed jobs and the number of retries is also configurable.

**4.2.1.3.1.4 SOP Specific Conformance for Storage Commitment SOP Class**

**4.2.1.3.1.4.1 Storage Commitment Operations (N-ACTION)**

The Storage AE will request storage commitment for the configured device for instances of the Ultrasound Image, Ultrasound Multi-frame Image and Structured Report Storage SOP Classes.

The Storage AE will consider Storage Commitment failed if no N-EVENT-REPORT is received for a Transaction UID within a configurable time period after receiving a successful N-ACTION response (duration of applicability for a Transaction UID).

The Storage AE does not send the optional Storage Media FileSet ID & UID Attributes or the Referenced Study Component Sequence Attribute in the N-ACTION

The Behavior of Storage AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 4.2-10  
STORAGE COMMITMENT N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Error Code</b>      | <b>Behavior</b>  |
|-----------------------|------------------------|------------------------|--|
| Success               | Success                | 0000                   | The request for storage comment is considered successfully sent. The system waits for the association of the N-Event-Report. |
| *                     | *                      | Any other status code. | The Association is aborted using A-Abort and the request for storage comment is marked as failed                             |

The behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-11  
STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>  |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the storage commitment job is marked as failed. |
| Association aborted by the SCP or network layers | The storage commitment job is marked as  |

|  |         |
|--|---------|
|  | failed. |
|--|---------|

#### 4.2.1.3.1.4.2 Storage Commitment Notification (N-EVENT-REPORT)

The Storage AE is capable of receiving an N-EVENT-REPORT notification if it has successfully negotiated a Presentation Context for the Storage Commitment Push Model.

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be cancelled.

The behavior of Storage AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below.

**Table 4.2-12  
STORAGE COMMITMENT N-EVENT-REPORT BEHAVIOR**

| Event Type Name                                       | Event Type ID | Behavior   |
|---|---------------|--|
| Storage Commitment Request Successful                 | 1             | The commit status is set to “Y” for each exam in the exam list.<br>Auto deletion for committed exam is not supported.  |
| Storage Commitment Request Complete – Failures Exists | 2             | The commit status is set to “N” for each exam in the exam list.<br>The Referenced SOP Instances under Failed SOP Sequence (0008, 1198) are logged. A send job that failed storage commitment will not be automatically restarted but can be restarted by user interaction. |

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the Table below.

**Table 4.2-13  
STORAGE COMMITMENT N-EVENT-REPORT RESPONSE STATUS REASONS**

| Service Status | Further Meaning        | Error Code | Behavior   |
|----------------|------------------------|------------|--|
| Success        | Success                | 0000       | The Storage commitment result has been successfully received.  |
| Failure        | Unrecognized Operation | 0211H      | The Transaction UID in the N_EVENT_REPORT request is not (was never issued within an N_ACTION request) |

|         |                    |       |   |
|---------|--------------------|-------|---|
| Failure | No Such Event Type | 0113H | An invalid Event Type ID was supplied in the N_EVENT_REPORT request |
| Failure | Processing Failure | 0110H | An internal error occurred during processing of the N_EVENT_REPORT  |

#### 4.2.1.3.1.5 SOP Specific Conformance for Verification

The Behavior when encountering status codes in a C-ECHO response is summarized in the Table below:

**Table 4.2-14  
VERIFICATION C-ECHO RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code            | Behavior                                |
|----------------|-----------------|-----------------------|---|
| Success        | Success         | 0000                  | Verification Status is set to 'Success' |
| *              | *               | Any other status code | Verification Status is set to 'Failed'  |

The Behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-15  
VERIFICATION COMMUNICATION FAILURE BEHAVIOR**


| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the verification job is marked as failed. |
| Association aborted by the SCP or network layers | The verification job is marked as failed.  |

#### 4.2.1.4 Association Acceptance Policy

##### 4.2.1.4.1 Activity – Receive Storage Commitment Response

##### 4.2.1.4.1.1 Description and Sequence of Activities

The Storage AE will accept associations in order to receive responses to a Storage Commitment Request.



**Figure 4.2-2**  
**SEQUENCING OF ACTIVITY - RECEIVE STORAGE COMMITMENT RESPONSE**

A possible sequence of interactions between the Storage AE and an Image Manager (e.g. a storage or archive device supporting Storage Commitment SOP Classes as an SCP) is illustrated in the Figure above:

1. The Image Manager opens a new association with the Storage AE.
2. The Image Manager sends an N-EVENT-REPORT request notifying the Storage AE of the status of a previous Storage Commitment Request. The Storage AE replies with an N-EVENT-REPORT response confirming receipt.
3. The Image Manager closes the association with the Storage AE.

**4.2.1.4.1.2 Accepted Presentation Contexts**

The Storage AE will accept Presentation Contexts as shown in the Table below.

**Table 4.2-16**  
**ACCEPTABLE PRESENTATION CONTEXTS FOR ACTIVITY**  
**RECEIVE STORAGE COMMITMENT RESPONSE**

| Presentation Context Table |     |                 |          |      |           |
|----------------------------|-----|-----------------|----------|------|-----------|
| Abstract Syntax            |     | Transfer Syntax |          | Role | Ext. Neg. |
| Name                       | UID | Name List       | UID List |      |           |
|                            |     |                 |          |      |           |

|                               |                      |  |  |     |      |
|-------------------------------|----------------------|--|--|-----|------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU | None |
| Verification                  | 1.2.840.10008.1.1    | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP | None |

#### 4.2.1.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

##### 4.2.1.4.1.3.1 Storage Commitment Notifications (N-EVENT-REPORT)

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be cancelled.

The behavior of Storage AE when receiving Event Types within the N-EVENT\_REPORT is summarized in Table 4.2-12.

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in Table 4.2-13.

#### 4.2.1.4.1.4 SOP Specific Conformance for Verification SOP Class

The Storage AE provides standard conformance to the Verification SOP Class as an SCP. If the C-ECHO request was successfully received, a 0000 (Success) status code will be returned in the C-ECHO response.

## 4.2.2 Workflow Application Entity Specification

### 4.2.2.1 SOP Classes

V7 provides Standard Conformance to the following SOP Classes:

**Table 4.2-17**  
**SOP CLASSES FOR AE WORKFLOW**

| SOP Classes                                | SOP Class UID           | SCU | SCP |
|--|-------------------------|-----|-----|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31  | Yes | No  |
| Modality Performed Procedure Step          | 1.2.840.10008.3.1.2.3.3 | Yes | No  |

## 4.2.2.2 Association Establishment Policy

### 4.2.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed.

**Table 4.2-18**

#### **DICOM APPLICATION CONTEXT FOR AE WORKFLOW**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

### 4.2.2.2.2 Number of Associations

V7 initiates one Association at a time for a Worklist request.

**Table 4.2-19**

#### **NUMBER OF ASSOCIATIONS INITIATED FOR AE WORKFLOW**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

### 4.2.2.2.3 Asynchronous Nature

V7 does not support asynchronous communications (multiple outstanding transactions over a single Association)

**Table 4.2-20**

#### **ASYNCHRONOUS NATURE AS A SCU FOR AE WORKFLOW**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

### 4.2.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-21**

#### **DICOM IMPLEMENTATION CLASS AND VERSION FOR AE WORKFLOW**

|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

### **4.2.2.3 Association Initiation Policy**

#### **4.2.2.3.1 Activity – Worklist Update**

##### **4.2.2.3.1.1 Description and Sequencing of Activities**


The request for a Worklist Update is initiated by user interaction.

The interactive Worklist Query will display a dialog for entering data as search criteria. When the Query is started on your request, only the data from the dialog will be inserted as matching keys into the query.

With automated worklist queries the V7 always requests all items for a Scheduled Procedure Step Start Date (actual date), Modality (US) and Scheduled Station AE Title.

Upon initiation of the request, the V7 will build an Identifier for the C-FIND request, will initiate an Association to send the request and will wait for Worklist responses. After retrieval of all responses, V7 will access the local database to add patient demographic data. The results will be displayed in a separate list, which will be cleared with the next worklist update.

V7 will initiate an Association in order to issue a C-FIND request according to the Modality Worklist Information Model.



**Figure 4.2-3**  
**SEQUENCING OF ACTIVITY - WORKLIST UPDATE**

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as an SCP) is illustrated in the figure above:

**4.2.2.3.1.2 Proposed Presentation Contexts**

V7 will propose Presentation Contexts as shown in the following table:

**Table 4.2-22**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY WORKLIST UPDATE**

| Presentation Context Table                 |                |                           |                      |      |           |
|--|----------------|---------------------------|----------------------|------|-----------|
| Abstract Syntax                            |                | Transfer Syntax           |                      | Role | Ext. Neg. |
| Name                                       | UID            | Name List                 | UID List             |      |           |
| Modality Worklist Information Model - FIND | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2    | SCU  | None      |
|  | 5.1.4.31       | Explicit VR Little Endian | 1.2.840.10008. 1.2.1 |      |           |

#### 4.2.2.3.1.3 SOP Specific Conformance for Modality Worklist

The behavior of V7 when encountering status codes in a Modality Worklist C-FIND response is summarized in the Table below. If any other SCP response status than “Success” or “Pending” is received by V7, a message “Query failed” will appear on the user interface.

**Table 4.2-23  
MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Error Code</b>            | <b>Behavior</b>  |
|-----------------------|---|------------------------------|--|
| Success               | Matching is complete  | 0000                         | The SCP has Completed the operation successfully.                                |
| Pending               | Matches are continuing  | FF00                         | Continue.  |
| Pending               | Matches are continuing -<br>Warning that one or more<br>Optional Keys were not<br>supported | FF01                         | Continue.  |
| *                     | *   | Any other<br>status<br>code. | The Association is aborted using A-Abort and the<br>Worklist is marked as failed |

The behavior of V7 during communication failure is summarized in the Table below.

**Table 4.2-24  
MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>   |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and<br>the worklist query is marked as failed. |
| Association aborted by the SCP or network layers | The Worklist query is marked as failed.   |

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The Table below provides a description of the V7 Worklist Request Identifier and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

Requested return attributes not supported by the SCP are set to have no value. Non-matching responses returned by the SCP due to unsupported optional matching keys are ignored. No attempt is made to filter out possible duplicate entries.

**Table 4.2-25**  
**WORKLIST REQUEST IDENTIFIER**

| Module Name<br>Attribute Name           | Tag       | VR | M   | R | Q | D | IOD |
|---|-----------|----|-----|---|---|---|-----|
| Scheduled Procedure Step                |           |    |     |   |   |   |     |
| Scheduled Procedure Step Sequence       | 0040,0100 | SQ |     | x |   |   |     |
| > Scheduled Station AET                 | 0040,0001 | AE | (S) | x | x |   |     |
| > Scheduled Procedure Step Start Date   | 0040,0002 | DA | S,R | x | x | x |     |
| > Scheduled Procedure Step Start Time   | 0040,0003 | TM |     | x |   |   |     |
| > Modality                              | 0008,0060 | CS | S   | x | x |   |     |
| > Scheduled Performing Physician's Name | 0040,0006 | PN |     | x |   | x |     |
| > Scheduled Procedure Step Description  | 0040,0007 | LO |     | x |   | x | x   |
| > Scheduled Station Name                | 0040,0010 | SH | S   | x | x |   |     |
| > Scheduled Procedure Step Location     | 0040,0011 | SH |     | x |   |   |     |
| > Requested Contrast Agent              | 0032,1070 | LO |     | x |   |   |     |
| > Scheduled Protocol Code Sequence      | 0040,0008 | SQ |     | x |   |   | x   |
| > Scheduled Procedure Step ID           | 0040,0009 | SH |     | x |   |   | x   |
| Requested Procedure                     |           |    |     |   |   |   |     |
| Requested Procedure ID                  | 0040,1001 | SH | S   | x | x |   | x   |
| Requested Procedure Description         | 0032,1060 | LO |     | x |   |   |     |
| Study Instance UID                      | 0020,000D | UI |     | x |   |   | x   |
| Requested Procedure Comments            | 0040,1400 | LT |     | x |   |   |     |
| Referenced Study Sequence               | 0008,1110 | SQ |     | x |   |   |     |
| Requested Procedure Code Sequence       | 0032,1064 | SQ |     | x |   |   |     |
| Names of Intended Recipients of Results | 0040,1010 | PN |     | x |   |   |     |
| Imaging Service Request                 |           |    |     |   |   |   |     |
| Accession Number                        | 0008,0050 | SH | S   | x | x | x | x   |
| Requesting Physician                    | 0032,1032 | PN |     | x |   |   |     |
| Referring Physician's Name              | 0008,0090 | PN |     | x |   | x | x   |
| Visit Status                            |           |    |     |   |   |   |     |
| Current Patient Location                | 0038,0300 | LO |     | x |   |   |     |
| Patient Identification                  |           |    |     |   |   |   |     |

|                            |           |    |   |   |   |   |   |
|----------------------------|-----------|----|---|---|---|---|---|
| Patient's Name             | 0010.0010 | PN | S | x | x | x | x |
| Patient ID                 | 0010,0020 | LO | S | x | x | x | x |
| Other Patient IDs          | 0010,1000 | LO |   | x |   |   | x |
| Other Patient IDs Sequence | 0010,1002 | SQ |   | x |   |   | x |
| > Patient ID               | 0010,0020 | LO |   | x |   | x | x |
| > Type of Patient ID       | 0010,0022 | CS |   | x |   |   | x |
| <b>Patient Demographic</b> |           |    |   |   |   |   |   |
| Patient's Birth Date       | 0010,0030 | DA |   | x |   | x | x |
| Patient's Sex              | 0010,0040 | CS |   | x |   | x | x |
| Patient's Size             | 0010,1020 | DS |   | x |   | x | x |
| Patient's Weight           | 0010,1030 | DS |   | x |   | x | x |
| Ethnic Group               | 0010,2160 | SH |   | x |   |   |   |
| Patient Comments           | 0010,4000 | LT |   | x |   |   |   |
| <b>Patient Medical</b>     |           |    |   |   |   |   |   |
| Medical Alerts             | 0010,2000 | LO |   | x |   |   |   |
| Additional Patient History | 0010,21B0 | LT |   | x |   | x |   |
| Pregnancy Status           | 0010,21C0 | US |   | x |   |   |   |
| Last Menstrual Date        | 0010,21D0 | DA |   | x |   | x | x |

The above table should read as follows:

- Module Name: The Name of the associated module for supported worklist attributes.
- Attribute Name: Attributes supported to build an V7 Worklist Request Identifier.
- Tag: DICOM tag for this attribute.
- VR: DICOM VR for this attribute.
- M: Matching keys for (automatic) Worklist Update. An "S" indicates that V7 supplies an attribute value for Single Value Matching or additional specific tags indicated by "(S)"; an "R" will indicate Range Matching.
- R: Return keys. An "X" will indicate that V7 will supply this attribute as Return Key with zero length for Universal Matching.
- Q: Interactive Query Key. An "X" will indicate that V7 will supply this attribute as matching key, if entered in the Setup Dialog.
- D: Displayed keys. An "X" indicates that this worklist attribute is displayed to the user during a patient registration dialog.
- IOD: An "X" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

#### **4.2.2.3.2 Activity – Acquire Images**

##### **4.2.2.3.2.1 Description and Sequencing of Activities**

An Association to the configured MPPS SCP system is established immediately after the first SOP Instance is acquired to send the MPPS N-Create message.


The “End Exam” button causes a message box in which a user can select “COMPLETED” or “DISCONTINUED” as a MPPS final state. An exam for which an MPPS instance is sent with a state of “COMPLETED” or “DISCONTINUED” can no longer be updated.

The V7 will support creation of “unscheduled cases” by allowing MPPS Instances to be communicated for locally registered Patients.

The V7 supports a 1-to-N relationship between Scheduled and Performed Procedure Steps.

V7 will initiate an Association to issue an:

- N-CREATE request according to the CREATE Modality Performed Procedure Step SOP Instance operation, or an:
- N-SET request to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.



**Figure 4.2-4  
SEQUENCING OF ACTIVITY - ACQUIRE IMAGES**

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the MPPS SOP Class as an SCP) is illustrated in the figure above:

#### 4.2.2.3.2.2 Proposed Presentation Contexts

V7 will propose Presentation Contexts as shown in the following table:

**Table 4.2-26**

**PROPOSED PRESENTATION CONTEXTS FOR REAL-WORLD ACTIVITY ACQUIRE IMAGES**

| Presentation Context Table |                |                           |                      |      |           |
|----------------------------|----------------|---------------------------|----------------------|------|-----------|
| Abstract Syntax            |                | Transfer Syntax           |                      | Role | Ext. Neg. |
| Name                       | UID            | Name List                 | UID List             |      |           |
| Modality Performed         | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2    | SCU  | None      |
| Procedure Step             | 3.1.2.3.3      | Explicit VR Little Endian | 1.2.840.10008. 1.2.1 |      |           |

**4.2.2.3.2.3 SOP Specific Conformance for MPPS**

The behavior of V7 when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in the Table below. If any other SCP response status than “Success” or “Warning” is received by V7, a message “MPPS failed” will appear on the user interface.

**Table 4.2-27**

**MPPS N-CREATE / N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning              | Error Code             | Behavior  |
|----------------|------------------------------|------------------------|---|
| Success        | Success                      | 0000                   | The SCP has Completed the operation successfully.                         |
| Warning        | Attribute Value Out of Range | 0116H                  | The MPPS Operation is considered successful.                              |
| *              | *                            | Any other status code. | The Association is aborted using A-Abort and the MPPS is marked as failed |

The behavior of V7 during communication failure is summarized in the table below:

**Table 4.2-28**

**MPPS COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the MPPS job is marked as failed. |
| Association aborted by the SCP or network layers | The MPPS job is marked as failed.  |

Table 4.2-29 provides a description of the MPPS N-CREATE and N-SET request identifiers send by V7. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent.

**Table 4.2-29  
MPPS N-CREATE / N-SET REQUEST IDENTIFIER**

| Attribute Name                               | Tag       | VR | N-CREATE                                 | N-SET |
|--|-----------|----|--|-------|
| Specific Character Set                       | 0008,0005 | CS | Ref. Section 6 SUPPORT OF CHARACTER SETS |       |
| <b>Performed Procedure Step Relationship</b> |           |    |  |       |
| Scheduled Step Attribute Sequence            | 0040,0270 | SQ |  |       |
| > Study Instance UID                         | 0020,000D | UI | From MWL or generated by device          |       |
| > Referenced Study Sequence                  | 0008,1110 | SQ | From MWL                                 |       |
| >> Referenced SOP Class UID                  | 0008.1150 | UI | From MWL                                 |       |
| >> Referenced SOP Instance UID               | 0008,1155 | UI | From MWL                                 |       |
| > Accession Number                           | 0008,0050 | SH | From MWL or user input                   |       |
| > Requested Procedure ID                     | 0040,1001 | SH | From MWL                                 |       |
| > Requested Procedure Description            | 0032,1060 | LO | From MWL                                 |       |
| > Scheduled Procedure Step ID                | 0040,0009 | SH | From MWL                                 |       |
| > Scheduled Procedure Step Description       | 0040,0007 | LO | From MWL                                 |       |
| > Scheduled Protocol Code Sequence           | 0040,0008 | SQ | From MWL                                 |       |
| >> Code Value                                | 0008,0100 | SH | From MWL                                 |       |
| >> Coding Scheme Designator                  | 0008,0102 | SH | From MWL                                 |       |
| >> Coding Scheme Version                     | 0008,0103 | SH | From MWL                                 |       |
| >> Code Meaning                              | 0008,0104 | LO | From MWL                                 |       |
| Patient's Name                               | 0010,0010 | PN | From MWL or user input                   |       |
| Patient ID                                   | 0010,0020 | LO | From MWL or user input                   |       |
| Patient's Birth Date                         | 0010,0030 | DA | From MWL or user input                   |       |
| Patient's Sex                                | 0010,0040 | CS | From MWL or user input                   |       |

|   |            |    |   |   |
|---|------------|----|---|---|
| Referenced Patient Sequence                                   | 0008,1120  | SQ | Zero length   |   |
| > Referenced SOP Class UID                                    | 0008,1150  | UI | Zero length   |   |
| > Referenced Instance UID                                     | 0008,1155  | UI | Zero length   |   |
| <b>Performed Procedure Step Information</b>                   |            |    |   |   |
| Series Number   | 0020, 0011 | IS | Generated by device                                   |   |
| Performed Procedure Step ID                                   | 0040,0253  | SH | Generated by device<br>(Study Date + Study Time)      |   |
| Performed Station AE Title                                    | 0040,0241  | AE | From Modality Setup                                   |   |
| Performed Station Name  | 0040,0242  | SH | From Modality Setup                                   |   |
| Performed Location  | 0040,0243  | SH | Zero length   |   |
| Performed Procedure Step Start Date                           | 0040,0244  | DA | Actual Start Date                                     |   |
| Performed Procedure Step Start Time                           | 0040,0245  | TM | Actual Start Time                                     |   |
| Performed Procedure Step Status                               | 0040,0252  | CS | "IN PROGRESS"   | "COMPLETED" or<br>"DISCONTINUED"                            |
| Performed Procedure Step Description                          | 0040,0254  | LO | From MWL or user input<br>(Same as Study Description) | From MWL or user input<br>(Same as Study Description)       |
| Performed Procedure Type Description                          | 0040,0255  | LO | Zero length   | Zero length   |
| Procedure Code Sequence                                       | 0008,1032  | SQ | From MWL  | From MWL  |
| > Code Value  | 0008,0100  | SH | From MWL  | From MWL  |
| > Coding Scheme Designator                                    | 0008,0102  | SH | From MWL  | From MWL  |
| > Coding Scheme Version                                       | 0008,0103  | SH | From MWL  | From MWL  |
| > Code Meaning  | 0008,0104  | LO | From MWL  | From MWL  |
| Performed Procedure Step End Date                             | 0040,0250  | DA | Zero length   | Actual End Date   |
| Performed Procedure Step End Time                             | 0040,0251  | TM | Zero length   | Actual End Time   |
| Performed Procedure Step Discontinuation Reason Code Sequence | 0040,0281  | SQ |   | Used when Performed Procedure Step Status is "DISCONTINUED" |

|  |           |    |  |                        |
|--|-----------|----|--|------------------------|
| > Code Value   | 0008,0100 | SH |  | From User Select       |
| > Coding Scheme Designator                                   | 0008,0102 | SH |  | From User Select       |
| > Coding Scheme Version                                      | 0008,0103 | SH |  |                        |
| > Code Meaning   | 0008,0104 | LO |  | From user select       |
| <b>Image Acquisition Results</b>                             |           |    |  |                        |
| Modality   | 0008,0060 | CS | "US"   |                        |
| Study ID   | 0020,0010 | SH | Requested Procedure ID<br>or Generated by device<br>(Study Date + Study<br>Time) |                        |
| Performed Protocol Code<br>Sequence                          | 0040,0260 | SQ | Zero length or Scheduled<br>Protocol Code Sequence                               |                        |
| Performed Series Sequence                                    | 0040,0340 | SQ | Zero length  | One or more items      |
| > Performing Physician's Name                                | 0008,1050 | PN |  | From MWL or user input |
| > Protocol Name  | 0018,1030 | LO |  | "FreeForm"             |
| > Operator's Name  | 0008,1070 | PN |  | From user input        |
| > Series Instance UID  | 0020,000E | UI |  | Generated by device    |
| > Series Description   | 0008,103E | LO |  | Zero length            |
| > Retrieve AE Title  | 0008,0054 | AE |  | Zero length            |
| > Referenced Image Sequence                                  | 0008,1140 | SQ |  | From Modality          |
| >> Referenced SOP Class UID                                  | 0008,1150 | UI |  | From Modality          |
| >> Referenced SOP Instance<br>UID                            | 0008,1155 | UI |  | From Modality          |
| > Referenced Non-Image<br>Composite SOP Instance<br>Sequence | 0040,0220 | SQ |  | From Modality          |
| >> Referenced SOP Class UID                                  | 0008,1150 | UI |  | From Modality          |
| >> Referenced SOP Instance<br>UID                            | 0008,1155 | UI |  | From Modality          |

#### 4.2.2.4 Association Acceptance Policy

The Workflow Application Entity does not accept Associations.

### 4.2.3 Hardcopy Application Entity Specification

#### 4.2.3.1 SOP Classes

V7 provides Standard Conformance to the following SOP Classes:

**Table 4.2-30**  
**SOP CLASSES FOR AE HARDCOPY**

| SOP Classes                           | SOP Class UID          | SCU | SCP |
|---------------------------------------|------------------------|-----|-----|
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9  | Yes | No  |
| Basic Color Print Management Meta     | 1.2.840.10008.5.1.1.18 | Yes | No  |

#### 4.2.3.2 Association Policies

##### 4.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-31**  
**DICOM APPLICATION CONTEXT FOR AE HARDCOPY**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 4.2.3.2.2 Number of Association

V7 can initiate one or more Associations at a time for each destination to which a transfer request is being processed in the active job queue list.

**Table 4.2-32**  
**NUMBER OF ASSOCIATIONS INITIATED FOR AE HARDCOPY**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | Unlimited (number of configured hardcopy devices) |
|---|---|

#### 4.2.3.2.3 Asynchronous Nature

V7 does not support asynchronous communications (multiple outstanding transactions over a single Association)

**Table 4.2-33**  
**ASYNCHRONOUS NATURE AS A SCU FOR AE HARDCOPY**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

#### 4.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-34**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE HARDCOPY**


|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

#### 4.2.3.3 Association Initiation Policy

##### 4.2.3.3.1 Activity – Film Images

##### 4.2.3.3.1.1 Description and Sequencing of Activities

A user composes images onto film sheets and requests them to be sent to a specific hardcopy device. The user can select the desired film format and number of copies. Each print-job is forwarded to the job queue and processed individually.



**Figure 4.2-5  
SEQUENCING OF ACTIVITY - FILM IMAGES**

A typical sequence of DIMSE messages sent over an association between Hardcopy AE and a Printer is illustrated in the Figure above:

Association Initiation Policies for “Send on end exam”, “Send after acquisition” and “Send manually” Mode are equal to the Sending images’ of the Storage Application Entity. (See 4.2.1.3.1.1)

Status of the print-job is reported through the job control interface. One or more job can be active at a time for each separate hardcopy device. If any response from the remote Application contains a status other than Success or Warning, the Association is aborted and the related job is switched to a failed state. It can be restarted any time by user interaction or, if configured, by automated retry.

#### **4.2.3.3.1.2 Proposed Presentation Contexts**

V7 is capable of proposing the Presentation Contexts shown in the Table below:

**Table 4.2-35**

**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY FILM IMAGES**

| Presentation Context Table |                |                           |                     |      |           |
|----------------------------|----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID            | Name List                 | UID List            |      |           |
| Basic Grayscale Print      | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
| Management Meta            | 5.1.1.9        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Basic Color Print          | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
| Management Meta            | 5.1.1.18       | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

**4.2.3.3.1.3 Common SOP Specific Conformance for all Print SOP Classes**

The general behavior of Hardcopy AE during communication failure is summarized in the table below. This behavior is common for all SOP Classes supported by Hardcopy AE.

**Table 4.2-36**

**HARDCOPY COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and the print job is marked as failed. |
| Association aborted by the SCP or network layers | The print job is marked as failed.  |

**4.2.3.3.1.4 SOP Specific Conformance for the Film Session SOP Class**

Hardcopy AE supports the following DIMSE operations for the Film Session SOP Class:

- N-CREATE
- N-DELETE

Details of the supported attributes and status handling behavior are described in the following subsections.

**4.2.3.3.1.4.1 Film Session SOP Class Operations (N-CREATE)**

The attributes supplied in an N-CREATE Request are listed in the Table below:

**Table 4.2-37**

**FILM SESSION SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| Attribute Name   | Tag       | VR | Value   | Presence of Value | Source |
|------------------|-----------|----|---|-------------------|--------|
| Number of Copies | 2000,0010 | IS | 1..99   | ALWAYS            | USER   |
| Print Priority   | 2000,0020 | CS | HIGH, MED or LOW  | ALWAYS            | USER   |
| Medium Type      | 2000,0030 | CS | PAPER, CLEAR FILM, BLUE FILM, MAMMO CLEAR FILM or MAMMO BLUE FILM | ALWAYS            | USER   |
| Film Destination | 2000,0040 | CS | MAGAZINE or PROCESSOR   | ALWAYS            | USER   |

The Behavior of Hardcopy AE when encountering status codes in an N-CREATE response is summarized in the table below:

**Table 4.2-38**

**FILM SESSION SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning              | Error Code             | Behavior   |
|----------------|------------------------------|------------------------|--|
| Success        | Success                      | 0000                   | The SCP has Completed the operation successfully.                              |
| Warning        | Attribute Value Out of Range | 0116H                  | System continues operations.   |
| Warning        | Attribute List Error         | 0107H                  | Same as above  |
| *              | *                            | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

**4.2.3.3.1.4.2 Film Session SOP Class Operations (N-DELETE)**

The behavior of Hardcopy AE when encountering status codes in an N-DELETE response is summarized in the Table below:

**Table 4.2-39**

**PRINTER SOP CLASS N-DELETE RESONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code |  |
|----------------|-----------------|------------|--|
|                |                 |            |  |

|         |         |                        |  |
|---------|---------|------------------------|--|
| Success | Success | 0000                   | The SCP has Completed the operation successfully.                              |
| *       | *       | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.5 SOP Specific Conformance for the Film Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.3.3.1.5.1 Film Box SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the table below:

**Table 4.2-40  
FILM BOX SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| Attribute Name                   | Tag       | VR | Value   | Presence of Value | Source |
|----------------------------------|-----------|----|---|-------------------|--------|
| Image Display Format             | 2010,0010 | ST | Enumerated values used (user configurable):<br>STANDARDX,Y,<br>where X and Y can take values from 1 to 5. . | ALWAYS            | USER   |
| Referenced Film Session Sequence | 2010.0500 | SQ |   | ALWAYS            | AUTO   |
| > Referenced SOP Class UID       | 0008,1150 | UI | 1.2.840.10008.5.1.1.1   | ALWAYS            | AUTO   |
| > Referenced SOP Instance UID    | 0008,1155 | UI | From created Film Session SOP Instance  | ALWAYS            | AUTO   |

|                           |           |    |   |        |      |
|---------------------------|-----------|----|---|--------|------|
| Film Orientation          | 2010,0040 | CS | PORTRAIT or LANDSCAPE   | ALWAYS | USER |
| Film Size ID              | 2010,0050 | CS | 8INX10IN, 8_5INX11IN,<br>10INX12IN, 10INX14IN,<br>11INX14IN, 11INX17IN,<br>14INX14IN, 14INX17IN,<br>24CMX24CM, 24CMX30CM,<br>A4, A3 | ALWAYS | USER |
| Magnification Type        | 2010,0060 | CS | REPLICATE, BILINEAR,<br>CUBIC, NONE   | ALWAYS | USER |
| Max Density               | 2010,0130 | US | 0 ~ 65535   | ANAP   | USER |
| Configuration Information | 2010,0150 | ST | Values are defined in Print<br>Conformance Statement  | ANAP   | USER |
| Smoothing Type            | 2010,0080 | CS | Values are defined in Print<br>Conformance Statement  | ANAP   | USER |
| Border Density            | 2010,0100 | CS | BLACK or WHITE  | ALWAYS | USER |
| Empty Image Density       | 2010,0110 | CS | BLACK or WHITE  | ALWAYS | USER |
| Min Density               | 2010,0120 | US | 0 ~ 65535   | ANAP   | USER |

The behavior of Hardcopy AE when encountering status codes in an N-CREATE responses is summarized in the table below:

**Table 4.2-41  
FILM BOX SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>       | <b>Error Code</b> | <b>Behavior</b>                                   |
|-----------------------|------------------------------|-------------------|---|
| Success               | Success                      | 0000              | The SCP has Completed the operation successfully. |
| Warning               | Attribute Value Out of Range | 0116H             | System continues operations.                      |
| Warning               | Attribute List Error         | 0107H             | Same as above                                     |

|         |   |                        |  |
|---------|---|------------------------|--|
| Warning | Requested Min Density or Max Density outside of printer's operating range | B605H                  | Same as above  |
| *       | *   | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.5.2 Film Box SOP Class Operations (N-ACTION)

An N-ACTION Request is issued to instruct the Print SCP to print the contents of the Film Box.

The behavior of Hardcopy AE when encountering status codes in an N-ACTION responses is summarized in the table below:

**Table 4.2-42**

#### **FILM BOX CLASS N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Error Code</b>      | <b>Behavior</b>  |
|-----------------------|------------------------|------------------------|--|
| Success               | Success                | 0000                   | The SCP has Completed the operation successfully.                              |
| *                     | *                      | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.6 SOP Specific Conformance for the Film Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Image Box SOP Class:

- N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.3.3.1.6.1 Image Box SOP Class Operations (N-SET)

The attributes supplied in an N-SET Request are listed in the Table below:

**Table 4.2-43**

**BASIC GRAYSCALE IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES**

| Attribute Name                 | Tag       | VR | Value                                 | Presence of Value | Source |
|--------------------------------|-----------|----|---------------------------------------|-------------------|--------|
| Image Position                 | 2020,0010 | US | 1 .. N (N = Row * Column of Film Box) | ALWAYS            | AUTO   |
| Basic Grayscale Image Sequence | 2020,0110 | SQ |                                       | ALWAYS            | AUTO   |
| > Samples Per Pixel            | 0028,0002 | US | 1                                     | ALWAYS            | AUTO   |
| > Photometric Interpretation   | 0028,0004 | CS | MONOCHROME2                           | ALWAYS            | AUTO   |
| > Rows                         | 0028,0010 | US | Number of Row Pixels of Image         | ALWAYS            | AUTO   |
| > Columns                      | 0028,0011 | US | Number of Column Pixels of Image      | ALWAYS            | AUTO   |
| > Bits Allocated               | 0028,0100 | US | 8                                     | ALWAYS            | AUTO   |
| > Bits Stored                  | 0028,0101 | US | 8                                     | ALWAYS            | AUTO   |
| > High Bit                     | 0028,0102 | US | 7                                     | ALWAYS            | AUTO   |
| > Pixel Representation         | 0028,0103 | US | 0                                     | ALWAYS            | AUTO   |
| > Pixel Data                   | 7FE0,0010 | OB | Pixels of Image                       | ALWAYS            | AUTO   |

**Table 4.2-44**

**BASIC COLOR IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES**

| Attribute Name               | Tag       | VR | Value                                 | Presence of Value | Source |
|------------------------------|-----------|----|---------------------------------------|-------------------|--------|
| Image Position               | 2020,0010 | US | 1 .. N (N = Row * Column of Film Box) | ALWAYS            | AUTO   |
| > Samples Per Pixel          | 0028,0002 | US | 3                                     | ALWAYS            | AUTO   |
| > Photometric Interpretation | 0028,0004 | CS | RGB                                   | ALWAYS            | AUTO   |
| > Planar Configuration       | 0028,0006 | US | 1                                     | ALWAYS            | AUTO   |
| > Rows                       | 0028,0010 | US | Number of Row Pixels of Image         | ALWAYS            | AUTO   |
| > Columns                    | 0028,0011 | US | Number of Column Pixels of            | ALWAYS            | AUTO   |

|                        |           |    |                 |        |      |
|------------------------|-----------|----|-----------------|--------|------|
|                        |           |    | Image           |        |      |
| > Bits Allocated       | 0028,0100 | US | 8               | ALWAYS | AUTO |
| > Bits Stored          | 0028,0101 | US | 8               | ALWAYS | AUTO |
| > High Bit             | 0028,0102 | US | 7               | ALWAYS | AUTO |
| > Pixel Representation | 0028,0103 | US | 0               | ALWAYS | AUTO |
| > Pixel Data           | 7FE0,0010 | OB | Pixels of Image | ALWAYS | AUTO |

The behavior of Hardcopy AE when encountering status codes in an N-SET response is summarized in the table below:

**Table 4.2-45  
IMAGE BOX SOP CLASS N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Error Code</b>      | <b>Behavior</b>  |
|-----------------------|------------------------|------------------------|--|
| Success               | Success                | 0000                   | The SCP has Completed the operation successfully.                              |
| *                     | *                      | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.4 Association Acceptance Policy

The Hardcopy Application Entity does not accept Associations.

## 4.2.4 Q/R Application Entity Specification

### 4.2.4.1 SOP Classes

V7 provides Standard Conformance to the following SOP Classes:

**Table 4.2-46**  
**SOP CLASSES FOR AE Q/R**

| SOP Classes                        | SOP Class UID               | SCU | SCP |
|------------------------------------|-----------------------------|-----|-----|
| Study Root Information Model- FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | No  |
| Study Root Information Model- MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Yes | No  |

### 4.2.4.2 Association Establishment Policy

#### 4.2.4.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed.

**Table 4.2-47**  
**DICOM APPLICATION CONTEXT FOR AE Q/R**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.4.2.2 Number of Associations

V7 initiates one Association at a time for a Q/R request.

**Table 4.2-48**  
**NUMBER OF ASSOCIATIONS INITIATED FOR AE Q/R**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.4.2.3 Asynchronous Nature

V7 does not support asynchronous communications (multiple outstanding transactions over a single Association)

**Table 4.2-49**  
**ASYNCHRONOUS NATURE AS A SCU FOR AE Q/R**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

#### 4.2.4.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-50**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE Q/R**


|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

#### 4.2.4.3 Association Initiation Policy

##### 4.2.4.3.1 Activity – Query study or series

##### 4.2.4.3.1.1 Description and Sequencing of Activities

The Query attempts to initiate a new association when the user selects Query from the user interface. When the Query is requested, the data from the user interface will be inserted as matching keys into the query form. When the request is initiated, the V7 will build an identifier for the C-FIND request, and it will initiate an association to send the request and will wait for Query responses. The results will be displayed in a study or series list.



**Figure 4.2-6**  
**SEQUENCING OF ACTIVITY - HANDLING QUERY STUDY OR SERIES**

**4.2.4.3.1.2 Proposed Presentation Contexts**

V7 will propose Presentation Contexts as shown in the following table:

**Table 4.2-51**  
**PROPOSED PRESENTATION CONTEXTS**  
**FOR REAL-WORLD ACTIVITY QUERY STUDY OR SERIES**

| Presentation Context Table |                |                           |                     |      |           |
|----------------------------|----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID            | Name List                 | UID List            |      |           |
| Study Root                 | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
| Information Model-FIND     | 5.1.4.1.2.2.1  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None      |
|                            |                | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None      |

**4.2.4.3.1.3 SOP Specific Conformance for Query SOP Classes**

The behavior of the V7 when encountering status codes in the Query C-FIND response is summarized in the table below. If any SCP response status other than “Successful” or “Pending” is received by V7, the message “Query failed” will appear in the user interface.

**Table 4.2-52**  
**QUERY C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning  | Error Code             | Behavior  |
|----------------|--|------------------------|---|
| Successful     | Matching is complete   | 0000                   | The SCP has completed the operation successfully.                           |
| Pending        | Matches are continuing   | FF00                   | The query is still ongoing.   |
| Pending        | Matches are continuing - Warning that one or more Optional Keys were not supported | FF01                   | The query is still ongoing.   |
| *              | *  | Any other status code. | The association is aborted using A-Abort and the Query is marked as failed. |

The behavior of the V7 during communication failure is summarized in the table below:

**Table 4.2-53  
QUERY COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Time Out   | The association is aborted using A-ABORT and the query is marked as failed. |
| Association aborted by the SCP or network layers | The study or series query is marked as failed.                              |

The system actually performs a number of C-FIND requests at multiple levels in the DICOM object hierarchy to get the data it requires to display studies or series. Table 4.2-54 provides a description of the query request identifiers.

**Table 4.2-54  
QUERY REQUEST IDENTIFIER FOR FIND-SCU**

| Attribute Name                 | Tag       | VR | M   | R | Q | D |
|--------------------------------|-----------|----|-----|---|---|---|
| <b>STUDY Level</b>             |           |    |     |   |   |   |
| Query/Retrieve Level           | 0008,0052 | CS | S   | x | x |   |
| Patient's ID                   | 0010,0020 | LO | S   | x | x | x |
| Patient's Name                 | 0010,0010 | PN | S   | x | x | x |
| Patient Birth Date             | 0010,0030 | DA |     | x |   |   |
| Patient Sex                    | 0010,0040 | CS |     | x |   |   |
| Patient Age                    | 0010,1010 | AS |     | x |   |   |
| Patient Comments               | 0010,4000 | LT |     | x |   |   |
| Retrieve AE Title              | 0008,0054 | AE |     | x |   |   |
| Study Description              | 0008,1030 | LO |     | x |   | x |
| Modalities In Study            | 0008,0061 | CS |     | x |   | x |
| Study Date                     | 0008,0020 | DA | S,R | x | x | x |
| Study Time                     | 0008,0030 | TM |     | x |   |   |
| Accession Number               | 0008,0050 | SH |     | x | x | x |
| Study Instance UID             | 0020,000D | UI |     | x |   |   |
| Study ID                       | 0020,0010 | SH |     | x |   |   |
| Referring Physician's Name     | 0008,0090 | PN |     | x |   |   |
| Performing Physician's Name    | 0008,1050 | PN |     | x |   |   |
| Number of Study Related Series | 0020,1206 | IS |     | x |   | x |
| <b>SERIES Level</b>            |           |    |     |   |   |   |
|                                |           | CS | S   | x | x |   |

|                                     |           |    |   |   |   |   |
|-------------------------------------|-----------|----|---|---|---|---|
| Query/Retrieve Level                | 0008,0052 |    |   |   |   |   |
| Series Number                       | 0020,0011 | IS |   | x |   | x |
| Series Description                  | 0008,103E | LO |   | x |   | x |
| Modality                            | 0008,0060 | CS |   | x |   | x |
| Series Date                         | 0008,0021 | DA |   | x |   | x |
| Series Time                         | 0008,0031 | TM |   | x |   |   |
| Manufacturer                        | 0008,0070 | LO |   | x |   |   |
| Operators Name                      | 0008,1070 | PN |   | x |   |   |
| Manufacturer Model Name             | 0008,1090 | LO |   | x |   |   |
| Body Part Examined                  | 0018,0015 | CS |   | x |   |   |
| StudyID                             | 0020,1110 | SH |   | x |   |   |
| Performed Procedure Step Start Date | 0040,0244 | DA |   | x |   |   |
| Series Instance UID                 | 0020,000E | UI |   | x |   |   |
| Study Instance UID                  | 0020,000D | UI | S |   | x |   |
| Number of Series Related Instances  | 0020,1209 | IS |   | x |   | x |

The table above should read as follows:

Attribute Name: Supported attributes that can build an V7 Query Request Identifier.

Tag: DICOM tag for this attribute.

VR: DICOM VR for this attribute.

M: Matching keys for (automatic) Query. An "S" indicates that the V7 can supply an attribute value for Single Value Matching or additional specific tags indicated by "(S)"; an "R" indicates Range Matching.

R: Return keys. An "X" indicates that the V7 will supply this attribute as the Return Key with zero length for Universal Matching.


Q: Interactive Query Key. An "X" indicates that the V7 will supply this attribute as a matching key, if entered in the Setup Dialog.

D: Displayed keys. An "X" indicates that this Query attribute is displayed to the user during a patient registration dialog.

### 4.2.4.3.2 Activity – Retrieve series

#### 4.2.4.3.2.1 Description and Sequencing of Activities

The retrieval function attempts to initiate a new association when the user selects Retrieve in the user interface. A single attempt will be made to retrieve the entity (series) from the selected Q/R AE. If retrieval fails, for whatever reason, no reattempt will be performed.



**Figure 4.2-7**  
**SEQUENCING OF ACTIVITY – HANDLING RETRIEVE SERIES**

#### 4.2.4.3.2.2 Proposed Presentation Contexts

V7 will propose Presentation Contexts as shown in the following table:

**Table 4.2-55**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY RETRIEVE SERIES**

| Presentation Context Table |                |                           |                     |      |           |
|----------------------------|----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID            | Name List                 | UID List            |      |           |
| Study Root                 | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCP  | None      |
| Information                | 5.1.4.1.2.2.2  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None      |
| Model- MOVE                |                | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCP  | None      |

#### 4.2.4.3.2.3 SOP Specific Conformance for Retrieval SOP Classes

The behavior of the V7 when encountering status codes in Retrieve C-MOVE response is summarized in the table below. If any SCP response status other than “Successful” or “Pending” is received by the V7, a “failed” message will appear on the user interface.

**Table 4.2-56  
RETRIEVE C-MOVE RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>                         | <b>Error Code</b>      | <b>Behavior</b>   |
|-----------------------|--|------------------------|---|
| Successful            | Sub-operations complete – no failures detected | 0000                   | All the Composite SOP Instances have been successfully sent to the C-MOVE Destination AE.   |
| Pending               | Sub-operations are still ongoing               | FF00                   | A response with this status code is sent every time a Composite SOP Instance has been successfully sent to the C-MOVE Destination AE. |
| *                     | *  | Any other status code. | The association is aborted using A-Abort and the retrieval is marked as failed  |

The behavior of the V7 during communication failure is summarized in the table below.

**Table 4.2-57  
RETRIEVE COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>   |
|--|---|
| Timeout  | The association is aborted using A-ABORT and the retrieval job is marked as failed. |
| Association aborted by the SCP or network layers | The retrieval is marked as failed.  |

**Table 4.2-58  
RETRIEVE REQUEST IDENTIFIER FOR MOVE-SCU**

| <b>Attribute Name</b> | <b>Tag</b> | <b>VR</b> | <b>M</b> | <b>R</b> | <b>Q</b> | <b>D</b> |
|-----------------------|------------|-----------|----------|----------|----------|----------|
| Query/Retrieve Level  | 0008,0052  | CS        | S        |          | x        |          |
| Study Instance UID    | 0020,000D  | UI        | S        |          | x        |          |
| Series Instance UID   | 0020,000E  | UI        | S        |          | x        |          |

#### 4.2.4.1 Association Acceptance Policy

The Q/R Application Entity does not accept associations.

#### 4.2.5 STORAGE-SCP Application Entity Specification

##### 4.2.5.1 SOP Classes

The V7 provides Standard Conformance to the following SOP Classes. Provide Storage SCP only Q/R service running.

**Table 4.2-59  
SOP CLASSES FOR AE STORAGE-SCP**

| SOP Classes                     | SOP Class UID                 | SCU | SCP |
|---------------------------------|-------------------------------|-----|-----|
| Verification                    | 1.2.840.10008.1.1             | Yes | Yes |
| US Image Storage                | 1.2.840.10008.5.1.4.1.1.6.1   | No  | Yes |
| US Multi-frame Storage          | 1.2.840.10008.5.1.4.1.1.3.1   | No  | Yes |
| CT Image Storage                | 1.2.840.10008.5.1.4.1.1.2     | No  | Yes |
| MR Image Storage                | 1.2.840.10008.5.1.4.1.1.4     | No  | Yes |
| MG present Image Storage        | 1.2.840.10008.5.1.4.1.1.1.2   | No  | Yes |
| MG Process Image Storage        | 1.2.840.10008.5.1.4.1.1.1.2.1 | No  | Yes |
| DX present Image Storage        | 1.2.840.10008.5.1.4.1.1.1.1   | No  | Yes |
| DX process Image Storage        | 1.2.840.10008.5.1.4.1.1.1.1.1 | No  | Yes |
| Standard PET Image Storage      | 1.2.840.10008.5.1.4.1.1.128   | No  | Yes |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7     | No  | Yes |

##### 4.2.5.1 Association Establishment Policy

###### 4.2.5.1.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed.

**Table 4.2-60  
DICOM APPLICATION CONTEXT FOR AE STORAGE-SCP**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.5.1.2 Number of Associations

The STORAGE-SCP AE can support multiple simultaneous associations requested by AEs. Each time the STORAGE-SCP AE receives an association request, a child process will be spawned to process the storage.

**Table 4.2-61**

**NUMBER OF ASSOCIATIONS INITIATED FOR AE STORAGE-SCP**

|   |           |
|---|-----------|
| Maximum number of simultaneous Associations | Unlimited |
|---|-----------|

#### 4.2.5.1.3 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-62**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE STORAGE-SCP**

|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

#### 4.2.5.2 Association Initiation Policy

The STORAGE-SCP Application Entity does not initiate associations.

#### 4.2.5.3 Association Acceptance Policy

##### 4.2.5.3.1 Activity – Receive Images

##### 4.2.5.3.1.1 Description and Sequencing of Activities

The STORAGE-SCP AE accepts associations only if they have valid Presentation Contexts. The STORAGE-SCP AE does not have a limit on the number of associations used to send images to it. Images belonging to more than one series can be sent over a single or multiple associations. Images belonging to a single Series can also be sent via different associations.

#### 4.2.5.3.1.2 Proposed Presentation Contexts

The V7 will propose Presentation Contexts as shown in the following table:

**Table 4.2-63  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY RECEIVE IMAGES**

| Presentation Context Table |                 |                           |                        |      |           |
|----------------------------|-----------------|---------------------------|------------------------|------|-----------|
| Abstract Syntax            |                 | Transfer Syntax           |                        | Role | Ext. Neg. |
| Name                       | UID             | Name List                 | UID List               |      |           |
| US Image Storage           | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.6.1   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| US Multi-frame Storage     | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.3.1   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| CT Image Storage           | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.2     | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| MR Image Storage           | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.4     | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| MG present Image Storage   | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.1.2   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| MG Process Image Storage   | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.1.2.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |
|                            |                 | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP  | None      |
| DX present Image Storage   | 1.2.840.10008.  | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP  | None      |
|                            | 5.1.4.1.1.1.1   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP  | None      |
|                            |                 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP  | None      |

|                                       |                                   |                           |                        |     |      |
|---------------------------------------|-----------------------------------|---------------------------|------------------------|-----|------|
|                                       |                                   | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP | None |
| DX process<br>Image Storage           | 1.2.840.10008.<br>5.1.4.1.1.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP | None |
|                                       |                                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP | None |
|                                       |                                   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP | None |
|                                       |                                   | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP | None |
| Standard PET<br>Image Storage         | 1.2.840.10008.<br>5.1.4.1.1.128   | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP | None |
|                                       |                                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP | None |
|                                       |                                   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP | None |
|                                       |                                   | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP | None |
| Secondary<br>Capture Image<br>Storage | 1.2.840.10008.<br>5.1.4.1.1.7     | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCP | None |
|                                       |                                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1    | SCP | None |
|                                       |                                   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2    | SCP | None |
|                                       |                                   | JPEG Baseline(Process 1)  | 1.2.840.10008.1.2.4.50 | SCP | None |

#### 4.2.5.3.1.3 SOP Specific Conformance for Storage SOP Classes

The behavior response of the V7 when encountering status codes in C-STORE is summarized in the table below. If any SCP response status other than "Successful" is received by the V7 it is marked as failed.

**Table 4.2-64  
C-STORE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning                       | Error Code             | Behavior   |
|----------------|---------------------------------------|------------------------|--|
| Successful     | Successfully stored the SOP instance. | 0000                   | The SCP has successfully stored the SOP Instance. If all SOP Instances succeed, the job is marked as complete. |
| *              | *                                     | Any other status code. | The association is aborted using A-Abort and the request to receive the image is marked as failed.             |

## 4.3 NETWORK INTERFACE

### 4.3.1 Physical Network Interface

V7 supports a single network interface. One of the following physical network interfaces will be available depending on hardware options installed:

**Table 4.3-1**  
**SUPPORTED PHYSICAL NETWORK INTERFACES**

|                   |
|-------------------|
| Ethernet 100baseT |
| Ethernet 10baseT  |

### 4.3.2 Additional Protocols

This product may be configured to get the local configuration via the DHCP. However it does not conform to other System Management Profiles as DNS nor LDAP.

### 4.3.3 IPv4 and IPv6 Support

This product supports IPv4 connections only.

## 4.4 CONFIGURATION

### 4.4.1 AE Title/Presentation Address Mapping

#### 4.4.1.1 Local AE Titles

All local applications use the AE Titles and TCP/IP Ports configured via the Setup/Connectivity/DICOM Menu. All local DICOM services use the same AE Title. The system listens for Verification requests and Commitment reports on the configured Port.

#### **4.4.1.2 Remote AE Title/Presentation Address Mapping**

The AE Title, host names and port numbers of remote applications are configured using the V7 Setup/Connectivity/DICOM Menu.

##### **4.4.1.2.1 Storage**

The Add button on the V7 Setup/Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Image Storage SCPs. Multiple remote Image Storage SCPs can be defined.

The Add button on the V7 Setup/ Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Structured Report Storage SCP. Only a single remote Structured Report Storage SCP can be defined.

The Add button on the V7 Setup/Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Storage Commitment SCP. Only a single remote Storage Commitment SCP can be defined and only one Image Storage SCP can be assigned for Storage Commitment.

##### **4.4.1.2.2 Workflow**

The Add button on the V7 Setup/Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Modality Worklist SCP. Only a single remote Modality Worklist SCP can be defined.

The Add button on the V7 Setup/Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote MPPS SCP. Only a single remote MPPS SCP can be defined.

##### **4.4.1.2.3 Hardcopy**

The Add button on the V7 Setup/Connectivity/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Print SCPs. Multiple remote Print SCPs can be defined.

## 4.4.2 Parameters

A number of parameters related to acquisition and general operation can be configured using the Setup/Connectivity/DICOM Menu. The Table below only shows those configuration parameters relevant to DICOM communications. See the V7 Manual for details on general configuration capabilities.

**Table 4.4-1  
CONFIGURATION PARAMETERS TABLE**


| Parameter                        | Configurable<br>(Yes/No) | Default Value       |
|----------------------------------|--------------------------|---------------------|
| <b>Local System Parameters</b>   |                          |                     |
| AE Title (Local System AE Title) | Yes                      | "MEDISON"           |
| Station Name                     | Yes                      | "Set Station Name"  |
| Port No. (Local Port Number)     | Yes                      | 1005                |
| SR Format                        | Yes                      | General Report      |
| Store SR at End of Exam          | Yes                      | UnChecked           |
| <b>Service Common Parameters</b> |                          |                     |
| Alias                            | Yes                      | Blank               |
| AE Title                         | Yes                      | Blank               |
| Host                             | Yes                      | Blank               |
| Port No.                         | Yes                      | 104                 |
| Retry Interval                   | Yes                      | 30 Sec.             |
| Connect Timeout                  | Yes                      | 30 Sec.             |
| Retry Interval                   | Yes                      | 30 Sec              |
| Maximum Retires                  | Yes                      | 1                   |
| Packet Size                      | Yes                      | 16000               |
| <b>Storage Parameters</b>        |                          |                     |
| Transfer Mode                    | Yes                      | "Send on end exam"  |
| Include 3D Volume                | Yes                      | Unchecked           |
| Include Multi Frame              | Yes                      | Checked             |
| Convert to Grayscale             | Yes                      | "No Conversion"     |
| Single Frame Compression         | Yes                      | JPEG Baseline / Low |
| Multi Frame Compression          | Yes                      | JPEG Baseline / Low |
| Multi Frame Frame Rate           | Yes                      | Full (30)           |
| Multi Frame Frame Image Size     | Yes                      | Original            |

|  |     |                                     |
|--|-----|-------------------------------------|
| Include Pixel Spacing                      | Yes | Unchecked                           |
| Window Center (VOI LUT)                    | Yes | 128                                 |
| Window Width (VOI LUT)                     | Yes | 256                                 |
| <b>Storage SR Paramater</b>                |     |                                     |
| Trasfer Mode                               | Yes | "Send on end exam"                  |
| <b>Performed Procedure Step Parameters</b> |     |                                     |
| Always complete exams                      | Yes | Checked                             |
| <b>Storage Commitment Parameters</b>       |     |                                     |
| Associated Storage Server                  | Yes | None                                |
| <b>Worklist Modality Parameters</b>        |     |                                     |
| Modality                                   | Yes | "US"                                |
| Exam Description                           | Yes | "Requested Procedure<br>Decription" |
| <b>Print Parameters</b>                    |     |                                     |
| Transfer Mode                              | Yes | "Send on end exam"                  |
| Color                                      | Yes | "Grayscale"                         |
| Medium Type                                | Yes | "BLUE FILM"                         |
| Image Display Format                       | Yes | Row:2, Col:3                        |
| Film Size                                  | Yes | 8 IN X 10 IN                        |
| Orientation                                | Yes | "PORTRAIT"                          |
| Destination                                | Yes | "MAGAZINE"                          |
| Magnification                              | Yes | "NONE"                              |
| Smoothing Type                             | Yes | Blank                               |
| Border Density                             | Yes | "BLACK"                             |
| Empty Density                              | Yes | "WHITE"                             |
| Priority                                   | Yes | "HIGH"                              |
| Min Density                                | Yes | 0                                   |
| Max Density                                | Yes | 0                                   |
| Copies                                     | Yes | 1                                   |
| Configuration Info                         | Yes | Blank                               |

## 5 MEDIA INTERCHANGE

### 5.1 IMPLEMENTATION MODEL

#### 5.1.1 Application Data Flow



**Figure 5.1-1**

#### **APPLICATION DATA FLOW DIAGRAM FOR MEDIA STORAGE**

- The Offline-Media Application Entity exports images and Structured Report to a Media Storage medium. It is associated with the local real-world activity “Export to Media”, “Export to Media” is performed upon user request for selected studies.

#### 5.1.2 Functional Definition of AEs

##### 5.1.2.1 Functional Definition of Offline-Media Application Entity

Activation of the “Export to Media” menu entry will pass the currently selected studies to the Offline-Media Application Entity. The SOP Instances associated with the selection will be collected into one or more export jobs. The contents of each export job will be written to a single media.

#### 5.1.3 Sequencing of Real-World Activities

At least one study must exist and be selected before the Offline-Media Application Entity can be invoked. The operator can insert a new media at any time before or after invocation of the Offline-Media Application Entity. If no media is available the export job can be cancelled immediately.

#### 5.1.4 File Meta Information Options

The implementation written to the File Meta Header in each file is:

**Table 5.1-1**

## DICOM IMPLEMENTATION CLASS AND VERSION FOR MEDIA STORAGE

|                             |                           |
|-----------------------------|---------------------------|
| Implementation Class UID    | 1.2.410.200001.101.11.801 |
| Implementation Version Name | V7                        |

## 5.2 AE SPECIFICATIONS

### 5.2.1 Offline-Media Application Entity Specification

The Offline-Media Application Entity provides standard conformance to the Media Storage Service Class. The Application Profiles and roles are listed below:

**Table 5.2-1**

#### APPLICATION PROFILES, ACTIVITIES AND ROLES FOR OFFLINE-MEDIA

| Application Profiles Supported | Real World Activity | Role     |
|--------------------------------|---------------------|----------|
| STD-US-SC-MF-CDR               | Export To Media     | FSC, FSU |
| STD-US-SC-MF-DVD               | Export To Media     | FSC, FSU |

#### 5.2.1.1 File Meta Information for the Application Entity

The File-Set Identifier included in the File Meta Header is "MED\_FSU".

#### 5.2.1.2 Real-World Activities

##### 5.2.1.2.1 Activity – Export to Media

The Offline-Media Application Entity acts as an FSC and FSU when requested to export SOP Instances from the local database to a media.

If the contents of the current selection do not fit on a single media, a separation into multiple export jobs which can be adapted by the user will be suggested.

The user will be prompted to insert a media for each export job. The contents of the export job will be written together with a corresponding DICOMDIR to a media. Writing in multi-session mode is supported.

### 5.2.1.2.1.1 Media Storage Application Profiles

The Offline-Media Application Entity supports the STD-US-SC-MF-CDR and STD-US-SC-MF-DVD Application Profile.

#### 5.2.1.2.1.1.1 Options

The Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the table below:

**Table 5.2-2**  
**IODS, SOP CLASSES AND TRANSFER SYNTAXES FOR OFFLINE MEDIA**

| <b>Information Object Definition</b>    | <b>SOP Class UID</b>         | <b>Transfer Syntax</b>          | <b>Transfer Syntax UID</b> |
|---|------------------------------|---------------------------------|----------------------------|
| Media Storage Directory Storage         | 1.2.840.10008.1.3.10         | Explicit VR Little Endian       | 1.2.840.10008.1.2.1        |
| US Image Storage                        | 1.2.840.10008.5.1.4.1.1.6.1  | Explicit VR Little Endian       | 1.2.840.10008.1.2.1        |
| US Multiframe Image Storage             | 1.2.840.10008.5.1.4.1.1.3.1  | JPEG Baseline Lossy Compression | 1.2.840.10008.1.2.4.50     |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.3 | Explicit VR Little Endian       | 1.2.840.10008.1.2.1        |

## **6 SUPPORT OF CHARACTER SETS**

All V7 DICOM applications support the

ISO\_IR 100 : Latin Alphabet No. 1

Supplementary set of ISO 8859

ISO 646

## 7 SECURITY

### 7.1 GENERAL

It is assumed that V7 is used within a secured environment. It is assumed that a secured environment includes as minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to V7.
- b. Firewall or router protections to ensure that V7 has only network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriately secure network channels (e.g. such as a Virtual Private Network (VPN)).

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

### 7.2 SECURE TRANSPORT CONNECTION PROFILES

This system supports the basic TLS Secure Transport Connection Profile and AES TLS Secure Transport Connection Profile. It supports the TLS 1.2 protocol and additional cipher suites for enhanced interoperability and security. Cipher Suites satisfied FIPS 140-2.

**Table 7.2-1  
SUPPORTED PROFILE AND CIPHER SUITES**

| Supported TLS Feature      | Mechanism   |
|----------------------------|---|
| Entity Authentication      | RSA based certificates  |
| Exchange of Master Secrets | RSA   |
| Data Integrity             | SHA   |
| Privacy (Cipher Suites)    | TLS_ECDHE_RSA_WITH_AES_256_GCM_SHA384 (0xc030)<br>TLS_ECDHE_RSA_WITH_AES_256_CBC_SHA384 (0xc028)<br>TLS_ECDHE_RSA_WITH_AES_256_CBC_SHA (0xc014)<br>TLS_DHE_RSA_WITH_AES_256_GCM_SHA384 (0x009f) |

|  |   |
|--|---|
|  | TLS_DHE_RSA_WITH_AES_256_CBC_SHA256 (0x006b)<br>TLS_DHE_RSA_WITH_AES_256_CBC_SHA (0x0039)<br>TLS_RSA_WITH_AES_256_GCM_SHA384 (0x009d)<br>TLS_RSA_WITH_AES_256_CBC_SHA256 (0x003d)<br>TLS_RSA_WITH_AES_256_CBC_SHA (0x0035)<br>TLS_ECDHE_RSA_WITH_AES_128_GCM_SHA256 (0xc02f)<br>TLS_ECDHE_RSA_WITH_AES_128_CBC_SHA256 (0xc027)<br>TLS_ECDHE_RSA_WITH_AES_128_CBC_SHA (0xc013)<br>TLS_DHE_RSA_WITH_AES_128_GCM_SHA256 (0x009e)<br>TLS_DHE_RSA_WITH_AES_128_CBC_SHA256 (0x006f)<br>TLS_DHE_RSA_WITH_AES_128_CBC_SHA (0x0033)<br>TLS_RSA_WITH_AES_128_GCM_SHA256 (0x009c)<br>TLS_RSA_WITH_AES_128_CBC_SHA256 (0x003c)<br>TLS_RSA_WITH_AES_128_CBC_SHA (0x002f)<br>TLS_ECDHE_RSA_WITH_3DES_EDE_CBC_SHA (0xc012)<br>TLS_DHE_RSA_WITH_3DES_EDE_CBC_SHA (0x0016)<br>TLS_RSA_WITH_3DES_EDE_CBC_SHA (0x000a) |
|--|---|

Certificates can be imported to the Trusted Certificate Authority for verifying incoming certificates.

The certificate for this system can be imported via media or created within the device. In case of importing, a private key is required along with the certificate. If the private key has a passphrase, then it must be entered as well.

In case of creating the device, self-signed certificate is generated and can be downloaded to media.

The following types of certificates are supported for import.

Base 64 Text: PEM format (.pem, .crt)

ASN.1 Binary: BER, DER, CER formats (.der, .cer, .crt)

## 8 ANNEXES

### 8.1 IOD CONTENTS

#### 8.1.1 Created SOP Instances

Table 8.1-1 specifies the attributes of an Ultrasound Image transmitted by the V7 storage applications.

8.1-3 specifies the attributes of a Comprehensive Structured Reports transmitted by the V7 storage applications.

The following tables use a number of abbreviations. The abbreviations used in the “Presence of...” column are:

- VNAP Value Not Always Present (attribute sends zero length if no value is present)
- ANAP Attribute Not Always Present
- ALWAYS Always Present
- EMPTY Attribute is sent without a value

The abbreviations used in the “Source” column:

- MWL the attribute value source Modality Worklist
- USER the attribute value source is from User input
- AUTO the attribute value is generated automatically
- MPPS the attribute value is the same as the Modality Performed Procedure Step service
- CONFIG the attribute value source is a configurable parameter

NOTE: All dates and times are encoded in the local configured calendar and time. Date, Time and Time zones are configured using the Setup Menu.

#### 8.1.1.1 US or US Multiframe Image IOD

**Table 8.1-1**

**IOD OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| IE      | Module        | Reference   | Presence of Module |
|---------|---------------|-------------|--------------------|
| Patient | Patient       | Table 8.1-4 | ALWAYS             |
| Study   | General Study | Table 8.1-5 | ALWAYS             |

|           |                       |              |                       |
|-----------|-----------------------|--------------|-----------------------|
|           | Patient Study         | Table 8.1-6  | ALWAYS                |
| Series    | General Series        | Table 8.1-7  | ALWAYS                |
| Equipment | General Equipment     | Table 8.1-8  | ALWAYS                |
| Image     | General Image         | Table 8.1-9  | ALWAYS                |
|           | Image Pixel           | Table 8.1-10 | ALWAYS                |
|           | Cine                  | Table 8.1-11 | Only if US Multiframe |
|           | Multi-Frame           | Table 8.1-12 | Only if US Multiframe |
|           | US Region Calibration | Table 8.1-13 | ANAP                  |
|           | US Image              | Table 8.1-14 | ALWAYS                |
|           | VOI LUT               | Table 8.1-15 | ALWAYS                |
|           | SOP Common            | Table 8.1-17 | ALWAYS                |

#### 8.1.1.1.1 Additional Module

**Table 8.1-2**  
**ADDITIONAL MODULES**

| Module      | Reference    | Presence of Module |
|-------------|--------------|--------------------|
| Image Plane | Table 8.1-16 | ANAP               |

#### 8.1.1.2 Comprehensive Structured Report IOD

**Table 8.1-3**  
**IOD OF CREATED COMPREHENSIVE STRUCTURED REPORT SOP INSTANCES**

| IE        | Module              | Reference    | Presence of Module |
|-----------|---------------------|--------------|--------------------|
| Patient   | Patient             | Table 8.1-4  | ALWAYS             |
| Study     | General Study       | Table 8.1-5  | ALWAYS             |
|           | Patient Study       | Table 8.1-6  | ALWAYS             |
| Series    | SR Document Series  | Table 8.1-18 | ALWAYS             |
| Equipment | General Equipment   | Table 8.1-8  | ALWAYS             |
| Document  | SR Document General | Table 8.1-19 | ALWAYS             |
|           | SR Document Content | Table 8.1-20 | ALWAYS             |
|           | SOP Common          | Table 8.1-21 | ALWAYS             |

### 8.1.1.3 Common Modules

**Table 8.1-4  
PATIENT MODULE OF CREATED SOP INSTANCES**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source        |
|----------------------------|-----------|----|---|-------------------|---------------|
| Patient's Name             | 0010,0010 | PN | From MWL or User Input. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain first 3 components (Last^First^Middle). Maximum 64 characters. | VNAP              | MWL/USER      |
| Patient ID                 | 0010,0020 | LO | From MWL, user input or generated by device. Maximum 64 characters.   | ALWAYS            | MWL/USER/AUTO |
| Patient's Birth Date       | 0010,0030 | DA | From MWL or user input  | VNAP              | MWL/USER      |
| Patient's Sex              | 0010,0040 | CS | From MWL or user input  | VNAP              | MWL/USER      |
| Other Patient IDs          | 0010,1000 | LO | From MWL or user input  | VNAP              | MWL/USER      |
| Other Patient IDs Sequence | 0010,1002 | SQ | From MWL or user input  | VNAP              | MWL/USER      |
| > Patient ID               | 0010,0020 | LO | From MWL or user input  | VNAP              | MWL/USER      |
| > Type of Patient ID       | 0010,0022 | CS | From MWL or user input  | VNAP              | MWL/USER      |

**Table 8.1-5  
GENERAL STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------|-----|----|-------|-------------------|--------|
|----------------|-----|----|-------|-------------------|--------|

|                               |           |    |  |        |          |
|-------------------------------|-----------|----|--|--------|----------|
| Study Instance UID            | 0020,000D | UI | From MWL or generated by device  | ALWAYS | MWL/AUTO |
| Study Date                    | 0008,0020 | DA | <yyyymmdd>   | ALWAYS | AUTO     |
| Study Time                    | 0008,0030 | TM | <hhmmss>   | ALWAYS | AUTO     |
| Referring Physician's Name    | 0008,0090 | PN | From MWL or user input   | VNAP   | MWL/USER |
| Study ID                      | 0020,0010 | SH | From Requested Procedure UID or System generate : Study Date + Study Time<br><yyyymmddhhmmss>  | ALWAYS | AUTO     |
| Accession Number              | 0008,0050 | SH | From MWL or user input   | VNAP   | MWL/USER |
| Study Description             | 0008,1030 | LO | From MWL (Scheduled procedure step description, Requested procedure description) or user input | ANAP   | MWL/USER |
| Referenced Study Sequence     | 0008,1110 | SQ | From MWL   | ANAP   | MWL      |
| > Referenced SOP Class UID    | 0008,1150 | UI | From MWL   | ANAP   | MWL      |
| > Referenced SOP Instance UID | 0008,1155 | UI | From MWL   | ANAP   | MWL      |
| Procedure Code Sequence       | 0008,1032 | SQ | From MWL   | ANAP   | MWL      |

**Table 8.1-6**

**PATIENT STUDY MODULE OF CREATED SOP INSTANCES**

| <b>Attribute Name</b> | <b>Tag</b> | <b>VR</b> | <b>Value</b>           | <b>Presence of Value</b> | <b>Source</b> |
|-----------------------|------------|-----------|------------------------|--------------------------|---------------|
| Patient's Size        | 0010,1020  | DS        | From MWL or user input | ANAP                     | MWL/USER      |

|                  |           |    |                        |      |          |
|------------------|-----------|----|------------------------|------|----------|
| Patient's Weight | 0010,1030 | DS | From MWL or user input | ANAP | MWL/USER |
|------------------|-----------|----|------------------------|------|----------|

**Table 8.1-7  
GENERAL SERIES MODULE OF CREATED SOP INSTANCES**

| Attribute Name                               | Tag       | VR | Value   | Presence of Value | Source   |
|--|-----------|----|---|-------------------|----------|
| Modality                                     | 0008,0060 | CS | US  | ALWAYS            | AUTO     |
| Series Instance UID                          | 0020,000E | UI | Generated by device   | ALWAYS            | AUTO     |
| Series Number                                | 0020,0011 | IS | Generated by device, increments from "1" in each study          | ALWAYS            | AUTO     |
| Series Date                                  | 0008,0021 | DA | <yyyymmdd>  | ALWAYS            | AUTO     |
| Series Time                                  | 0008,0031 | TM | <hhmmss>  | ALWAYS            | AUTO     |
| Performing Physician's Name                  | 0008,1050 | PN | From MWL or user input  | ANAP              | MWL/USER |
| Operators' Name                              | 0008,1070 | PN | From user input   | ANAP              | USER     |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ALWAYS            | MPPS     |
| > Referenced SOP Class UID                   | 0008,1150 | UI | MPPS SOP Class UID "1.2.840.10008.3.1.2.3.3"                    | ALWAYS            | MPPS     |
| > Referenced SOP Instance UID                | 0008,1155 | UI | MPPS SOP Instance UID   | ALWAYS            | MPPS     |
| Body Part Examined                           | 0018,0015 | CS | From user input   | ANAP              | USER     |
| Request Attributes Sequence                  | 0040,0275 | SQ | Zero or 1 item will be present                                  | ANAP              | AUTO     |
| > Requested Procedure ID                     | 0040,1001 | SH | From MWL  | ANAP              | MWL      |

|  |           |    |                           |        |              |
|--|-----------|----|---------------------------|--------|--------------|
| > Scheduled<br>Procedure Step<br>ID          | 0040,0009 | SH | From MWL                  | ANAP   | MWL          |
| > Scheduled<br>Procedure Step<br>Description | 0040,0007 | LO | From MWL                  | ANAP   | MWL          |
| > Scheduled<br>Protocol Code<br>Sequence     | 0040.0008 | SQ | From MWL                  | ANAP   | MWL          |
| Performed<br>Procedure Step<br>ID            | 0040,0253 | SH | Same as MPPS              | ALWAYS | MPPS         |
| Performed<br>Procedure Step<br>Start Date    | 0040,0244 | DA | Same as Study Date        | ALWAYS | AUTO         |
| Performed<br>Procedure Step<br>Start Time    | 0040,0245 | TM | Same as Study Time        | ALWAYS | AUTO         |
| Performed<br>Procedure Step<br>Description   | 0040,0254 | LO | Same as Study Description | ANAP   | MWL/US<br>ER |

**Table 8.1-8**

**GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES**

| <b>Attribute Name</b>        | <b>Tag</b> | <b>VR</b> | <b>Value</b>                      | <b>Presence<br/>of Value</b> | <b>Source</b> |
|------------------------------|------------|-----------|-----------------------------------|------------------------------|---------------|
| Manufacturer                 | 0008,0070  | LO        | SAMSUNG MEDISON CO., LTD,<br>Ltd. | ALWAYS                       | AUTO          |
| Institution Name             | 0008,0080  | LO        | From user input                   | ANAP                         | CONFIG        |
| Station Name                 | 0008,1010  | SH        | From user input                   | ANAP                         | CONFIG        |
| Manufacturer's<br>Model Name | 0008,1090  | LO        | V7                                | ALWAYS                       | AUTO          |
| Device Serial<br>Number      | 0018,1000  | LO        | Generated by device               | ALWAYS                       | AUTO          |
| Software                     | 0018,1020  | LO        | Generated by device               | ALWAYS                       | AUTO          |

|          |  |  |  |  |  |
|----------|--|--|--|--|--|
| Versions |  |  |  |  |  |
|----------|--|--|--|--|--|

#### 8.1.1.4 US or US Multiframe Image Module

**Table 8.1-9**

**GENERAL IMAGE MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name                 | Tag       | VR | Value   | Presence of Value | Source |
|--------------------------------|-----------|----|---|-------------------|--------|
| Instance Number                | 0020,0013 | IS | Generated by device, increments from "1" in each series                                   | ALWAYS            | AUTO   |
| Patient Orientation            | 0020,0020 | CS | NULL  |                   |        |
| Content Date                   | 0008,0023 | DA | <yyyymmdd>  | ALWAYS            | AUTO   |
| Content Time                   | 0008,0033 | TM | <hhmmss>  | ALWAYS            | AUTO   |
| Image Type                     | 0008,0008 | CS | "ORIGINAL" and "PRIMARY"  | ALWAYS            | AUTO   |
| Acquisition Date               | 0008,0022 | DA | <yyyymmdd>  | ALWAYS            | AUTO   |
| Acquisition Time               | 0008,0032 | TM | <hhmmss>  | ALWAYS            | AUTO   |
| Acquisition DateTime           | 0008,002A | DT | <yyyymmddhhmmss>  | ALWAYS            | AUTO   |
| Ultrasound Color Data Present  | 0028,0014 | US | Color data not present = "00"<br>Color data is present = "01"                             | ALWAYS            | AUTO   |
| Lossy Image Compression        | 0028,2110 | CS | US = "00" (uncompressed) or<br>"01" (lossy compressed)<br>US-MF = "01" (lossy compressed) | ALWAYS            | AUTO   |
| Lossy Image Compression Ratio  | 0028,2112 | DS | Used if (0028, 2110) = "01",<br>Calculated by device                                      | ANAP              | AUTO   |
| Lossy Image Compression Method | 0028,2114 | CS | "ISO_10918_1", used if<br>(0028,2110) = "01"  | ANAP              | AUTO   |

**Table 8.1-10**

**IMAGE PIXEL MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name             | Tag       | VR             | Value   | Presence of Value | Source |
|----------------------------|-----------|----------------|---|-------------------|--------|
| Samples per Pixel          | 0028,0002 | US             | "3" for RGB or YBR_FULL_422<br>"1" for MONOCHROME2                      | ALWAYS            | AUTO   |
| Photometric Interpretation | 0028,0004 | CS             | Uncompressed = "RGB" or<br>"MONOCHROME2"<br>Compressed = "YBR_FULL_422" | ALWAYS            | AUTO   |
| Rows                       | 0028,0010 | US             | US = "924", US-MF = CONFIG<br>(Default 924)                             | ALWAYS            | AUTO   |
| Columns                    | 0028,0011 | US             | US = "1232", US-MF = CONFIG<br>(Default 1232)                           | ALWAYS            | AUTO   |
| Bits Allocated             | 0028,0100 | US             | "8"   | ALWAYS            | AUTO   |
| Bits Stored                | 0028,0101 | US             | "8"   | ALWAYS            | AUTO   |
| High Bit                   | 0028,0102 | US             | "7"   | ALWAYS            | AUTO   |
| Pixel Representation       | 0028,0103 | US             | "0"   | ALWAYS            | AUTO   |
| Pixel Data                 | 7FE0,0010 | OW<br>or<br>OB | Generated by device   | ALWAYS            | AUTO   |
| Planar Configuration       | 0028,0006 | US             | "0"   | ALWAYS            | AUTO   |
| Private Creator            | 7FE1,0010 | LO             | "MEDISON_US"  | ANAP              | AUTO   |
| 3D Volume                  | 7FE1,1002 | OB             | 3D Volume Data  | ANAP              | AUTO   |

**Table 8.1-11**

**CINE MODULE OF CREATED US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value             | Presence of Value | Source |
|----------------|-----------|----|-------------------|-------------------|--------|
| Frame Time     | 0018,1063 | DS | Milliseconds      | ANAP              | AUTO   |
| Cine Rate      | 0018,0040 | IS | Frames per second | ANAP              | AUTO   |

**Table 8.1-12**

**MULTI-FRAME MODULE OF CREATED US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------|-----|----|-------|-------------------|--------|
|----------------|-----|----|-------|-------------------|--------|

|                         |           |    |                          |      |      |
|-------------------------|-----------|----|--------------------------|------|------|
| Number of Frames        | 0028,0008 | IS | Numbers of Frames        | ANAP | AUTO |
| Frame Increment Pointer | 0028,0009 | AT | "1577059" : (0018, 1063) | ANAP | AUTO |

**Table 8.1-13**

**US REGION CALIBRATION MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name                 | Tag       | VR | Value  | Presence of Value | Source |
|--------------------------------|-----------|----|--|-------------------|--------|
| Sequence of Ultrasound Regions | 0018,6011 | SQ | Generated by device. A sequence is present for each region in the system display.            | ANAP              | AUTO   |
| > Region Location Min x0       | 0018,6018 | UL | Left position of region  | ALWAYS            | AUTO   |
| > Region Location Min y0       | 0018,601A | UL | Top position of region   | ALWAYS            | AUTO   |
| > Region Location Max x1       | 0018,601C | UL | Right position of region   | ALWAYS            | AUTO   |
| > Region Location Max y1       | 0018,601E | UL | Bottom position of region  | ALWAYS            | AUTO   |
| > Physical Units X Direction   | 0018,6024 | US | 2D Image : 0003H = cm<br>M-Mode : 0004H = seconds<br>Doppler : 0004H = seconds               | ALWAYS            | AUTO   |
| > Physical Units Y Direction   | 0018,6026 | US | 2D Image : 0003H = cm<br>M-Mode : 0003H = cm<br>Doppler : 0005H = hertz or<br>0007H = cm/sec | ALWAYS            | AUTO   |
| > Physical Delta X             | 0018,602C | FD | The physical value per pixel increment   | ALWAYS            | AUTO   |
| > Physical Delta Y             | 0018,602E | FD | The physical value per pixel increment   | ALWAYS            | AUTO   |

|                         |           |    |  |        |      |
|-------------------------|-----------|----|--|--------|------|
| > Region Spatial Format | 0018,6012 | US | 2D Tissue : 0001H<br>M-Mode Tissue or flow : 0002H<br>Spectral (CW or PW Doppler) :<br>0003H       | ALWAYS | AUTO |
| > Region Data Type      | 0018,6014 | US | Tissue : 0001H<br>Color Flow : 0002H<br>PW Spectral Doppler : 0003H<br>CW Spectral Doppler : 0004H | ALWAYS | AUTO |
| > Region Flags          | 0018,6016 | UL | See DICOM PS 3.3 C.8.5.5.1.3   | ALWAYS | AUTO |

**Table 8.1-14**

**US IMAGE MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source |
|----------------------------|-----------|----|---|-------------------|--------|
| Samples Per Pixel          | 0028,0002 | US | "3" for RGB or YBR_FULL_422<br>"1" for MONOCHROME2  | ALWAYS            | AUTO   |
| Photometric Interpretation | 0028,0004 | CS | Uncompressed = "RGB" or<br>"MONOCHROME2"<br>Compressed = "YBR_FULL_422"                   | ALWAYS            | AUTO   |
| Bits Allocated             | 0028,0100 | US | "8"   | ALWAYS            | AUTO   |
| Bits Stored                | 0028,0101 | US | "8"   | ALWAYS            | AUTO   |
| High Bit                   | 0028,0102 | US | "7"   | ALWAYS            | AUTO   |
| Planar Configuration       | 0028,0006 | US | "0"   | ALWAYS            | AUTO   |
| Pixel Representation       | 0028,0103 | US | "0"   | ALWAYS            | AUTO   |
| Image Type                 | 0008,0008 | CS | "ORIGINAL" and "PRIMARY"  | ALWAYS            | AUTO   |
| Lossy Image Compression    | 0028,2110 | CS | US = "00" (uncompressed) or<br>"01" (lossy compressed)<br>US-MF = "01" (lossy compressed) | ALWAYS            | AUTO   |

**Table 8.1-15**

**VOI LUT MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value           | Presence of Value | Source |
|----------------|-----------|----|-----------------|-------------------|--------|
| Window Center  | 0028,1050 | DS | default : "128" | ALWAYS            | CONFIG |
| Window Width   | 0028,1051 | DS | default : "256" | ALWAYS            | CONFIG |

**Table 8.1-16**

**IMAGE PLANE MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value  | Presence of Value | Source |
|----------------|-----------|----|--|-------------------|--------|
| Pixel Spacing  | 0028,0030 | DS | <p>In case that All following conditions are satisfied, This attribute is inserted.</p> <ol style="list-style-type: none"> <li>1. User shall select the option activating Pixel Spacing at the DICOM Setup.</li> <li>2. Image shall have regions consisting of only tissue and color</li> <li>3. For all regions, Units for X and Y direction shall be "cm"</li> <li>4. For all regions, Delta X of US Region calibration module shall have the same value.</li> <li>5. For all regions, Delta Y of US Region calibration module shall have the same value.</li> </ol> | ANAP              | AUTO   |

**Table 8.1-17**

**SOP COMMON MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value   | Presence of Value | Source |
|----------------|-----------|----|---|-------------------|--------|
| SOP Class UID  | 0008,0016 | UI | US =<br>"1.2.840.10008.5.1.4.1.1.6.1"<br>US-MF =<br>"1.2.840.10008.5.1.4.1.1.3.1" | ALWAYS            | AUTO   |

|                        |           |    |  |        |      |
|------------------------|-----------|----|--|--------|------|
| SOP Instance UID       | 0008,0018 | UI | Generated by device                      | ALWAYS | AUTO |
| Specific Character Set | 0008,0005 | CS | Ref. Section 6 SUPPORT OF CHARACTER SETS | ALWAYS | AUTO |

### 8.1.1.5 Comprehensive Structured Report Modules

**Table 8.1-18**

#### SR DOCUMENT SERIES MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES

| Attribute Name                               | Tag       | VR | Value   | Presence of Value | Source |
|--|-----------|----|---|-------------------|--------|
| Modality                                     | 0008,0060 | CS | SR  | ALWAYS            | AUTO   |
| Series Instance UID                          | 0020,000E | UI | Generated by device   | ALWAYS            | AUTO   |
| Series Number                                | 0020,0011 | IS | "2"   | ALWAYS            | AUTO   |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ALWAYS            | MPPS   |
| > Referenced SOP Class UID                   | 0008,1150 | UI | MPPS SOP Class UID "1.2.840.10008.3.1.2.3.3"                    | ALWAYS            | MPPS   |
| > Referenced SOP Instance UID                | 0008,1155 | UI | MPPS SOP Instance UID   | ALWAYS            | MPPS   |

**Table 8.1-19**

#### SR DOCUMENT GENERAL MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES

| Attribute Name    | Tag       | VR | Value   | Presence of Value | Source |
|-------------------|-----------|----|---|-------------------|--------|
| Instance Number   | 0020,0013 | IS | Generated by device, increments from "1" in each series | ALWAYS            | AUTO   |
| Completion Flag   | 0040,A491 | CS | "PARTIAL"   | ALWAYS            | AUTO   |
| Verification Flag | 0040,A493 | CS | "UNVERIFIED"  | ALWAYS            | AUTO   |
| Content Date      | 0008,0023 | DA | <yyyymmdd>  | ALWAYS            | AUTO   |

|   |           |    |                                 |        |          |
|---|-----------|----|---------------------------------|--------|----------|
| Content Time                                  | 0008,0033 | TM | <hhmmss>                        | ALWAYS | AUTO     |
| Referenced Request Sequence                   | 0040,A370 | SQ | 1 item will be present          | ANAP   | AUTO     |
| > Study Instance UID                          | 0020,000D | UI | From MWL or generated by device | ANAP   | MWL/AUTO |
| > Referenced Study Sequence                   | 0008,1110 | SQ | From MWL                        | ANAP   | MWL      |
| >> Referenced SOP Class UID                   | 0008,1150 | UI | From MWL                        | ANAP   | MWL      |
| >> Referenced SOP Instance UID                | 0008,1155 | UI | From MWL                        | ANAP   | MWL      |
| > Accession Number                            | 0008,0050 | SH | From MWL or user input          | VNAP   | MWL/USER |
| > Placer Order Number/Imaging Service Request | 0040,2016 | LO | NULL                            | VNAP   | AUTO     |
| > Filler Order Number/Imaging Service Request | 0040,2017 | LO | NULL                            | VNAP   | AUTO     |
| > Requested Procedure ID                      | 0040,1001 | SH | From MWL                        | VNAP   | MWL      |
| > Requested Procedure Description             | 0032,1060 | LO | From MWL                        | VNAP   | MWL      |
| > Requested Procedure Code Sequence           | 0032,1064 | SQ | From MWL                        | VNAP   | MWL      |
| Performed Procedure Code Sequence             | 0040,A372 | SQ | NULL                            | VNAP   | AUTO     |

**Table 8.1-20**

**SR DOCUMENT CONTENT MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| Attribute Name                        | Tag       | VR | Value   | Presence of Value | Source |
|---------------------------------------|-----------|----|---|-------------------|--------|
| Value Type                            | 0040,A040 | CS | "CONTAINER"   | ALWAYS            | AUTO   |
| Concept Name Code Sequence            | 0040,A043 | SQ | 1 item will be present  | ALWAYS            | AUTO   |
| > Include 'Code Sequence Macro'       |           |    | "EV(125000, DCM, "OB-GYN Ultrasound Procedure Report") for OB-GYN | ALWAYS            | AUTO   |
| Include 'Container Macro'             |           |    |   | ALWAYS            | AUTO   |
| Content Sequence                      | 0040,A730 | SQ | One or more items may be included in this sequence                | ALWAYS            | AUTO   |
| > Relationship Type                   | 0040,A010 | CS | Ref. Section 9 STRUCTURED REPORT TEMPLATES                        | ALWAYS            | AUTO   |
| > Include Document Relationship Macro |           |    | Ref. Section 9 STRUCTURED REPORT TEMPLATES                        | ALWAYS            | AUTO   |
| > Include Document Content Macro      |           |    | Ref. Section 9 STRUCTURED REPORT TEMPLATES                        | ALWAYS            | AUTO   |

**Table 8.1-21**

**SOP COMMON MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| Attribute Name         | Tag       | VR | Value                                    | Presence of Value | Source |
|------------------------|-----------|----|--|-------------------|--------|
| SOP Class UID          | 0008,0016 | UI | "1.2.840.10008.5.1.4.1.1.88.33"          | ALWAYS            | AUTO   |
| SOP Instance UID       | 0008,0018 | UI | Generated by device                      | ALWAYS            | AUTO   |
| Specific Character Set | 0008,0005 | CS | Ref. Section 6 SUPPORT OF CHARACTER SETS | ALWAYS            | AUTO   |

**8.1.2 Used Fields in received IOD by application**

The V7 storage application does not receive SOP Instances. The usage of attributes received via Modality Worklist is described in section 4.2.2.3.1.3.

### 8.1.3 Attribute mapping

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in the Table below. The format and conversions used in Table are the same as the corresponding table in IHE Technical Framework.

**Table 8.1-22  
ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST, IMAGE AND MPPS**

| <b>Modality Worklist</b>              | <b>Image IOD</b>                       | <b>MPPS IOD</b>                        |
|---------------------------------------|--|--|
| Patient's Name                        | Patient's Name                         | Patient's Name                         |
| Patient ID                            | Patient ID                             | Patient ID                             |
| Patient's Birth Date                  | Patient's Birth Date                   | Patient's Birth Date                   |
| Patient's Sex                         | Patient's Sex                          | Patient's Sex                          |
| Patient's Size                        | Patient's Size                         | _____                                  |
| Patient's Weight                      | Patient's Weight                       | _____                                  |
| Referring Physician's Name            | Referring Physician's Name             | _____                                  |
| Scheduled Performing Physician's Name | Performing Physician's Name            | Performing Physician's Name            |
|                                       | _____                                  | Scheduled Step Attributes Sequence     |
| Study Instance UID                    | Study Instance UID                     | > Study Instance UID                   |
| Referenced Study Sequence             | Referenced Study Sequence              | > Referenced Study Sequence            |
| Accession Number                      | Accession Number                       | > Accession Number                     |
|                                       | Request Attributes Sequence            | _____                                  |
| Requested Procedure ID                | > Requested Procedure ID               | > Requested Procedure ID               |
| Requested Procedure Description       | _____                                  | > Requested Procedure Description      |
| Scheduled Procedure Step ID           | > Scheduled Procedure Step ID          | > Scheduled Procedure Step ID          |
| Scheduled Procedure Step Description  | > Scheduled Procedure Step Description | > Scheduled Procedure Step Description |
| Scheduled Protocol Code Sequence      | > Scheduled Protocol Code Sequence     | > Scheduled Protocol Code Sequence     |
| Requested Procedure ID                | Study ID                               | Study ID                               |
| _____                                 | Performed Procedure Step ID            | Performed Procedure Step ID            |

|                                   |  |                                      |
|-----------------------------------|--|--------------------------------------|
| _____                             | Performed Procedure Step Start Date          | Performed Procedure Step Start Date  |
| _____                             | Performed Procedure Step Start Time          | Performed Procedure Step Start Time  |
| _____                             | Performed Procedure Step Description         | Performed Procedure Step Description |
| _____                             | _____  | Performed Series Sequence            |
| Requested Procedure Code Sequence | Procedure Code Sequence                      | Procedure Code Sequence              |
| _____                             | Referenced Performed Procedure Step Sequence | _____                                |
| _____                             | > Referenced SOP Class UID                   | SOP Class UID                        |
| _____                             | > Referenced SOP Instance UID                | SOP Instance UID                     |

### 8.1.4 Coerced/Modified Fields

The Modality Worklist AE will truncate attribute values received in the response to a Modality Worklist Query if the value length is longer than the maximum length permitted by the attribute's VR.

## 8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

The Private Attributes added to create SOP Instances are listed in the Table below. V7 reserves blocks of private attributes in groups 7FE1. Further details on usage of these private attributes are contained in Section 8.1

**Table 8.2-1  
DATA DICTIONARY OF PRIVATE ATTRIBUTES**

| Tag          | Attribute Name  | VR | VM |
|--------------|-----------------|----|----|
| (7FE1, 0010) | Private Creator | LO | 1  |
| (7FE1, 1002) | 3D Volume       | OB | 1  |

## 8.3 CODED TERMINOLOGY AND TEMPLATES

The Workflow AE is capable of supporting arbitrary coding schemes for Procedure and Protocol Codes. The contents of Requested Procedure Code Sequence (0032, 1064) and Scheduled Protocol Code Sequence (0040, 0008) supplied in Worklist Items will be mapped to Image IOD and MPPS attributes as described in Section 8.1.3

## **8.4 STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES**

No Specialized or Private SOP Classes are supported.

### **8.4.1 US OR US MULTIFRAME IMAGE STORAGE SOP CLASS**

The US or US Multiframe Image Storage SOP Classes are extended to create a Standard Extended SOP Class by addition of standard and private attributes to the created SOP Instances as documented in section 8.1

3D Volume Data is transferred to the configured Storage Server, if “Send 3D Volume” option is enabled in the Setup Dialog.

## **8.5 PRIVATE TRANSFER SYNTAXES**

No Private Transfer Syntaxes are supported.

## 9 STRUCTURED REPORT TEMPLATES

This Section uses the following forms for describing Structured Report Templates used in V7.

|   | Rel with Parent | VT | Concept Name | Presence of Value | Comments |
|---|-----------------|----|--------------|-------------------|----------|
| 1 |                 |    |              |                   |          |
| 2 |                 |    |              |                   |          |

|     | NL | REL | VT | Concept Name | Unit / CODE Value | Ref TID | Ref CID | Comments or Label |
|-----|----|-----|----|--------------|-------------------|---------|---------|-------------------|
| A-1 |    |     |    |              |                   |         |         |                   |
| A-2 |    |     |    |              |                   |         |         |                   |

|                   |   |
|-------------------|---|
| Rel with Parent   | Relationship  |
| VT                | Value Type  |
| Concept Name      | Any constraints on Concept Name are specified in this field as defined or enumerated coded entries, or as baseline or defined context groups.   |
| Presence of Value | Ref. Section 8.1.1  |
| Comments          | Description about Reference section or used values.   |
| Label             | Name which is indicated in the system   |
| NL                | The nesting level of Content Items is denoted by ">" symbols  |
| REL               | Relationship  |
| Unit/Code, Value  | Applied unit, enumerated coded entries, or the reference of Context Group.  |
| Ref TID           | Referenced Template ID Number   |
| Ref CID           | Referenced Context ID Number. The left side of "/" shows a CID value applied in "Concept Name" column and the right side shows a CID value applied in "Unit/Code, Value" column. (e.g. 228/12012) |

## 9.1 OB-GYN REPORT TEMPLATES

### 9.1.1 OB-GYN Ultrasound Procedure Report (TID 5000)

**Table 9.1-1  
OB-GYN ULTRASOUND PROCEDURE REPORT TEMPLATE**

|    | Rel with Parent | VT        | Concept Name   | Presence of Value | Comments              |
|----|-----------------|-----------|--|-------------------|-----------------------|
| 1  |                 | CONTAINER | EV (125000, DCM, "OB-GYN Ultrasound Procedure Report") | ALWAYS            |                       |
| 2  | HAS CONCEPT MOD | INCLUDE   | DTID (1204) Language of Content Item and Descendants   |                   |                       |
| 3  | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observation Context                        | ALWAYS            | Ref. Section 9.1.1.1  |
| 4  | CONTAINS        | INCLUDE   | DTID (5001) Patient Characteristics                    | ANAP              | Ref. Section 9.1.1.2  |
| 5  | CONTAINS        | CONTAINER | DT (111028, DCM, "Image Library")                      |                   |                       |
| 6  | CONTAINS        | IMAGE     | No Purpose of reference                                |                   |                       |
| 7  | CONTAINS        | INCLUDE   | DTID (5002) OB-GYN Procedure Summary Section           | ANAP              | Ref. Section 9.1.1.3  |
| 8  | CONTAINS        | INCLUDE   | DTID (5004) Fetal Biometry Ratio Section               | ANAP              | Ref. Section 9.1.1.4  |
| 9  | CONTAINS        | INCLUDE   | DTID (5005) Fetal Biometry Section                     | ANAP              | Ref. Section 9.1.1.5  |
| 10 | CONTAINS        | INCLUDE   | DTID (5006) Fetal Long Bones Section                   | ANAP              | Ref. Section 9.1.1.6  |
| 11 | CONTAINS        | INCLUDE   | DTID (5007) Fetal Cranium Section                      | ANAP              | Ref. Section 9.1.1.7  |
| 12 | CONTAINS        | INCLUDE   | DTID (5011) Early Gestation Section                    | ANAP              | Ref. Section 9.1.1.8  |
| 13 | CONTAINS        | INCLUDE   | DTID (5009) Fetal Biophysical Profile Section          | ANAP              | Ref. Section 9.1.1.9  |
| 14 | CONTAINS        | INCLUDE   | DTID (5010) Amniotic Sac Section                       | ANAP              | Ref. Section 9.1.1.10 |
| 15 | CONTAINS        | INCLUDE   | DTID (5015) Pelvis and Uterus Section                  | ANAP              | Ref. Section 9.1.1.11 |
| 16 | CONTAINS        | INCLUDE   | DTID (5012) Ovaries Section                            | ANAP              | Ref. Section 9.1.1.12 |
| 17 | CONTAINS        | INCLUDE   | DTID (5013) Follicles Section                          | ANAP              | Ref. Section 9.1.1.13 |
| 18 | CONTAINS        | INCLUDE   | DTID (SM99003) Cyst Section                            | ANAP              | Ref. Section 9.1.1.14 |
| 19 | CONTAINS        | INCLUDE   | DTID (5025) OB-GYN Fetal Vascular Measurement Group    | ANAP              | Ref. Section 9.1.1.15 |
| 20 | CONTAINS        | INCLUDE   | DTID (5026) OB-GYN Pelvic Vascular                     | ANAP              | Ref. Section 9.1.1.16 |
| 21 | CONTAINS        | INCLUDE   | DTID (SM99005) OB-GYN Mass and Flow Section            | ANAP              | Ref. Section 9.1.1.17 |

|    |          |         |   |      |                       |
|----|----------|---------|---|------|-----------------------|
| 22 | CONTAINS | INCLUDE | DTID (5025) OB-GYN Fetal Vascular Measurement Group | ANAP | Ref. Section 9.1.1.18 |
|----|----------|---------|---|------|-----------------------|

### 9.1.1.1 Observation ConText (TID 1001)

**Table 9.1-2  
OBSERVATION CONTEXT IN OB-GYN SR**

|     | REL             | VT    | Concept Name                          | Unit / CODE Value   | Label                |
|-----|-----------------|-------|---------------------------------------|---|----------------------|
| A-1 | HAS OBS CONTEXT | CODE  | (121005, DCM, "Observer Type")        | (121006, DCM, "Person")   |                      |
| A-2 | HAS OBS CONTEXT | PNAME | (121008, DCM, "Person Observer Name") |   | Ref. Physician       |
| A-3 | HAS OBS CONTEXT | CODE  | (121024, DCM, "Subject Class")        | (121025 ,DCM,"Patient")   |                      |
| A-4 | HAS OBS CONTEXT | PNAME | (121029,DCM, "Subject Name")          |   | Last Name,First Name |
| A-5 | HAS OBS CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")    |   | BirthDate            |
| A-6 | HAS OBS CONTEXT | CODE  | (121032,DCM, "Subject Sex")           | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown sex") | M<br>F<br>O          |
| A-7 | HAS OBS CONTEXT | NUM   | (121033,DCM, "Subject Age")           | (mo, UCUM, "month")   | Not Used             |

### 9.1.1.2 Patient Characteristics (TID 5001)

**Table 9.1-3  
PATIENT CHARACTERISTICS IN OB-GYN SR**

|       | REL      | VT        | Concept Name                           | Unit / CODE Value                                    | Label       |
|-------|----------|-----------|--|--|-------------|
| A-8   | CONTAINS | CONTAINER | (121118,DCM "Patient Characteristics") |  |             |
| A-8-1 | CONTAINS | TEXT      | (121106,DCM, "Comment")                |  | Description |
| A-8-2 | CONTAINS | NUM       | (8302-2, LN, "Patient Height")         | (cm, UCUM, "centimeter")<br>(mm, UCUM, "millimeter") | Height      |

|       |          |     |                                      |                         |         |
|-------|----------|-----|--------------------------------------|-------------------------|---------|
| A-8-3 | CONTAINS | NUM | (29463-7, LN, "Patient Weight")      | (kg, UCUM, "kilograms") | Weight  |
| A-8-4 | CONTAINS | NUM | (11996-6, LN "Gravida")              | (1, UCUM, "no units")   | Gravida |
| A-8-5 | CONTAINS | NUM | (11977-6, LN, "Para")                | (1, UCUM, "no units")   | Para    |
| A-8-6 | CONTAINS | NUM | (11612-9, LN, "Aborta")              | (1, UCUM, "no units")   | Aborta  |
| A-8-7 | CONTAINS | NUM | (33065-4, LN, "Ectopic Pregnancies") | (1, UCUM, "no units")   | Ectopic |

### 9.1.1.3 OB-GYN Summary Section (TID 5002)

**Table 9.1-4**

#### **OB-GYN Procedure Summary Section**

|         | REL             | VT        | Concept Name                            | Unit / CODE Value     | Label           | Comments     |
|---------|-----------------|-----------|---|-----------------------|-----------------|--------------|
| A-9     | CONTAINS        | CONTAINER | (121111, DCM, "Summary")                |                       |                 |              |
| A-9-1   | CONTAINS        | DATE      | Context ID 12003 Extended OB-GYN Dates  | yyyymmdd              | Estab. Due Date | Table 9.1-20 |
| A-9-2   | CONTAINS        | NUM       | (11878-6, LN, "Number of Fetuses")      | (1, UCUM, "no units") |                 |              |
| A-9-3   | CONTAINS        | TEXT      | (121106, DCM, "Comment")                |                       | Comment         |              |
| A-9-4   | CONTAINS        | CONTAINER | (125008, DCM, "Fetus Summary")          |                       |                 |              |
| A-9-4-1 | HAS OBS CONTEXT | TEXT      | (11951-1, LN, "Fetus ID")               |                       |                 |              |
| A-9-4-2 | CONTAINS        | NUM       | (11878-6, LN, "Number of Fetuses")      |                       | Gestations      |              |
| A-9-4-3 | CONTAINS        | NUM       | (18185-9, LN, "Gestational Age")        |                       | GA              |              |
|         |                 |           | (11885-1, LN, "Gestational Age by LMP") |                       | GA(LMP)         |              |
| A-9-4-4 | CONTAINS        | NUM       | (11727-5, LN, "Estimated Weight")       | (kg, UCUM, "kg")      | EFW             |              |

|           |               |      |   |   |            |                   |
|-----------|---------------|------|---|---|------------|-------------------|
| A-9-4-4-1 | INFERRED FROM | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values") | (Context ID 12014) OB Fetal Body Weight Equations and Tables              |            | Ref. Table 9.1-38 |
| A-9-4-5   | CONTAINS      | NUM  | (11767-1, LN, "EFW percentile rank")                          | (%, UCUM, "Percent")  | Pctl.(EFW) |                   |
| A-9-4-5-1 | INFERRED FROM | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values") | (Context ID 12016) Estimated Fetal Weight Percentile Equations and Tables |            | Ref. Table 9.1-40 |
| A-9-4-6   | CONTAINS      | NUM  | (11948-7, LN, "Fetal Heart Rate")                             | (bpm, UCUM "bpm")   | FHR        |                   |

**9.1.1.4 OB-GYN Fetal Biometry Ratio Section (TID 5004)**

**Table 9.1-5**

**Fetal Biometry Ratio Section in OB-GYN SR**

|        | REL             | VT        | Concept Name   | Unit / CODE Value     | Comments          |
|--------|-----------------|-----------|--|-----------------------|-------------------|
| A-10   | CONTAINS        | CONTAINER | (125001, DCM, "Fetal Biometry Ratios")                       |                       |                   |
| A-10-1 | HAS OBS CONTEXT | TEXT      | (11951-1, LN, "Fetus ID")                                    |                       |                   |
| A-10-2 | CONTAINS        | NUM       | (12004, CID, "Fetal Biometry Ratios")                        | (1, UCUM, "no units") | Ref. Table 9.1-21 |
|        | CONTAINS        | NUM       | Context ID 12004 Extended Fetal Biometry Ratios Measurements | (%, UCUM, "%")        |                   |

### 9.1.1.5 OB-GYN Fetal Biometry Section (TID 5005)

**Table 9.1-6**

#### **Fetal Biometry Section in OB-GYN SR**

|            | REL                   | VT        | Concept Name   | Unit / CODE Value   | Comments             |
|------------|-----------------------|-----------|--|---|----------------------|
| A-11       | CONTAINS              | CONTAINER | (125002, DCM, "Fetal Biometry")                                |   |                      |
| A-11-1     | HAS OBS<br>CONTEXT    | TEXT      | (11951-1, LN, "Fetus ID")                                      |   |                      |
| A-11-2     | CONTAINS              | CONTAINER | (125005, DCM, "Biometry Group")                                |   |                      |
| A-11-2-1   | CONTAINS              | NUM       | Context ID 12005 Extended Fetal<br>Biometry Measurements       | (cm, UCUM,<br>"centimeter")                                     | Ref. Table<br>9.1-22 |
|            |                       |           |  | (mm, UCUM,<br>"millimeter")                                     |                      |
|            |                       |           |  | (cm <sup>2</sup> , UCUM,<br>"Square<br>centimeter")             |                      |
| A-11-2-1-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM, "Derivation")                                    | Common CID-<br>Derivation                                       |                      |
| A-11-2-2   | CONTAINS              | NUM       | (18185-9, LN, "Gestational Age")                               | (d, UCUM, "days")   |                      |
| A-11-2-2-1 | INFERRED<br>FROM      | CODE      | (121420 , DCM, "Equation")<br>(121424, DCM, "Table of Values") | (ContextID 12013)<br>Gestational Age<br>Equations and<br>Tables | Ref. Table<br>9.1-38 |

|            |               |      |  |  |                   |
|------------|---------------|------|--|--|-------------------|
| A-11-2-3   | CONTAINS      | NUM  | (125012, DCM, "Growth Percentile Rank")<br>(125013, DCM, "Growth Z-score") | (percentile, UCUM, "percentile")                       |                   |
| A-11-2-3-1 | INFERRED FROM | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values")              | (ContextID 12015)<br>Fetal Growth Equations and Tables | Ref. Table 9.1-39 |

### 9.1.1.6 OB-GYN Fetal Long Bones Section (TID 5006)

**Table 9.1-7**

#### **Fetal Long Bones Sections in OB-GYN SR**

|            | REL             | VT        | Concept Name   | Unit / CODE Value                          | Comments                                |
|------------|-----------------|-----------|--|--|---|
| A-12       | CONTAINS        | CONTAINER | (125003, DCM, "Fetal Long Bones")                                |  |   |
| A-12-1     | HAS OBS CONTEXT | TEXT      | (11951-1, LN, "FetusID")   |  | Will be present if more than one fetus. |
| A-12-2     | CONTAINS        | CONTAINER | (125005, DCM, "Biometry Group")                                  |  |   |
| A-12-2-1   | CONTAINS        | NUM       | Context ID 12006 Extended Fetal Long Bones Biometry Measurements | (cm, UCUM, "centimeter")                   | Ref. Table 9.1-23                       |
| A-12-2-1-1 | HAS CONCEPT MOD | CODE      | (121401, DCM, "Derivation")                                      | Common CID-Derivation                      |   |
| A-12-2-1-2 | HAS CONCEPT MOD | CODE      | (G-C0E3, SRT, "Finding Site")                                    | (ContextID 7304)<br>Implant Target Anatomy | Ref. Table 9.1-34                       |

|              |                       |      |   |   |                      |
|--------------|-----------------------|------|---|---|----------------------|
| A-12-2-1-2-1 | HAS<br>CONCEPT<br>MOD | CODE | (G-C171, SRT, "Laterality")                                   | (G-A100, SRT, "Right")<br>(G-A101, SRT, "Left")                 |                      |
| A-12-2-2     | CONTAINS              | NUM  | (18185-9, LN, "Gestational Age")                              | (d, UCUM, "day")  |                      |
| A-12-2-2-1   | INFERRED<br>FROM      | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values") | (ContextID 12013)<br>Gestational Age<br>Equations and<br>Tables | Ref. Table<br>9.1-37 |
| A-12-2-3     | CONTAINS              | NUM  | (125012, DCM, "Growth Percentile Rank")                       | (percentile, UCUM, "percentile")                                |                      |
| A-12-2-3-1   | INFERRED<br>FROM      | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values") | (ContextID 12015)<br>Fetal Growth<br>Equations and<br>Tables    | Ref. Table<br>9.1-39 |

### 9.1.1.7 OB-GYN Fetal Cranium Section (TID 5007)

**Table 9.1-8**

#### **Fetal Cranium Sections in OB-GYN SR**

|          | REL                | VT        | Concept Name                            | Unit / CODE Value  | Comments                                |
|----------|--------------------|-----------|---|--|---|
| A-13     | CONTAINS           | CONTAINER | (125004, DCM, "Fetal Cranium")          |  |   |
| A-13-1   | HAS OBS<br>CONTEXT | TEXT      | (11951-1, LN, "FetusID")                |  | Will be present if more than one fetus. |
| A-13-2   | CONTAINS           | CONTAINER | (125005, DCM, "Biometry Group")         |  |   |
| A-13-2-1 | CONTAINS           | NUM       | Context ID 12007 Extended Fetal Cranium | (cm, UCUM, "centimeter")<br>(mm, UCUM, "millimeter")<br>(cm2, UCUM, "Square centimeter") | Ref. Table 9.1-24                       |

|                  |                       |      |  |  |                   |
|------------------|-----------------------|------|--|--|-------------------|
| A-13-2-1-1       | HAS<br>CONCEPT<br>MOD | CODE | (121401, DCM, "Derivation")                                      | Common CID-<br>Derivation                                    |                   |
| A-13-2-1-2       | HAS<br>CONCEPT<br>MOD | CODE | (G-C0E3, SRT, "Finding Site")                                    | (ContextID 12022)<br>Fetal Cranium<br>Anatomic Sites         | Table 9.1-35      |
| A-13-2-1-2-<br>1 | HAS<br>CONCEPT<br>MOD | CODE | (G-C171, SRT, "Laterality")                                      | (G-A100, SRT,<br>"Right")<br>(G-A101, SRT, "Left")           |                   |
| A-13-2-2         | CONTAINS              | NUM  | (18185-9, LN, "Gestational<br>Age")                              | (d, UCUM, "day")   |                   |
| A-13-2-2-1       | INFERRED<br>FROM      | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of<br>Values") | (ContextID 12013)<br>Gestational Age<br>Equations and Tables | Ref. Table 9.1-36 |
| A-13-2-3         | CONTAINS              | NUM  | (125012, DCM, "Growth<br>Percentile Rank")                       | (percentile, UCUM,<br>"percentile")                          |                   |
| A-13-2-3-1       | INFERRED<br>FROM      | CODE | (121420, DCM, "Equation")<br>(121424, DCM, "Table of<br>Values") | (ContextID 12015)<br>Fetal Growth<br>Equations and Tables    | Ref. Table 9.1-39 |

### 9.1.1.8 OB-GYN Early Gestation Section (TID 5011)

**Table 9.1-9**

#### Early Gestation Section in OB-GYN SR

|        | REL                | VT        | Concept Name                        | Unit / CODE Value | Label | Comments |
|--------|--------------------|-----------|-------------------------------------|-------------------|-------|----------|
| A-14   | CONTAINS           | CONTAINER | (125009, DCM, "Early<br>Gestation") |                   |       |          |
| A-14-1 | HAS OBS<br>CONTEXT | TEXT      | (11951-1, LN, "Fetus<br>ID")        |                   |       |          |

|            |                       |           |  |  |  |                   |
|------------|-----------------------|-----------|--|--|--|-------------------|
| A-14-2     | CONTAINS              | CONTAINER | (125005, DCM, "Biometry Group")  |  |  |                   |
| A-14-2-1   | CONTAINS              | NUM       | Context ID 12009<br>Extended Early<br>Gestation Biometry<br>Measurements | (cm, UCUM, "centimeter")<br>(mm, UCUM, "millimeter")<br>(cm2, UCUM, "Square centimeter") |  | Ref. Table 9.1-26 |
| A-14-2-1-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM, "Derivation")  | Common CID-Derivation  |  |                   |
| A-14-2-2   | CONTAINS              | NUM       | (18185-9, LN Gestational Age   | (d, UCUM, "day")   |  |                   |
| A-14-2-2-1 | INFERRED<br>FROM      | CODE      | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values")            | (ContextID 12013)Gestational Age Equations and Tables                                    |  | Ref. Table 9.1-36 |
| A-14-2-3   | CONTAINS              | NUM       | (125012, DCM, "Growth Percentile Rank")                                  | (percentile, UCUM, "percentile")   |  |                   |
| A-14-2-3-1 | INFERRED<br>FROM      | CODE      | (121420, DCM, "Equation")<br>(121424, DCM, "Table of Values")            | (ContextID 12015)Fetal Growth Equations and Tables                                       |  | Ref. Table 9.1-39 |

### 9.1.1.9 OB-GYN Fetal Biophysical Profile Section (TID 5009)

Table 9.1-10

#### Fetal Biophysical Profile Section in OB-GYN SR

|      | REL      | VT        | Concept Name                         | Unit / CODE Value | Label | Comments |
|------|----------|-----------|--------------------------------------|-------------------|-------|----------|
| A-15 | CONTAINS | CONTAINER | (125006, DCM, "Biophysical Profile") |                   |       |          |

|        |                    |      |  |                               |                                 |                                       |
|--------|--------------------|------|--|-------------------------------|---------------------------------|---------------------------------------|
| A-15-1 | HAS OBS<br>CONTEXT | TEXT | (11951-1, LN, "Fetus<br>ID")                         |                               |                                 |                                       |
| A-15-2 | CONTAINS           | NUM  | (11631-9, LN, "Gross<br>Body Movement")              | ({0:2}, UCUM, "range<br>0:2") | Fetal<br>Movements              | value as<br>entered in<br>the Report. |
|        |                    |      | (11632-7, LN, "Fetal<br>Breathing")                  |                               | Fetal<br>Breathing<br>Movements |                                       |
|        |                    |      | (11635-0, LN, "Fetal<br>Tone")                       |                               | Fetal Tone                      |                                       |
|        |                    |      | (11635-5, LN, "Fetal<br>Heart Reactivity")           |                               | Nonstress<br>Test               |                                       |
|        |                    |      | (11630-1, LN,<br>"Amniotic Fluid<br>Volume")         |                               | Amniotic<br>Fluid<br>Volume     |                                       |
|        |                    |      | (11634-3, LN,<br>"Biophysical Profile<br>Sum Score") |                               | (1, UCUM, "no units")           | Total                                 |

#### 9.1.1.10 OB-GYN Amniotic Sac Section (TID 5010)

Table 9.1-11

#### Amniotic Sac Section in OB-GYN SR

|        | REL                   | VT        | Concept Name                                | Unit / CODE<br>Value                 | Label | Comments |
|--------|-----------------------|-----------|---|--------------------------------------|-------|----------|
| A-16   | CONTAINS              | CONTAINER | (121070, DCM,<br>"Findings")                |                                      | AFI   |          |
| A-16-1 | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT,<br>"Finding Site")            | (T-F1300,<br>SRT, "Amniotic<br>Sac") |       |          |
| A-16-2 | CONTAINS              | NUM       | (11627-7, LN,<br>"Amniotic Fluid<br>Index") | (cm, UCUM,<br>"centimeter")          | AFI   |          |

|          |                       |      |  |                           |                        |                      |
|----------|-----------------------|------|--|---------------------------|------------------------|----------------------|
|          |                       |      | Context ID 12008<br>Extended OB-<br>GYN Amniotic Sac |                           |                        | Ref. Table<br>9.1-25 |
|          |                       |      | (99004-01, MDSN,<br>"MVP")                           |                           | Max Vertical<br>Pocket |                      |
| A-16-2-1 | HAS<br>CONCEPT<br>MOD | CODE | (121401, DCM,<br>"Derivation")                       | Common CID-<br>Derivation |                        |                      |

### 9.1.1.11 OB-GYN Pelvis and Uterus Section (TID 5015)

**Table 9.1-12**

#### **Pelvis and Uterus Section in OB-GYN SR**

|                | REL                   | VT        | Concept Name                             | Unit / CODE<br>Value                               | Label       | Comments |
|----------------|-----------------------|-----------|--|--|-------------|----------|
| A-17           | CONTAINS              | CONTAINER | (125011, DCM,<br>"Pelvis and<br>Uterus") |  |             |          |
| A-17-1         | CONTAINS              | CONTAINER | (T-83000, SRT,<br>"Uterus")              |  |             |          |
| A-17-1-1       | CONTAINS              | NUM       | (11865-3, LN,<br>"Uterus Width")         | (cm, UCUM,<br>"centimeter")                        | Uterus W    |          |
|                |                       |           | (11842-2, LN,<br>"Uterus Length")        |  | Uterus L    |          |
|                |                       |           | (11859-6, LN,<br>"Uterus Height")        |  | Uterus H    |          |
| A-17-1-1-<br>1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM,<br>"Derivation")           | Common CID-<br>Derivation                          |             |          |
| A-17-1-2       | CONTAINS              | NUM       | (33192-6, LN,<br>"Uterus Volume")        | (cm <sup>3</sup> , UCUM,<br>"Cubic<br>centimeter") | Uterus Vol. |          |

|          |                       |      |   |  |             |                      |
|----------|-----------------------|------|---|--|-------------|----------------------|
| A-17-2   | CONTAINS              | NUM  | Context ID 12011<br>Extended<br>Ultrasound Pelvis<br>and Uterus | (cm, UCUM,<br>"centimeter")<br>(mm, UCUM,<br>"millimeter") |             | Ref. Table<br>9.1-27 |
| A-17-2-1 | HAS<br>CONCEPT<br>MOD | CODE | (121401, DCM,<br>"Derivation")                                  | Common CID-<br>Derivation                                  |             |                      |
| A-17-3   | CONTAINS              | NUM  | Cervix Volume   | (ml, UCUM,<br>"milliliter")                                | Cervix Vol. |                      |

### 9.1.1.12OB-GYN Ovaries Section (TID 5012)

**Table 9.1-13**  
**Ovary in OB-GYN SR**

|          | REL                   | VT        | Concept Name                             | Unit / CODE<br>Value        | Label               | Comments |
|----------|-----------------------|-----------|--|-----------------------------|---------------------|----------|
| A-18     | CONTAINS              | CONTAINER | (121070, DCM,<br>"Findings")             |                             |                     |          |
| A-18-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT,<br>"Finding Site")         | (T-87000,<br>SRT, "Ovary")  |                     |          |
| A-18-2   | CONTAINS              | CONTAINER | (T-87000, SRT,<br>"Ovary")               |                             |                     |          |
| A-18-2-1 | CONTAINS              | NUM       | (11829-9, LN,<br>"Left Ovary<br>Width")  | (cm, UCUM,<br>"centimeter") | Lt. Ovary<br>Width  |          |
|          |                       |           | (11840-6, LN,<br>"Left Ovary<br>Length") |                             | Lt. Ovary<br>Length |          |
|          |                       |           | (11857-0, LN,<br>"Left Ovary<br>Height") |                             | Lt. Ovary<br>Height |          |
|          |                       |           | (99005-29, MDSN,<br>"Left Ovary Area")   | Lt. Ovary<br>Area           |                     |          |

|            |                       |           |  |  |                     |  |
|------------|-----------------------|-----------|--|--|---------------------|--|
| A-18-2-1-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM,<br>"Derivation")             | Common CID-<br>Derivation                                  |                     |  |
| A-18-2-2   | CONTAINS              | NUM       | (12164-0, LN,<br>"Left Ovary<br>Volume")   | (cm3, UCUM,<br>"Cubic<br>centimeter")                      | Lt. Ovary<br>Volume |  |
| A-18-3     | CONTAINS              | CONTAINER | (T-87000, SRT,<br>"Ovary")                 |  | Right Ovary         |  |
| A-18-3-1   | CONTAINS              | NUM       | (11830-7, LN,<br>"Right Ovary<br>Width")   | (cm, UCUM,<br>"centimeter")<br>(mm, UCUM,<br>"millimeter") | Rt. Ovary W         |  |
|            |                       |           | (11841-4, LN,<br>"Right Ovary<br>Length")  |  | Rt. Ovary L         |  |
|            |                       |           | (11858-8, LN,<br>"Right Ovary<br>Height")  |  | Rt. Ovary H         |  |
|            |                       |           | (99005-28, MDSN,<br>"Right Ovary<br>Area") | (cm2,<br>UCUM, "Square<br>centimeter")                     | Rt. Ovary<br>Area   |  |
| A-18-3-1-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM,<br>"Derivation")             | Common CID-<br>Derivation                                  |                     |  |
| A-18-3-2   | CONTAINS              | NUM       | (12165-7, LN,<br>"Right Ovary<br>Volume")  | (cm3, UCUM,<br>"Cubic<br>centimeter")                      | Rt. Ovary<br>Volume |  |

### 9.1.1.13 OB-GYN Follicles Section (TID 5013)

**Table 9.1-14**

#### **Follicle Section in OB-GYN SR**

|        | REL            | VT        | Concept Name                     | Unit / CODE Value                     | Label | Comments |
|--------|----------------|-----------|----------------------------------|---------------------------------------|-------|----------|
| A-19   | CONTAINS       | CONTAINER | (121070, DCM,<br>"Findings")     |                                       |       |          |
| A-19-1 | HAS<br>CONCEPT | CODE      | (G-C0E3, SRT,<br>"Finding Site") | (T-87600, SRT,<br>"Ovarian Follicle") |       |          |

|            |                       |           |  |  |              |  |
|------------|-----------------------|-----------|--|--|--------------|--|
|            | MOD                   |           |  |  |              |  |
| A-19-2     | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT,<br>"Laterality")           | (G-A100, SRT,<br>"Right")<br>(G-A101, SRT, "Left")         |              |  |
| A-19-3     | CONTAINS              | CONTAINER | (125007, DCM,<br>"Measurement<br>Group") |  |              |  |
| A-19-3-1   | HAS OBS<br>CONTEXT    | TEXT      | (12510, DCM,<br>"Identifier")            |  | "1", "2" ... |  |
| A-19-3-2   | CONTAINS              | NUM       | (GD705, SRT,<br>"Volume")                | (ml, UCUM, "milliliter")                                   | Vol.         |  |
| A-19-3-3   | CONTAINS              | NUM       | (11793-7, LN, "Follicle<br>Diameter")    | (cm, UCUM,<br>"centimeter")<br>(mm, UCUM,<br>"millimeter") | "1", "2" ... |  |
| A-19-3-3-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM,<br>"Derivation")           | Common CID-<br>Derivation                                  |              |  |

### 9.1.1.14OB-GYN Cyst Section Section (TID SM99003)

Table 9.1-15

#### Cyst Section in OB-GYN SR

|          | REL                   | VT        | Concept Name                     | Unit / CODE Value                                     | Label       | Comments |
|----------|-----------------------|-----------|----------------------------------|---|-------------|----------|
| A-20     | CONTAINS              | CONTAINER | (121070, DCM,<br>"Findings")     |   |             |          |
| A-20-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding<br>Site") | (99009-01 , MDSN,<br>"Cyst")                          |             |          |
| A-20-2   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT,<br>"Laterality")   | (G-A100, SRT,<br>"Right")<br>(G-A101, SRT,<br>"Left") |             |          |
| A-20-3   | CONTAINS              | CONTAINER | (99009-01 , MDSN,<br>"Cyst")     |   |             |          |
| A-20-3-1 | CONTAINS              | NUM       | (99005-26, MSDN,                 | (cm, UCUM,  | Ovary Width |          |

|            |                       |      |                                 |                          |              |  |
|------------|-----------------------|------|---------------------------------|--------------------------|--------------|--|
|            |                       |      | "Cyst Width")                   | "centimeter")            |              |  |
|            |                       |      | (99005-24, MSDN, "Cyst Length") | (mm, UCUM, "millimeter") | Cyst Length  |  |
|            |                       |      | (99005-25, MSDN, "Cyst Height") |                          | Ovary Height |  |
| A-20-3-1-1 | HAS<br>CONCEPT<br>MOD | CODE | (121401, DCM, "Derivation")     | Common CID-Derivation    |              |  |
| A-20-3-2   | CONTAINS              | NUM  | (99005-27, MSDN, "Cyst Volume") | (ml, UCUM, "milliliter") | Ovary Volume |  |

### 9.1.1.15 OB-GYN Fetal Vascular Measurement Group (TID 5025)

Table 9.1-16

#### OB-GYN Fetal Vascular Measurement Group in OB-GYN SR

|          | REL                   | VT        | Concept Name  | Unit / CODE Value                              | Label | Comments          |
|----------|-----------------------|-----------|---|--|-------|-------------------|
| A-21     | CONTAINS              | CONTAINER | (121070, DCM, "Findings")   |  |       |                   |
| A-21-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding Site")   | (T-F6800, SRT, "Embryonic Vascular Structure") |       |                   |
| A-21-2   | CONTAINS              | CONTAINER | Context ID 12141<br>Extended Fetal Vasculature<br>Anatomical Location |  |       | Ref. Table 9.1-28 |
| A-21-2-1 | HAS OBS<br>CONTEXT    | TEXT      | (11951-1, LN, "Fetus ID")   |  |       |                   |
| A-21-2-2 | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT, "Laterality")   | (G-A100, SRT, "Right")                         |       |                   |
|          |                       |           |   | (G-A101, SRT, "Left")                          |       |                   |
|          |                       |           |   | (G-A102, SRT, "Unilateral")                    |       |                   |
| A-21-2-3 | CONTAINS              | NUM       | Context ID 12119<br>Vascular Ultrasound<br>Property                   |  |       | Ref. Table 9.1-30 |

|            |                       |      |                                |                           |  |  |
|------------|-----------------------|------|--------------------------------|---------------------------|--|--|
| A-21-2-3-1 | HAS<br>CONCEPT<br>MOD | CODE | (121401, DCM,<br>"Derivation") | Common CID-<br>Derivation |  |  |
|------------|-----------------------|------|--------------------------------|---------------------------|--|--|

### 9.1.1.16 OB-GYN Pelvic Vascular Measurement Group (TID 5026)

Table 9.1-17

#### OB-GYN Pelvic Vascular Measurement Group in OB-GYN SR

|            | REL                   | VT        | Concept Name  | Unit / CODE Value                                 | Label        | Comments             |
|------------|-----------------------|-----------|---|---|--------------|----------------------|
| A-22       | CONTAINS              | CONTAINER | (121070, DCM,<br>"Findings")  |   |              |                      |
| A-22-1     | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding<br>Site")  | (T-D6007, SRT,<br>"Pelvic Vascular<br>Structure") |              |                      |
| A-22-2     | CONTAINS              | CONTAINER | Context ID 12140<br>Extended Pelvic<br>Vasculature Anatomical<br>Location |   |              | Ref. Table<br>9.1-29 |
| A-22-2-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT,<br>"Laterality")  | (SRT, G-A100,<br>"Right")                         |              |                      |
|            |                       |           |   | (SRT, G-A101,<br>"Left")                          |              |                      |
|            |                       |           |   | (SRT, G-A103,<br>"Unilateral")                    |              |                      |
| A-22-2-2   | HAS<br>CONCEPT<br>MOD | TEXT      | (112050, DCM,<br>"Anatomic Identifier")                                   |   | "1", "2" ... |                      |
| A-22-2-3   | CONTAINS              | NUM       | Context ID 12119<br>Vascular Ultrasound<br>Property                       |   |              | Ref. Table<br>9.1-30 |
| A-22-2-3-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM,<br>"Derivation")  | Common CID-<br>Derivation                         |              |                      |

### 9.1.1.17 OB-GYN Mass and Flow Section (TID SM99005)

Table 9.1-18

### Mass and Flow Section in OB-GYN SR

|            | REL                   | VT        | Concept Name                                     | Unit / CODE Value           | Label                    | Comments                        |
|------------|-----------------------|-----------|--|-----------------------------|--------------------------|---------------------------------|
| A-23       | CONTAINS              | CONTAINER | (121070, DCM, "Findings")                        |                             |                          |                                 |
| A-23-1     | HAS<br>CONCEPT<br>MOD | CODE      | Context ID 99100<br>Gynecology Finding Site      |                             |                          | <b>Ref<br/>Table<br/>9.1-36</b> |
| A-23-2     | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT, "Laterality")                      | (G-A100, SRT, "Right")      |                          |                                 |
|            |                       |           |  | (G-A101, SRT, "Left")       |                          |                                 |
|            |                       |           |  | (G-A102, SRT, "Unilateral") |                          |                                 |
| A-23-3     | CONTAINS              | CONTAINER | Context ID 99103<br>Gynecology Mass and Flow     |                             |                          | Ref Table 9.1-36                |
| A-23-3-1   | HAS OBS<br>CONTEXT    | TEXT      | (12510, DCM, "Identifier")                       |                             | "1", "2" ...             |                                 |
| A-23-3-2   | CONTAINS              | NUM       | (11840-6, LN, "Length")                          | (cm, UCUM, "centimeter")    |                          |                                 |
|            | CONTAINS              | NUM       | (11857-0, LN, "Height")                          |                             |                          |                                 |
|            | CONTAINS              | NUM       | (11829-9, LN, "Width")                           |                             | (mm, UCUM, "millimeter") |                                 |
| A-23-3-2-1 | HAS<br>CONCEPT<br>MOD | CODE      | (121401, DCM, "Derivation")                      | Common CID-Derivation       |                          |                                 |
| A-23-3-3   | CONTAINS              | NUM       | (G-D705, SRT, "Volume")                          | (ml, UCUM, "milliliter")    | Vol.                     |                                 |
| A-23-3-4   | CONTAINS              | NUM       | Context ID 12119<br>Vascular Ultrasound Property |                             |                          | Ref. Table 9.1-30               |
| A-23-3-4-1 | HAS<br>CONCEPT        | CODE      | (121401, DCM, "Derivation")                      | Common CID-Derivation       |                          |                                 |

|  |     |  |  |  |  |  |
|--|-----|--|--|--|--|--|
|  | MOD |  |  |  |  |  |
|--|-----|--|--|--|--|--|

**9.1.1.18 OB-GYN User Creation Group Section (TID SM99010)**

**Table 9.1-19**

**User Creation Group Section in OB-GYN SR**

|          | REL             | VT        | Concept Name  | Unit / CODE Value  | Label        | Comments  |
|----------|-----------------|-----------|---|--|--------------|---|
| A-24     | CONTAINS        | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |  |              | *Creation ID: Randomly generated 7-digit unique ID          |
| A-24-1   | HAS OBS CONTEXT | TEXT      | (11951-1, LN, "Fetus ID")                             |  | "1", "2" ... | Will be present if the creation group has fetus             |
|          | HAS CONCEPT MOD | CODE      | (G-C171, SRT, "Laterality")                           | (G-A100, SRT, "Right")<br>(G-A101, SRT, "Left")<br>(G-A102, SRT, "Unilateral") |              | Will be present if the creation group has laterality option |
| A-24-2   | CONTAINS        | CONTAINER | (99900-Creation ID, MDSN, "User Creation Item Name")  |  |              |   |
| A-24-2-1 | HAS CONCEPT MOD | CODE      | (121401, DCM, "Derivation")                           | Common CID-Derivation  |              |   |

## 9.1.2 DCMR Context Groups used in V7

### 9.1.2.1 Standard Extended Context Groups in OB-GYN SR

**Table 9.1-20**  
**Context ID 12003 Extended OB-GYN Dates**

| CSD | CV      | CM                              | Label         |
|-----|---------|---------------------------------|---------------|
| LN  | 11778-8 | EDD                             | Estab.DueDate |
| LN  | 11779-6 | EDD from LMP                    | EDD(LMP)      |
| LN  | 11781-2 | EDD from average ultrasound age | EDD(AUA)      |
| LN  | 11955-2 | LMP                             | LMP           |
| LN  | 11976-8 | Ovulation date                  | Exp.Ovul.     |

**Table 9.1-21**  
**Context ID 12004 Extended Fetal Biometry Ratios Measurements**

| CSD  | CV       | CM  | Laterality | Label         |
|------|----------|---|------------|---------------|
| LN   | 11871-1  | FL/AC   | N/A        | FL/AC         |
| LN   | 11872-9  | FL/BPD  | N/A        | FL/BPD        |
| LN   | 11873-7  | FL/HC   | N/A        | FL/HC         |
| LN   | 11823-2  | Cephalic Index  | N/A        | CI(BPD/OFD)   |
| LN   | 11947-9  | HC/AC   | N/A        | HC/AC         |
| MDSN | 99000-01 | FL/FOOT   | N/A        | FL/FOOT       |
| MDSN | 99000-02 | ThC/AC  | N/A        | ThC/AC        |
| MDSN | 99000-08 | Anterior Horn Lateral ventricular to Hemispheric Width Ratio  | N/A        | Va/Hem        |
| MDSN | 99000-09 | Posterior Horn Lateral ventricular to Hemispheric Width Ratio | N/A        | Vp/Hem        |
| MDSN | 99000-14 | Lateral Ventricle width to Hemispheric Width Ratio            | N/A        | Ventricle/Hem |

**Table 9.1-22**  
**Context ID 12005 Extended Fetal Biometry Measurements**

| CSD | CV | CM | Label |
|-----|----|----|-------|
|-----|----|----|-------|

|      |          |  |              |
|------|----------|--|--------------|
| LN   | 11820-8  | Biparietal Diameter                    | BPD          |
| LN   | 11851-3  | Occipital-Frontal Diameter             | OFD          |
| LN   | 11984-2  | Head Circumference                     | HC           |
| LN   | 11818-2  | Anterior-Posterior Abdominal Diameter  | APD          |
| LN   | 11862-0  | Tranverse Abdominal Diameter           | TAD          |
| LN   | 11979-2  | Abdominal Circumference                | AC           |
| LN   | 11963-6  | Femur Length                           | FL           |
| MDSN | 99000-10 | Right Femur Length                     | Rt. FL       |
| MDSN | 99000-11 | Left Femur Length                      | Lt. FL       |
| LN   | 11819-0  | Anterior-Posterior Trunk Diameter      | APTD         |
| LN   | 11864-6  | Transverse Thoracic Diameter           | TTD          |
| LN   | 11988-3  | Thoracic Circumference                 | ThC          |
| MDSN | 99001-20 | AVol                                   | AVol         |
| MDSN | 99001-21 | TVol                                   | TVol         |
| LN   | 11965-1  | Foot length                            | Foot         |
| MDSN | 99001-18 | Kidney length                          | Renal L      |
| LN   | 11834-9  | Left Kidney length                     | Lt. Renal L  |
| LN   | 11825-7  | Left Kidney width                      | Lt. Renal AP |
| MDSN | 99001-19 | Kidney width                           | Renal AP     |
| LN   | 11836-4  | Right Kidney length                    | Rt. Renal L  |
| LN   | 11827-3  | Right Kidney width                     | Rt. Renal AP |
| LN   | 33068-8  | Thoracic Area                          | ThA          |
| MDSN | 99001-01 | Middle Abdominal Diameter              | MAD          |
| MDSN | 99005-13 | Right Pelvis                           | Rt. Pelvis   |
| MDSN | 99005-14 | Left Pelvis                            | Lt. Pelvis   |
| MDSN | 99005-01 | Pelvis                                 | Pelvis       |
| MDSN | 99001-02 | Fetal Trunk Area                       | FTA          |
| MDSN | 99001-03 | APTDxTTD                               | APTDxTTD     |
| MDSN | 99001-04 | Ear Length                             | Ear          |
| MDSN | 99001-05 | Middle Phalanx                         | MP           |
| MDSN | 99001-06 | Thoracic Anteriorposterior Diameter    | ThD ap       |
| MDSN | 99001-07 | Thoracic Transverse Diameter           | ThD trans    |
| MDSN | 99001-08 | Heart Anteriorposterior Diameter       | HrtD ap      |
| MDSN | 99001-09 | Heart Transverse Diameter              | HrtD trans   |
| MDSN | 99001-11 | Cardio-Thoracic Area Ratio by Distance | CTAR(D)      |

|      |          |   |         |
|------|----------|---|---------|
| MDSN | 99001-12 | Heart Area                                  | HrtA    |
| MDSN | 99001-13 | Cardio-Thoracic Area Ratio by Area          | CTAR(A) |
| MDSN | 99001-23 | Cardio-Thoracic Area Ratio by Circumference | CTAR(C) |

**Table 9.1-23**

**Context ID 12006 Extended Fetal Long Bones Biometry Measurements**

| CSD  | CV       | CM              | Laterality       | Label     |
|------|----------|-----------------|------------------|-----------|
| LN   | 11966-9  | Humerus length  | N/A, Right, Left | HUM       |
| LN   | 11969-3  | Ulna length     | N/A, Right, Left | ULNA      |
| LN   | 11968-5  | Tibia length    | N/A, Right, Left | TIB       |
| LN   | 11967-7  | Radius length   | N/A, Right, Left | RAD       |
| LN   | 11964-4  | Fibula length   | N/A, Right, Left | FIB       |
| LN   | 11962-8  | Clavicle length | N/A, Right, Left | CLAV      |
| MDSN | 99002-01 | Vertebral       | N/A              | Vertebral |

**Table 9.1-24**

**Context ID 12007 Extended Fetal Cranium**

| CSD  | CV       | CM                                       | Laterality       | Label        |
|------|----------|--|------------------|--------------|
| LN   | 11863-8  | Trans Cerebellar Diameter                | N/A              | CEREB        |
| LN   | 11860-4  | Cisterna Magna length                    | N/A              | CM           |
| LN   | 12146-7  | Nuchal Fold thickness                    | N/A              | NF           |
| LN   | 33069-6  | Nuchal Translucency                      | N/A              | NT           |
| LN   | 11629-3  | Outer Orbital Diameter                   | N/A              | OOD          |
| LN   | 33070-4  | Inner Orbital Diameter                   | N/A              | IOD          |
| LN   | 33197-5  | Anterior Horn Lateral ventricular width  | N/A, Right, Left | Va           |
| LN   | 33196-7  | Posterior Horn Lateral ventricular width | N/A, Right, Left | Vp           |
| LN   | 12171-5  | Lateral Ventricle width                  | N/A, Right, Left | Ventricle    |
| LN   | 12170-7  | Width of Hemisphere                      | N/A, Right, Left | HEM          |
| MDSN | 99004-02 | Frontomaxillary facial angle             | N/A              | FMF angle    |
| MDSN | 99004-03 | Corpus callosum length                   | N/A              | CC length    |
| MDSN | 99004-04 | Corpus callosum width                    | N/A              | CC width     |
| MDSN | 99004-05 | Corpus callosum thickness                | N/A              | CC thickness |

**Table 9.1-25**  
**Context ID 12008 Extended OB-GYN Amniotic Sac**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                | <b>Label</b> |
|------------|-----------|--------------------------|--------------|
| LN         | 11624-4   | First Quadrant Diameter  | Q1           |
| LN         | 11626-9   | Second Quadrant Diameter | Q2           |
| LN         | 11625-1   | Third Quadrant Diameter  | Q3           |
| LN         | 11623-6   | Fourth Quadrant Diameter | Q4           |
| LN         | 11627-7   | Amniotic Fluid Index     | AFI          |
| SRT        | M-02550   | Diameter                 | MVP          |

**Table 9.1-26**  
**Context ID 12009 Extended Early Gestation Biometry Measurements**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                | <b>Label</b> |
|------------|-----------|--------------------------|--------------|
| LN         | 11850-5   | Gestational Sac Diameter | GS           |
| LN         | 11957-8   | Crown Rump Length        | CRL          |
| LN         | 11816-6   | Yolk Sac length          | YS           |
| LN         | 33071-2   | Spine Length             | SL           |

**Table 9.1-27**  
**Context ID 12011 Extended Ultrasound Pelvis and Uterus**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>             | <b>Label</b> |
|------------|-----------|-----------------------|--------------|
| LN         | 11961-0   | Cervix Length         | Cervix L     |
| LN         | 12145-9   | Endometrium Thickness |              |
| MDSN       | 99005-02  | Cervix Height         | Cervix H     |
| MDSN       | 99005-03  | Cervix Width          | Cervix W     |
| MDSN       | 99005-04  | Cervix Volume         | Cervix Vol.  |

**Table 9.1-28**  
**Context ID 12141 Extended Fetal Vasculature Anatomical Location**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>              | <b>Label</b>               |
|------------|-----------|------------------------|----------------------------|
| SRT        | T-42000   | Aorta                  | Fetal Aorta                |
| SRT        | T-D0765   | Descending Aorta       | Dsc Aorta (in Fetal Heart) |
| SRT        | T-45600   | Middle Cerebral Artery | Mid Cereb A                |
| SRT        | T-44000   | Pulmonary Artery       | MPA (in Fetal Heart)       |

|      |          |                    |  |
|------|----------|--------------------|--|
| SNM3 | T-45010  | Carotid artery     | Lt. Fetal Carotid<br>Rt. Fetal Carotid |
| MDSN | 99008-02 | Ductus Venosus     | Ductus Venosus                         |
| MDSN | 99008-03 | Renal Artery       | Lt. Renal A<br>Rt. Renal A             |
| SRT  | T-48710  | Inferior vena cava | IVC (in Fetal Heart)                   |
| MDSN | 99008-07 | Ductus Atriosus    | Duct A                                 |
| MDSN | 99008-09 | Ascending Aorta    | Asc Aorta                              |
| SRT  | T-F1810  | Umbilical Artery   | Umbilical A                            |
| SRT  | T-F1820  | Umbilical Vein     | Umbilical V                            |

**Table 9.1-29**

**Context ID 12140 Extended Pelvic Vasculature Anatomical Location**

| CSD  | CV       | CM                           | Label  |
|------|----------|------------------------------|--|
| SRT  | T-46980  | Ovarian Artery               | Lt. Ovarian A<br>Rt. Ovarian A   |
| SRT  | T-46820  | Uterine Artery               | Lt. Uterine A (in OB or Gynecology)<br>Rt. Uterine A (in OB or Gynecology) |
| SRT  | T-F1412  | Vitelline Artery of Placenta | Placenta A   |
| MDSN | 99007-01 | Perisystic Flow              | Perisystic Flow  |
| MDSN | 99007-02 | Endometrial Flow             | Endometrial Flow   |

**Table 9.1-30**

**Context ID 12119 Vascular Ultrasound Property**

| CSD     | CV   | CM | Label |
|---------|--|----|-------|
| INCLUDE | CID 12120 Extended Blood Velocity Measurements |    |       |
| INCLUDE | CID 12121 Vascular Indices and Ratios          |    |       |
| INCLUDE | CID 12122 Other Vascular Properties            |    |       |

**Table 9.1-31**

**Context ID 12120 Extended Blood Velocity Measurement**

| CSD | CV      | CM                         | Label |
|-----|---------|----------------------------|-------|
| LN  | 11653-3 | End Diastolic Velocity     | ED    |
| LN  | 11726-7 | Peak Systolic Velocity     | PS    |
| LN  | 11665-7 | Minimum Diastolic Velocity | MD    |

|      |          |                             |        |
|------|----------|-----------------------------|--------|
| LN   | 20352-1  | Time averaged mean velocity | TAmean |
| LN   | 11692-1  | Time averaged peak velocity | TAmx   |
| LN   | 11726-7  | Peak Systolic Velocity      | S      |
| MDSN | 99008-05 | Peak Diastolic Velocity     | D      |
| LN   | 11665-7  | Minimum Diastolic Velocity  | a      |
| MDSN | 99200-03 | Max Velocity                | Vmax   |

**Table 9.1-32**  
**Context ID 12121 Vascular Indices and Ratios**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                            | <b>Label</b> |
|------------|-----------|--------------------------------------|--------------|
| LN         | 20167-3   | Acceleration Index                   | Acc          |
| SRT        | R-101BA   | Lumen Area Stenosis                  | %StA         |
| SRT        | R-101BB   | Lumen Diameter Stenosis              | %StD         |
| LN         | 12008-9   | Pulsatility Index                    | PI           |
| LN         | 12023-8   | Resistivity Index                    | RI           |
| LN         | 12144-2   | Systolic to Diastolic Velocity Ratio | S/D          |
| MDSN       | 99200-01  | Diastolic to Systolic Velocity Ratio | D/S          |
| MDSN       | 99200-00  | Systolic to Atrial Velocity Ratio    | S/a          |
| MDSN       | 99200-09  | Preload Index                        | PLI          |
| MDSN       | 99200-10  | Peak velocity index for veins        | PVIV         |

**Table 9.1-33**  
**Context ID 12122 Other Vascular Properties**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                           | <b>Label</b> |
|------------|-----------|-------------------------------------|--------------|
| LN         | 20168-1   | Acceleration Time                   | AccT         |
| LN         | 20217-6   | Deceleration Time                   | DecT         |
| MDSN       | 99200-07  | Deceleration                        | Dec          |
| SRT        | G-0364    | Vessel lumen diameter               | Dout         |
| SRT        | R-1025C   | Vessel Intimal Diameter             | Din          |
| SRT        | R-1025D   | Vessel Intimal Cross-Sectional Area | Ain          |
| SRT        | G-0365    | Vessel outside diameter             | Vesl. Dist.  |
| SRT        | G-0366    | Vessel lumen cross-sectional area   | Aout         |
| LN         | 33878-0   | Volume flow                         | Vol. Flow    |
| LN         | 20247-3   | Peak Gradient                       | PGmax        |

|      |          |               |        |
|------|----------|---------------|--------|
| LN   | 20256-4  | Mean Gradient | PGmean |
| MDSN | 99200-04 | Duration Time | Dur T  |
| LN   | 8867-4   | Heart rate    | HR     |

**Table 9.1-34**  
**Context ID 7304 Implant Target Anatomy**

| <b>CSD</b> | <b>CV</b> | <b>CM</b> | <b>Label</b> |
|------------|-----------|-----------|--------------|
| SRT        | T-12410   | Humerus   | HUM          |
| SRT        | T-12420   | Radius    | RAD          |
| SRT        | T-12430   | Ulna      | ULNA         |
| SRT        | T-12440   | Tibia     | TIB          |
| SRT        | T-12450   | Fibula    | FIB          |
| SRT        | T-12310   | Clavicle  | CLAV         |

**Table 9.1-35**  
**Context ID 12022 Fetal Cranium Anatomic Sites**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                        | <b>Label</b> |
|------------|-----------|----------------------------------|--------------|
| SRT        | T-A1700   | Anterior Horn Lateral Ventricle  | Va           |
| SRT        | T-A1710   | Posterior Horn Lateral Ventricle | Vp           |
| SRT        | T-A010F   | Cerebral hemisphere              | Hem          |

**Table 9.1-36**  
**Context ID 99100 Gynecology Finding Site**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>         | <b>HERA W10 Label</b> |
|------------|-----------|-------------------|-----------------------|
| SRT        | M-03000   | Mass              | Mass                  |
| MDSN       | 99009-04  | Endometrial Polyp | Endo. Polyp           |
| MDSN       | 99009-05  | Ovarian Mass      | Ovarian Mass          |
| MDSN       | 99009-08  | Ectopic Pregnancy | Ectopic Pregnancy     |
| MDSN       | 99009-09  | Uterine Fibroid   | Uterine Fibroid       |
| MDSN       | 99009-10  | Cervix            | Cervix Flow           |

**Table 9.1-37**  
**Context ID 99103 Gynecology Mass and Flow**

| <b>CSD</b> | <b>CV</b> | <b>CM</b> | <b>HERA W10 Label</b> |
|------------|-----------|-----------|-----------------------|
| SRT        | M-03000   | Mass      | Mass                  |

|      |                 |                        |                 |
|------|-----------------|------------------------|-----------------|
| MDSN | <b>99007-03</b> | Endometrial Polyp Flow | Endo. Polyp     |
| MDSN | 99007-04        | Ovarian Mass Flow      | Ovarian Mass    |
| MDSN | 99007-07        | Ectopic Flow           | Ectopic Flow    |
| MDSN | 99007-08        | Uterine Fibroid Flow   | Uterine Fibroid |
| MDSN | 99007-09        | Cervical Flow          | Cervix Flow     |

### 9.1.2.2 Gestational Age Equations and Tables (Context Group 12013)

**Table 9.1-37**

#### **Gestational Age Equations and Tables**

| <b>Coding Scheme Designator</b> | <b>Code Value</b> | <b>Code Meaning</b>           |
|---------------------------------|-------------------|-------------------------------|
| LN                              | 11889-3           | AC, Campbell 1975             |
| LN                              | 11892-7           | AC, Hadlock 1984              |
| LN                              | 33076-1           | AC, Shinozuka 1996            |
| LN                              | 11902-4           | BPD, Hadlock 1984             |
| LN                              | 33538-0           | BPD, Hansmann 1986            |
| LN                              | 11905-7           | BPD, Jeanty 1984              |
| LN                              | 11906-5           | BPD, Kurtz 1980               |
| LN                              | 33082-9           | BPD, Osaka 1989               |
| LN                              | 11907-3           | BPD, Sabbagha 1978            |
| LN                              | 33084-5           | BPD, Shinozuka 1996           |
| LN                              | 33086-0           | BPD-oi, Chitty 1997           |
| LN                              | 33087-8           | BPD-oo, Chitty 1997           |
| LN                              | 33088-6           | Clavical length, Yarkoni 1985 |
| LN                              | 11910-7           | CRL, Hadlock 1992             |
| LN                              | 33540-6           | CRL, Hansmann 1986            |
| LN                              | 11913-1           | CRL, Nelson 1981              |
| LN                              | 33093-6           | CRL, Osaka 1989               |
| LN                              | 33094-4           | CRL, Rempen 1991              |
| LN                              | 11914-9           | CRL, Robinson 1975            |
| LN                              | 33095-1           | CRL, Shinozuka 1996           |
| LN                              | 33096-9           | CRL, Tokyo 1986               |
| LN                              | 33098-5           | FL, Chitty 1997               |
| LN                              | 11920-6           | FL, Hadlock 1984              |

|    |         |  |
|----|---------|--|
| LN | 33541-4 | FL, Hansmann 1986                            |
| LN | 11922-2 | FL, Hohler 1982                              |
| LN | 11923-0 | FL, Jeanty 1984                              |
| LN | 33101-7 | FL, Osaka 1989                               |
| LN | 33102-5 | FL, Shinozuka 1996                           |
| LN | 33103-3 | FL, Tokyo 1986                               |
| LN | 11928-9 | GS, Hellman 1969                             |
| LN | 33107-4 | GS, Nyberg 1992                              |
| LN | 33108-2 | GS, Tokyo 1986                               |
| LN | 33110-8 | HC measured, Chitty 1997                     |
| LN | 33111-6 | HC derived, Chitty 1997                      |
| LN | 11932-1 | HC, Hadlock 1984                             |
| LN | 33543-0 | HC, Hansmann 1986                            |
| LN | 11936-2 | Humerus, Jeanty 1984                         |
| LN | 11937-0 | Humerus, Merz 1987                           |
| LN | 33117-3 | Humerus Length, Osaka 1989                   |
| LN | 33118-1 | Length of Vertebra, Tokyo 1986               |
| LN | 33120-7 | OFD, Hansmann 1986                           |
| LN | 11941-2 | Tibia, Jeanty 1984                           |
| LN | 33134-8 | TCD, Hill 1990                               |
| LN | 33138-9 | Fetal Trunk Cross-Sectional Area, Osaka 1989 |
| LN | 11944-6 | Ulna, Jeanty 1984                            |
| LN | 11929-7 | GS, Rempen 1991                              |
| LN | 33083-7 | BPD, Rempen 1991                             |

### 9.1.2.3 OB Fetal Body Weight Equations and Tables (Context ID 12014)

**Table 9.1-38**

#### **OB Fetal Body Weight Equations and Tables**

| <b>Coding Scheme Designator</b> | <b>Code Value</b> | <b>Code Meaning</b>                  |
|---------------------------------|-------------------|--------------------------------------|
| LN                              | 11756-4           | EFW by AC, Campbell 1975             |
| LN                              | 11738-2           | EFW by AC, BPD, Hadlock 1984         |
| LN                              | 11735-8           | EFW by AC, BPD, FL, Hadlock 1985     |
| LN                              | 11732-5           | EFW by AC, BPD, FL, HC, Hadlock 1985 |

|    |         |                                       |
|----|---------|---------------------------------------|
| LN | 11751-5 | EFW by AC, FL, Hadlock 1985           |
| LN | 11746-5 | EFW by AC, FL, HC, Hadlock 1985       |
| LN | 33139-7 | EFW by BPD, TTD, Hansmann 1986        |
| LN | 11739-0 | EFW by AC and BPD, Shepard 1982       |
| LN | 33140-5 | EFW by BPD, FTA, FL, Osaka 1990       |
| LN | 33144-7 | EFW by BPD, APAD, TAD, FL, Tokyo 1987 |

#### 9.1.2.4 Fetal Growth Equations and Tables (Context ID 12015)

**Table 9.1-39**

#### **Fetal Growth Equations and Tables**

| <b>Coding Scheme Designator</b> | <b>Code Value</b> | <b>Code Meaning</b>                |
|---------------------------------|-------------------|------------------------------------|
| LN                              | 33145-4           | AC by GA, ASUM 2000                |
| LN                              | 33146-2           | AC by GA, Hadlock 1984             |
| LN                              | 33147-0           | AC (measured) by GA, Chitty 1994   |
| LN                              | 33546-3           | AC (derived) by GA, Chitty 1994    |
| LN                              | 33149-6           | AC by GA, Shinozuka 1996           |
| LN                              | 33150-4           | AxT by GA, Shinozuka 1996          |
| LN                              | 33151-2           | BPD by GA, ASUM 2000               |
| LN                              | 33198-3           | BPD by GA, Hadlock 1984            |
| LN                              | 33556-2           | BPD outer-inner by GA, Chitty 1994 |
| LN                              | 33152-0           | BPD outer-outer by GA, Chitty 1994 |
| LN                              | 33156-1           | BPD by GA, Shinozuka 1996          |
| LN                              | 33155-3           | BPD by GA, Rempen 1991             |
| LN                              | 33159-5           | CRL by GA ASUM 2000                |
| LN                              | 33161-1           | CRL by GA, Rempen1991              |
| LN                              | 33160-3           | CRL by GA, Shinozuka 1996          |
| LN                              | 33164-5           | Fibula by GA, Jeanty 1983          |
| LN                              | 33165-2           | FL by GA, ASUM 2000                |
| LN                              | 33166-0           | FL by GA, Hadlock 1984             |
| LN                              | 33167-8           | FL by GA, Chitty 1994              |
| LN                              | 33171-0           | GS by GA, Rempen 1991              |
| LN                              | 33170-2           | FL by GA, Shinozuka 1996           |

|    |         |                                 |
|----|---------|---------------------------------|
| LN | 33172-8 | HC by GA, ASUM 2000             |
| LN | 33173-6 | HC by GA, Hadlock 1984          |
| LN | 33174-4 | HC derived by GA, Chitty 1994   |
| LN | 33177-7 | Humerus Length by GA, ASUM 2000 |
| LN | 33178-5 | OFD by GA, ASUM 2000            |
| LN | 33180-1 | Radius by GA, Jeanty 1983       |
| LN | 33181-9 | TCD by GA Goldstein 1987        |

### 9.1.2.5 Estimated Fetal Weight Percentile Equations and Tables (Context ID 12016)

**Table 9.1-40**

#### **Estimated Fetal Weight Percentile Equations and Tables**

| <b>Coding Scheme Designator</b> | <b>Code Value</b> | <b>Code Meaning</b>       |
|---------------------------------|-------------------|---------------------------|
| LN                              | 33183-5           | FWP by GA, Hadlock 1991   |
| LN                              | 33184-3           | FWP by GA, Williams, 1982 |
| LN                              | 33189-2           | FWP by GA, Brenner 1976   |

## 9.2 VASCULAR ULTRASOUND REPORT TEMPLATES

### 9.2.1 Vascular Ultrasound Report (TID 5100)

**Table 9.2-1  
VASCULAR ULTRASOUND REPORT TEMPLATE**

|   | Rel with Parent | VT        | Concept Name  | Comments  | Label            |
|---|-----------------|-----------|---|---|------------------|
| 1 |                 | CONTAINER | (125100, DCM, "Vascular Ultrasound Procedure Report") |   |                  |
| 2 | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observation Context                       |   |                  |
| 3 | CONTAINS        | INCLUDE   | DTID (5101) Vascular Patient Characteristics          |   |                  |
| 4 | CONTAINS        | INCLUDE   | DTID (5102) Vascular Procedure Summary Section        |   |                  |
| 5 | CONTAINS        | INCLUDE   | DTID (5103) Vascular Ultrasound Section               | \$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head")          | TCD (Unilateral) |
|   |                 |           |   | \$SectionLaterality = EV (G-A103, SRT, "Unilateral")                |                  |
|   |                 |           |   | \$Anatomy = DCID (12106) Intracranial Cerebral Vessels (unilateral) |                  |
| 6 | CONTAINS        | INCLUDE   | DTID (5103) Vascular Ultrasound Section               | \$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head")          | TCD              |
|   |                 |           |   | \$SectionLaterality = EV (G-A100, SRT, "Right")                     |                  |
|   |                 |           |   | \$Anatomy = DCID (12105) Intracranial Cerebral Vessels              |                  |
| 7 | CONTAINS        | INCLUDE   | DTID (5103) Vascular Ultrasound Section               | \$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head")          | TCD              |
|   |                 |           |   | \$SectionLaterality = EV (G-A101, SRT, "Left")                      |                  |
|   |                 |           |   | \$Anatomy = DCID (12105)  |                  |

|    |          |         |  |  |           |
|----|----------|---------|--|--|-----------|
|    |          |         |  | Intracranial Cerebral Vessels  |           |
| 8  | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-40501,<br>SRT, "Blood Vessel of Head")<br>\$Anatomy = DCID (12105)<br>Intracranial Cerebral Vessels   | TCD       |
| 9  | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-45005,<br>SRT, "Artery of neck")<br>\$SectionLaterality = EV (G-A101,<br>SRT, "Left")<br>\$Anatomy = DCID (12104)<br>Extracranial Arteries                | Carotid   |
| 10 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-45005,<br>SRT, "Artery of neck")<br>\$SectionLaterality = EV (G-A100,<br>SRT, "Right")<br>\$Anatomy = DCID (12104)<br>Extracranial Arteries               | Carotid   |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-45005,<br>SRT, "Artery of neck")<br>\$Anatomy = DCID (12104)<br>Extracranial Arteries   | Carotid   |
| 11 |          |         |  | \$SectionScope = DT (T-47040,<br>SRT, "Artery of Lower Extremity")<br>\$SectionLaterality = EV (G-A101,<br>SRT, "Left")<br>\$Anatomy = DCID (12109) Lower<br>Extremity Arteries  | LE Artery |
| 12 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-47040,<br>SRT, "Artery of Lower Extremity")<br>\$SectionLaterality = EV (G-A100,<br>SRT, "Right")<br>\$Anatomy = DCID (12109) Lower<br>Extremity Arteries | LE Artery |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-47040,<br>SRT, "Artery of Lower Extremity")   | LE Artery |

|    |          |         |   |   |           |
|----|----------|---------|---|---|-----------|
|    |          |         |   | \$Anatomy = DCID (12109) Lower Extremity Arteries               |           |
| 13 |          |         |   | \$SectionScope = DT (T-49403, SRT, "Vein of Lower Extremity")   | LE Vein   |
|    |          |         |   | \$SectionLaterality = EV (G-A101, SRT, "Left")                  |           |
|    |          |         |   | \$Anatomy = DCID (12110) Lower Extremity Veins                  |           |
| 14 | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-49403, SRT, "Vein of Lower Extremity")   | LE Vein   |
|    |          |         |   | \$SectionLaterality = EV (G-A100, SRT, "Right")                 |           |
|    |          |         |   | \$Anatomy = DCID (12110) Lower Extremity Veins                  |           |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-49403, SRT, "Vein of Lower Extremity")   | LE Vein   |
|    |          |         |   | \$Anatomy = DCID (12110) Lower Extremity Veins                  |           |
| 15 |          |         |   | \$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") | UE Artery |
|    |          |         |   | \$SectionLaterality = EV (G-A101, SRT, "Left")                  |           |
|    |          |         |   | \$Anatomy = DCID (12107) Upper Extremity Arteries               |           |
| 16 | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") | UE Artery |
|    |          |         |   | \$SectionLaterality = EV (G-A100, SRT, "Right")                 |           |
|    |          |         |   | \$Anatomy = DCID (12107) Upper Extremity Arteries               |           |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") | UE Artery |
|    |          |         |   | \$Anatomy = DCID (12107) Upper Extremity Arteries               |           |
| 17 |          |         |   | \$SectionScope = DT (T-49103, SRT, "Vein Of Upper Extremity")   | UE Vein   |

|    |          |         |   |  |                             |
|----|----------|---------|---|--|-----------------------------|
|    |          |         |   | \$SectionLaterality = EV (G-A101, SRT, "Left")                     |                             |
|    |          |         |   | \$Anatomy = DCID (12108) Upper Extremity Veins                     |                             |
| 18 | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-49103, SRT, "Vein Of Upper Extremity")      | UE Vein                     |
|    |          |         |   | \$SectionLaterality = EV (G-A100, SRT, "Right")                    |                             |
|    |          |         |   | \$Anatomy = DCID (12108) Upper Extremity Veins                     |                             |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-49103, SRT, "Vein Of Upper Extremity")      | UE Vein                     |
|    |          |         |   | \$Anatomy = DCID (12108) Upper Extremity Veins                     |                             |
| 19 |          |         |   | \$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") | Abdomen Renal               |
|    |          |         |   | \$SectionLaterality = EV (G-A100, SRT, "Right")                    |                             |
|    |          |         |   | \$Anatomy = DCID (12115) Renal Vessels                             |                             |
| 20 | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") | Abdomen Renal               |
|    |          |         |   | \$SectionLaterality = EV (G-A101, SRT, "Left")                     |                             |
|    |          |         |   | \$Anatomy = DCID (12115) Renal Vessels                             |                             |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular Ultrasound Section | \$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") | Abdomen Renal               |
|    |          |         |   | \$Anatomy = DCID (12115) Renal Vessels                             |                             |
| 21 |          |         |   | \$SectionScope = DT (T-46002, SRT, "Artery of Abdomen")            | Abdomen Artery (Unilateral) |
|    |          |         |   | \$SectionLaterality = EV (G-A103, SRT, "Unilateral")               |                             |
|    |          |         |   | \$Anatomy = DCID (12112) Abdominal Arteries (unilateral)           |                             |

|    |          |         |  |  |                              |
|----|----------|---------|--|--|------------------------------|
| 22 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-46002,<br>SRT, "Artery of Abdomen") | Abdomen Artery               |
|    |          |         |  | \$SectionLaterality = EV (G-A100,<br>SRT, "Right")         |                              |
|    |          |         |  | \$Anatomy = DCID (12111)<br>Abdominal Arteries (lateral)   |                              |
| 23 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-46002,<br>SRT, "Artery of Abdomen") | Abdomen Artery               |
|    |          |         |  | \$SectionLaterality = EV (G-A101,<br>SRT, "Left")          |                              |
|    |          |         |  | \$Anatomy = DCID (12111)<br>Abdominal Arteries (lateral)   |                              |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-46002,<br>SRT, "Artery of Abdomen") | Abdomen Artery               |
|    |          |         |  | \$Anatomy = DCID (12111)<br>Abdominal Arteries (lateral)   |                              |
| 25 |          |         |  | \$SectionScope = DT (T-487A0,<br>SRT, "Vein of Abdomen")   | Abdomen Vein<br>(Unilateral) |
|    |          |         |  | \$SectionLaterality = EV (G-A103,<br>SRT, "Unilateral")    |                              |
|    |          |         |  | \$Anatomy = DCID (12114)<br>Abdominal Veins (unilateral)   |                              |
| 26 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-487A0,<br>SRT, "Vein of Abdomen")   | Abdomen Vein                 |
|    |          |         |  | \$SectionLaterality = EV (G-A100,<br>SRT, "Right")         |                              |
|    |          |         |  | \$Anatomy = DCID (12113)<br>Abdominal Veins (lateral)      |                              |
|    | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-487A0,<br>SRT, "Vein of Abdomen")   | Abdomen Vein                 |
|    |          |         |  | \$SectionLaterality = EV (G-A101,<br>SRT, "Left")          |                              |
|    |          |         |  | \$Anatomy = DCID (12113)<br>Abdominal Veins (lateral)      |                              |
| 28 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-487A0,<br>SRT, "Vein of Abdomen")   | Abdomen Vein                 |

|    |          |         |  |   |                      |
|----|----------|---------|--|---|----------------------|
|    |          |         |  | \$Anatomy = DCID (12113)<br>Abdominal Veins (lateral)   |                      |
| 29 | CONTAINS | INCLUDE | DTID (5103) Vascular<br>Ultrasound Section | \$SectionScope = DT (T-D4000,<br>SRT,"Abdomen")<br>\$Anatomy = DCID (6204) Anatomic<br>Non-Colon Findings | Abdomen<br>(B- MODE) |

### 9.2.1.1 Observation ConText (TID 1001)

**Table 9.2-2  
OBSERVATION CONTEXT IN VASCULAR SR**

|     | REL             | VT    | Concept Name                             | Unit / CODE Value  | Label                    |
|-----|-----------------|-------|--|--|--------------------------|
| B-1 | HAS OBS CONTEXT | CODE  | (121005, DCM, "Observer<br>Type")        | (121006, DCM,<br>"Person")   |                          |
| B-2 | HAS OBS CONTEXT | PNAME | (121008, DCM, "Person<br>Observer Name") |  | Ref. Physician           |
| B-3 | HAS OBS CONTEXT | CODE  | (121024, DCM, "Subject Class")           | (121025 ,DCM,"Patient")  |                          |
| B-4 | HAS OBS CONTEXT | PNAME | (121029, DCM, "Subject Name")            |  | Last Name,<br>First Name |
| B-5 | HAS OBS CONTEXT | DATE  | (121031, DCM, "Subject Birth<br>Date")   | DCID (7456) Units of<br>Measure for Age                              |                          |
| B-6 | HAS OBS CONTEXT | CODE  | (121032, DCM, "Subject Sex")             | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown<br>sex") | Gender                   |
| B-7 | HAS OBS CONTEXT | NUM   | (121033, DCM, "Subject Age")             | (mo, UCUM, "month")  | Not Used                 |

### 9.2.1.2 Vascular Patient Characteristics (TID 5101)

**Table 9.2-3  
VASCULAR PATIENT CHARACTERISTICS IN VASCULAR SR**

|     | REL      | VT        | Concept Name         | Unit / CODE Value | Label |
|-----|----------|-----------|----------------------|-------------------|-------|
| B-8 | CONTAINS | CONTAINER | (121118,DCM "Patient |                   |       |

|         |               |      |   |  |  |
|---------|---------------|------|---|--|--|
|         |               |      | Characteristics")                         |  |  |
| B-8-1   | CONTAINS      | NUM  | (121033, DCM, "Subject Age")              | Units = DCID (7456) Units of Measure for Age             |  |
| B-8-2   | CONTAINS      | CODE | (121032, DCM, "Subject Sex")              | DCID (7455) Sex  |  |
| B-8-3   | CONTAINS      | NUM  | (82776, LN, " Body Surface Area)          |  |  |
| B-8-3-1 | INFERRED FROM | CODE | (8278-4, LN, "Body Surface Area Formula") | (122241, DCM, "BSA = $0.007184*WT^{0.425}*HT^{0.725}$ ") |  |

### 9.2.1.3 Vascular Procedure Summary Section (TID 5102)

Table 9.2-4

#### VASCULAR PROCEDURE SUMMARY SECTION IN VASCULAR SR

|       | REL      | VT        | Concept Name             | Unit / CODE Value | Label |
|-------|----------|-----------|--------------------------|-------------------|-------|
| B-9   | CONTAINS | CONTAINER | (121111,DCM "Summary")   |                   |       |
| B-9-1 | CONTAINS | TEXT      | (121106, DCM, "Comment") |                   |       |

### 9.2.1.4 Vascular Ultrasound Section (TID 5103)

Table 9.2-5

#### VASCULAR ULTRASOUND SECTION IN VASCULAR SR

|          | REL                   | VT        | Concept Name                               | Unit / CODE Value                     | Label |
|----------|-----------------------|-----------|--|---------------------------------------|-------|
| B-10     |                       | CONTAINER | DT (121070, DCM, "Findings")               |                                       |       |
| B-10-1   | HAS<br>CONCEPT<br>MOD | CODE      | EV (G-C0E3, SRT, "Finding Site")           | \$SectionScope                        |       |
| B-10-2   | HAS<br>CONCEPT<br>MOD | CODE      | EV (G-C171, SRT, "Laterality")             | \$SectionLaterality                   |       |
|          | CONTAINS              | INCLUDE   | DTID (5104) Vascular Measurement Group     | \$AnatomyGroup = \$Anatomy            |       |
| B-10-3   | CONTAINS              | CONTAINER | \$AnatomyGroup                             |                                       |       |
| B-10-3-1 | HAS<br>CONCEPT<br>MOD | CODE      | EV (G-A1F8, SRT, "Topographical Modifier") | DCID (12116) Vessel Segment Modifiers |       |

|            |                       |               |  |   |         |
|------------|-----------------------|---------------|--|---|---------|
| B-10-3-2   | HAS<br>CONCEPT<br>MOD | CODE          | EV (125101, DCM, "Vessel Branch")                        | DCID (12117) Vessel Branch Modifiers                      |         |
|            | CONTAINS              | INCLUDE       | DTID (300) Measurement                                   | \$Measurement = DCID (12119) Vascular Ultrasound Property |         |
|            |                       |               |  | \$Derivation = DCID (3627) Measurement Type               |         |
| B-10-3-3   |                       | NUM           | \$Measurement  | Units = \$Units   |         |
| B-10-3-3-1 | HAS<br>CONCEPT<br>MOD | CODE          | EV (121401, DCM, "Derivation")                           | \$Derivation  |         |
|            | CONTAINS              | INCLUDE       | DTID (300) Measurement                                   | \$Measurement = \$AnatomyRatio                            |         |
| B-10-3-4   | NUM                   | \$Measurement | Units = \$Units  | B-10-3-4  |         |
|            | CONTAINS              | INCLUDE       | DTID (SM99110) Vascular User Creation Group Section      | CONTAINS  | INCLUDE |
| B-10-4     | CONTAINS              | CONTAINER     | EV (99900-Creation ID, MDSN, "User Creation Group Name") | *Creation ID:<br>Randomly generated 7-digit unique ID     |         |
| B-10-4-1   | HAS<br>CONCEPT<br>MOD | CODE          | EV (G-A1F8, SRT, "Topographical Modifier")               | DCID (12116) Vessel Segment Modifiers                     |         |
| B-10-4-2   | CONTAINS              | NUM           | EV (99900-Creation ID, MDSN, "User Creation Item Name")  |   |         |
| B-10-4-2-1 | HAS<br>CONCEPT<br>MOD | CODE          | EV (121401, DCM, "Derivation")                           | \$Derivation  |         |

## 9.2.2 Vascular Measurement and Calculation used in Vascular SR

**Table 9.2-6  
Carotid**

| Label | DICOM SR Concept Name (CDS<br>CV CM) | Laterality | Topographical Modifier | Item<br>Configuration |
|-------|--------------------------------------|------------|------------------------|-----------------------|
|-------|--------------------------------------|------------|------------------------|-----------------------|

|               |   |                                   |  |  |
|---------------|---|-----------------------------------|--|--|
| Subclavian A  | SRTT-46100\Subclavian Artery  | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| CCA           | SRTT-45100\Common Carotid Artery  | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 2<br>Ref. Table 9.2<br>17 |
| Bulb          | SRTT-45170\Carotid Bulb   | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 2<br>Ref. Table 9.2<br>17 |
| ICA           | SRTT-45300\Internal Carotid Artery  | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 2<br>Ref. Table 9.2<br>17 |
| ECA           | SRTT-45200\External Carotid Artery  | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 2<br>Ref. Table 9.2<br>17 |
| Vertebral A   | SRTT-45700\Vertebral Artery   | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| Vertebral 0~4 | (99201-11, MDSN, "Vertebral0")<br>(99201-10, MDSN, "Vertebral1")<br>(99201-9, MDSN, "Vertebral2")<br>(99201-8, MDSN, "Vertebral3")<br>(99201-7, MDSN, "Vertebral4") | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid-<br>longitudinal<br>SRTG-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |

**Table 9.2-7  
LE Artery**

| Label | DICOM SR Concept Name (CDS<br>CV CM) | Laterality                        | Topographical Modifier               | Item<br>Configuration        |
|-------|--------------------------------------|-----------------------------------|--------------------------------------|------------------------------|
| CIA   | SRTT-46710\Common Iliac Artery       | SRTG-A100\Right<br>SRTG-A101\Left | SRTG-A118\Proximal<br>SRTG-A188\Mid- | Calculation<br>Items Table 1 |

|             |  |                                     |   |                      |
|-------------|--|-------------------------------------|---|----------------------|
|             |  |                                     | longitudinal<br>SRT\G-A119\Distal   | Ref. Table 9.2<br>16 |
| IIA         | SRT\T-46740\Internal Iliac Artery      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| EIA         | SRT\T-46910\External Iliac Artery      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| CFA         | SRT\T-47400\Common Femoral Artery      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| SFA         | SRT\T-47403\Superficial Femoral Artery | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| DFA         | SRT\T-47440\Profunda Femoris Artery    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| Popliteal A | SRT\T-47500\Popliteal Artery           | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| ATA         | SRT\T-47700\Anterior Tibial Artery     | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| PTA         | SRT\T-47600\Posterior Tibial Artery    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |                      |
| Peroneal A  | SRT\T-47630\Peroneal Artery            | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-                                      |                      |

|              |                                   |                                     |   |
|--------------|-----------------------------------|-------------------------------------|---|
|              |                                   |                                     | longitudinal<br>SRT\G-A119\Distal   |
| DPA          | SRT\T-47741\Dorsalis Pedis Artery | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| MPA          | SRT\T-47690\Plantar Arterial Arch | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| LPA          | SRT\T-47690\Plantar Arterial Arch | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Metatarsal A | MDSN\99201-1\Metatarsal Artery    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Digital A    | MDSN\99201-2\Digitgal Artery      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-8**  
**LE Vein**

| Label | DICOM SR Concept Name (CDS<br>CV CM) | Laterality                          | Topographical Modifier  | Item<br>Configuration                                |
|-------|--------------------------------------|-------------------------------------|---|--|
| CIV   | SRT\T-48920\Common Iliac Vein        | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 3<br>Ref. Table 9.2<br>18 |
| IIV   | SRT\T-48940\Internal iliac vein      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| EIV   | SRT\T-48930\External Iliac Vein      | SRT\G-A100\Right                    | SRT\G-A118\Proximal   |  |

|             |                                     |                                     |   |
|-------------|-------------------------------------|-------------------------------------|---|
|             |                                     | SRT\G-A101\Left                     | SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal                        |
| CFV         | SRT\G-035B\Common Femoral Vein      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| PFV         | SRT\T-49660\Profunda Femoris Vein   | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| SFV         | SRT\G-035A\Superficial Femoral Vein | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| GSV         | SRT\T-49530\Great Saphenous Vein    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Popliteal V | SRT\T-49640\Popliteal Vein          | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| LSV         | SRT\T-49550\Lesser Saphenous Vein   | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| ATV         | SRT\T-49630\Anterior Tibial Vein    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| PTV         | SRT\T-49620\Posterior Tibial Vein   | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Peroneal V  | SRT\T-49650\Peroneal Vein           | SRT\G-A100\Right                    | SRT\G-A118\Proximal   |

|              |                                    |                                     |   |
|--------------|------------------------------------|-------------------------------------|---|
|              |                                    | SRT\G-A101\Left                     | SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal                        |
| MPV          | MDSN\99203-01\Medial Plantar Vein  | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| LPV          | MDSN\99203-02\Lateral Plantar Vein | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Metatarsal V | MDSN\99203-03\Metatarsal Vein      | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Digital V    | MDSN\99203-04\Digital Vein         | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-9  
UE Artery**

| Label        | DICOM SR Concept Name (CDS CV CM) | Laterality                          | Topographical Modifier  | Item Configuration                                   |
|--------------|-----------------------------------|-------------------------------------|---|--|
| Subclavian A | SRT\T-46100\Subclavian Artery     | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| Axillary A   | SRT\T-47100\Axillary Artery       | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |
| Brachial A   | SRT\T-47160\Brachial Artery       | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |

|          |                                     |                                     |   |
|----------|-------------------------------------|-------------------------------------|---|
| Radial A | SRT\T-47300\Radial Artery           | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Ulnar A  | SRT\T-47200\Ulnar Artery            | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| SPA      | SRT\T-47240\Superficial Palmar Arch | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-10**  
**UE Vein**

| Label              | DICOM SR Concept Name (CDS CV CM) | Laterality                          | Topographical Modifier  | Item Configuration                                   |
|--------------------|-----------------------------------|-------------------------------------|---|--|
| Internal Jugular V | SRT\T-48170\Internal Jugular vein | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 3<br>Ref. Table 9.2<br>18 |
| Innominate V       | SRT\T-48620\Innominate vein       | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| Subclavian V       | SRT\T-48330\Subclavian vein       | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| Axillary V         | SRT\T-49110\Axillary vein         | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| Brachial V         | SRT\T-49350\Brachial vein         | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal                      |  |
|                    |                                   |                                     |   |  |

|            |                           |                                     |   |
|------------|---------------------------|-------------------------------------|---|
|            |                           |                                     | SRT\G-A119\Distal   |
| Cephalic V | SRT\T-49240\Cephalic vein | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Basilic V  | SRT\T-48052\Basilic vein  | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Radial V   | SRT\T-49340\Radial vein   | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| Ulnar V    | SRT\T-49330\Ulnar vein    | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-11**  
**TCD**

| Label     | DICOM SR Concept Name (CDS CV CM)                   | Laterality                          | Topographical Modifier  | Item Configuration                                   |
|-----------|---|-------------------------------------|---|--|
| Basilar A | SRT\T-45800\Basilar Artery                          | SRT\ G-<br>A103\Unilateral          | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| ACA       | SRT\T-45540\Anterior Cerebral<br>Artery             | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| MCA       | SRT\T-45600\Middle Cerebral<br>Artery               | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| PCA (P1)  | SRT\R-10253\Posterior Cerebral<br>Artery P1 Segment | SRT\G-A100\Right<br>SRT\G-A101\Left |   |  |

|               |   |                                     |   |
|---------------|---|-------------------------------------|---|
| PCA (P2)      | SRT\R-10255\Posterior Cerebral Artery P2 Segment  | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| ACommA        | SRT\T-45530\Anterior communicating artery         | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| PCommA        | SRT\T-45520\Posterior communicating artery        | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| ExtICA        | MDSN\99201-6\Exterior intracranial carotid artery | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| TICA          | SRT\R-102BD\Terminal internal carotid artery      | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| Siphon        | SRT\T-45308\Carotid Siphon                        | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| Ophthalmic A. | SRT\T-45400\Ophthalmic Artery                     | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| Vertebral A.  | SRT\T-45700\Vertebral artery                      | SRT\G-A100\Right<br>SRT\G-A101\Left |   |
| Vertebral 4   | MDSN\99201-7\Vertebral 4                          | SRT\G-A100\Right<br>SRT\G-A101\Left |   |

**Table 9.2-12  
Abdomen Artery**

| Label     | DICOM SR Concept Name (CDS CV CM) | Laterality                 | Topographical Modifier  | Item Configuration                                   |
|-----------|-----------------------------------|----------------------------|---|--|
| Aorta     | SRT\T-42000\Aorta                 | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 4<br>Ref. Table 9.2<br>19 |
| Celiac A  | SRT\T-46400\Celiac Axis           | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| Splenic A | SRT\T-46460\Splenic Artery        | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal                      |  |

|                    |  |                                     |   |
|--------------------|--|-------------------------------------|---|
|                    |  |                                     | SRT\G-A119\Distal   |
| SMA                | SRT\T-46510\Superior Mesenteric Artery       | SRT\ G-A103\Unilateral              | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| IMA                | SRT\T-46520\Inferior Mesenteric Artery       | SRT\ G-A103\Unilateral              | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| C Hepatic A        | SRT\T-46421\Common Hepatic Artery            | SRT\ G-A103\Unilateral              | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| R Hepatic A        | SRT\T-46423\Right Branch of Hepatic Artery   | SRT\ G-A103\Unilateral              | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| L Hepatic A        | SRT\T-46427\Left Branch of Hepatic Artery    | SRT\ G-A103\Unilateral              | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Renal Aortic Ratio | LN\33869-9\Renal Artery/Aorta velocity ratio | SRT\ G-A103\Unilateral              |   |
| Common Iliac A     | SRT\T-46710\Common Iliac Artery              | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| External Iliac A   | SRT\T-46910\External Iliac Artery            | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Internal Iliac A   | SRT\T-46740\Internal Iliac Artery            | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-13  
Abdomen Vein**

| Label       | DICOM SR Concept Name (CDS CV CM)               | Laterality                 | Topographical Modifier  | Item Configuration                                   |
|-------------|---|----------------------------|---|--|
| M Hepatic V | SRT\T-48726\Middle Hepatic Vein                 | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 3<br>Ref. Table 9.2<br>18 |
| R Hepatic V | SRT\T-48725\Right Hepatic Vein                  | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| L Hepatic V | SRT\T-48727\Left Hepatic Vein                   | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| IVC         | SRT\T-48710\Inferior Vena Cava                  | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| M Portal V  | SRT\T-48810\Portal Vein                         | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| R Portal V  | SRT\T-48813\Right Main Branch of<br>Portal Vein | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| L Portal V  | SRT\T-48814\Left Main Branch of<br>Portal Vein  | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |
| SMV         | SRT\T-48840\Superior Mesenteric<br>Vein         | SRT\ G-<br>A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |  |

|           |                                      |                        |   |
|-----------|--------------------------------------|------------------------|---|
| IMV       | SRT\T-48910\Inferior Mesenteric Vein | SRT\ G-A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |
| Splenic V | SRT\T-48890\Splenic Vein             | SRT\ G-A103\Unilateral | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |

**Table 9.2-14  
Abdomen Renal**

| Label           | DICOM SR Concept Name (CDS CV CM)        | Laterality                          | Topographical Modifier  | Item Configuration                                   |
|-----------------|--|-------------------------------------|---|--|
| Renal A         | SRT\T-46600\Renal Artery                 | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal | Calculation<br>Items Table 1<br>Ref. Table 9.2<br>16 |
| Upper Renal A   | SRT\T-46600\Renal Artery                 | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |
| Lower Renal A   | SRT\T-46600\Renal Artery                 | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |
| Arcuate Renal A | SRT\T-4668A\Arcuate Artery of the Kidney | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |
| Lobular Renal A | SRT\T-4667D\Interlobar Artery of Kidney  | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal<br>SRT\G-A119\Distal |  |
| Arcuate A       | MDSN\99201-3\Arcuate Artery              | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-longitudinal                      |  |

|             |                             |                                     |   |
|-------------|-----------------------------|-------------------------------------|---|
|             |                             |                                     | SRT\G-A119\Distal   |
| Segmental A | SRTT-46659\Segmental Artery | SRT\G-A100\Right<br>SRT\G-A101\Left | SRT\G-A118\Proximal<br>SRT\G-A188\Mid-<br>longitudinal<br>SRT\G-A119\Distal |
| RenalVein   | SRTT-48740\Renal Vein       | SRT\G-A100\Right<br>SRT\G-A101\Left |   |

**Table 9.2-15  
Abdomen (2D)**

| Label         | DICOM SR Concept Name (CDS CV CM) | Laterality                          | Topographical Modifier |
|---------------|-----------------------------------|-------------------------------------|------------------------|
| Liver         | (T-62000, SRT, "Liver")           | SRT\ G-A103\Unilateral              |                        |
| Spleen        | (T-C3000, SRT, "Spleen")          | SRT\ G-A103\Unilateral              |                        |
| Gall bladder  | (T-63000, SRT, "Gall bladder")    | SRT\ G-A103\Unilateral              |                        |
| Pancreas      | (T-65000, SRT, "Pancreas")        | SRT\ G-A103\Unilateral              |                        |
| Pancreas Head | (99016-2, MDSN, "Pancreas Head")  |                                     |                        |
| Pancreas Body | (99016-3, MDSN, "Pancreas Body")  |                                     |                        |
| Pancreas Tail | (99016-4, MDSN, "Pancreas Tail")  |                                     |                        |
| Bowel         | (99016-6, MDSN, "Bowel")          |                                     |                        |
| Kidney        | (T-71000, SRT, "Kidney")          | SRT\G-A100\Right<br>SRT\G-A101\Left |                        |
| Kidney Vol.   | 99016-19, MDSN, "Kidney Volume")  |                                     |                        |
| Kidney L      | (99016-20 MDSN, "Kidney Length")  |                                     |                        |
| Kidney W      | (99016-21 MDSN, "Kidney Width")   |                                     |                        |
| Kidney H      | (99016-22 MDSN, "Kidney Height")  |                                     |                        |

**Table 9.2-116  
Calculation Items Table 1**

| Label | DICOM SR Concept Name (CDS CV CM)     |
|-------|---------------------------------------|
| PS    | LN\11726-7\Peak Systolic Velocity     |
| ED    | LN\11653-3\End Diastolic Velocity     |
| MD    | LN\11665-7\Minimum Diastolic Velocity |

|              |  |
|--------------|--|
| TAmax        | LN\11692-1\Time averaged peak velocity             |
| TAmean       | LN\20352-1\Time averaged mean velocity             |
| PGmax        | LN\20247-3\Peak Gradient                           |
| PGmean       | LN\20256-4\Mean Gradient                           |
| PS/ED        | LN\12144-2\Systolic to Diastolic Velocity Ratio    |
| ED/PS        | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |
| RI           | LN\12023-8\Resistivity Index                       |
| PI           | LN\12008-9\Pulsatility Index                       |
| AccT         | LN\20168-1\Acceleration Time                       |
| Acc          | MDSN\99200-06\Acceleration                         |
| DecT         | LN\20217-6\Deceleration Time                       |
| Dec          | MDSN\99200-07\Deceleration                         |
| Aout         | SRT\G-0366\Vessel lumen cross-sectional area       |
| Ain          | SRT\R-1025D\Vessel Intimal Cross-Sectional Area    |
| %StA         | SRT\R-101BA\Lumen Area Stenosis                    |
| Dout         | SRT\G-0364\Vessel lumen diameter                   |
| Din          | SRT\R-1025C\Vessel Intimal Diameter                |
| %StD         | SRT\R-101BB\Lumen Diameter Stenosis                |
| Vesl. Area   | MDSN\99200-02\Vessel Area                          |
| Vol. Flow(A) | LN\33878-0\Volume flow                             |
| Vesl. Dist   | SRT\G-0365\Vessel outside diameter                 |
| Vol. Flow(D) | LN\33878-0\Volume flow                             |

**Table 9.2-17**  
**Calculation Items Table 2**

| <b>Label</b> | <b>DICOM SR Concept Name (CDS CV CM)</b>           |
|--------------|--|
| PS           | LN\11726-7\Peak Systolic Velocity                  |
| ED           | LN\11653-3\End Diastolic Velocity                  |
| MD           | LN\11665-7\Minimum Diastolic Velocity              |
| TAmax        | LN\11692-1\Time averaged peak velocity             |
| TAmean       | LN\20352-1\Time averaged mean velocity             |
| PGmax        | LN\20247-3\Peak Gradient                           |
| PGmean       | LN\20256-4\Mean Gradient                           |
| PS/ED        | LN\12144-2\Systolic to Diastolic Velocity Ratio    |
| ED/PS        | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |

|              |   |
|--------------|---|
| RI           | LN\12023-8\Resistivity Index                    |
| PI           | LN\12008-9\Pulsatility Index                    |
| AccT         | LN\20168-1\Acceleration Time                    |
| Acc          | MDSN\99200-06\Acceleration                      |
| DecT         | LN\20217-6\Deceleration Time                    |
| Dec          | MDSN\99200-07\Deceleration                      |
| Aout         | SRT\G-0366\Vessel lumen cross-sectional area    |
| Ain          | SRT\R-1025D\Vessel Intimal Cross-Sectional Area |
| %StA         | SRT\R-101BA\Lumen Area Stenosis                 |
| Dout         | SRT\G-0364\Vessel lumen diameter                |
| Din          | SRT\R-1025C\Vessel Intimal Diameter             |
| %StD         | SRT\R-101BB\Lumen Diameter Stenosis             |
| Vesl. Area   | MDSN\99200-02\Vessel Area                       |
| Vol. Flow(A) | LN\33878-0\Volume flow                          |
| Vesl. Dist   | SRT\G-0365\Vessel outside diameter              |
| Vol. Flow(D) | LN\33878-0\Volume flow                          |
| IMT          | MDSN\99200-05\Intima-media thickness            |

**Table 9.2-18**  
**Calculation Items Table 3**

| Label      | DICOM SR Concept Name (CDS CV CM)  |
|------------|------------------------------------|
| Vmax       | MDSN\99200-03\Max Velocity         |
| Dur T      | MDSN\99200-04\Duration Time        |
| Vesl. Dist | SRT\G-0365\Vessel outside diameter |
| Vesl. Area | MDSN\99200-02\Vessel Area          |

**Table 9.2-19**  
**Calculation Items Table 4**

| Label       | DICOM SR Concept Name (CDS CV CM) |
|-------------|-----------------------------------|
| SAG AP      | SRT\G-A22A\Length                 |
| Trans Width | SRT\G-A220\Width                  |
| TRANS AP    | DCM\121207\Height                 |
| Volume      | SRT\G-D705\Volume                 |
| PS          | LN\11726-7\Peak Systolic Velocity |
| ED          | LN\11653-3\End Diastolic Velocity |

|                    |  |
|--------------------|--|
| MD                 | LN\11665-7\Minimum Diastolic Velocity              |
| TAm <sub>ax</sub>  | LN\11692-1\Time averaged peak velocity             |
| TAm <sub>ean</sub> | LN\20352-1\Time averaged mean velocity             |
| PG <sub>max</sub>  | LN\20247-3\Peak Gradient                           |
| PG <sub>mean</sub> | LN\20256-4\Mean Gradient                           |
| PS/ED              | LN\12144-2\Systolic to Diastolic Velocity Ratio    |
| ED/PS              | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |
| RI                 | LN\12023-8\Resistivity Index                       |
| PI                 | LN\12008-9\Pulsatility Index                       |
| AccT               | LN\20168-1\Acceleration Time                       |
| Acc                | MDSN\99200-06\Acceleration                         |
| DecT               | LN\20217-6\Deceleration Time                       |
| Dec                | MDSN\99200-07\Deceleration                         |
| A <sub>out</sub>   | SRT\G-0366\Vessel lumen cross-sectional area       |
| A <sub>in</sub>    | SRT\R-1025D\Vessel Intimal Cross-Sectional Area    |
| %StA               | SRT\R-101BA\Lumen Area Stenosis                    |
| D <sub>out</sub>   | SRT\G-0364\Vessel lumen diameter                   |
| D <sub>in</sub>    | SRT\R-1025C\Vessel Intimal Diameter                |
| %StD               | SRT\R-101BB\Lumen Diameter Stenosis                |
| Vesl. Area         | MDSN\99200-02\Vessel Area                          |
| Vol. Flow(A)       | LN\33878-0\Volume flow                             |
| Vesl. Dist         | SRT\G-0365\Vessel outside diameter                 |
| Vol. Flow(D)       | LN\33878-0\Volume flow                             |

### 9.3 ADULT ECHOCARDIOGRAPHY PROCEDURE REPORT TEMPLATES

#### 9.3.1 Adult Echocardiography Procedure Report (TID 5200)

Table 9.3-1

Adult Echocardiography Ultrasound Procedure Report Tempalte

| No | Rel With Parent | VT        | Concept Name  | Comments |
|----|-----------------|-----------|---|----------|
| 1  |                 | CONTAINER | EV (125200, DCM, "Adult Echocardiography Procedure Report") |          |

|    |                       |           |  |   |
|----|-----------------------|-----------|--|---|
| 2  | HAS<br>CONCEPT<br>MOD | INCLUDE   | DTID (1204) Language of<br>Content Item and<br>Descendants |   |
| 3  | HAS OBS<br>CONTEXT    | INCLUDE   | DTID (1001) Observation<br>Context                         |   |
| 4  | CONTAINS              | INCLUDE   | DTID (5201)<br>Echocardiography Patient<br>Characteristics |   |
| 5  | CONTAINS              | CONTAINER | (111028, DCM, "Image<br>Library")                          |   |
| 6  | CONTAINS              | IMAGE     | No purpose of reference                                    |   |
| 7  | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-32600, SRT, "Left Ventricle")        |
|    |                       |           |  | \$MeasType = DCID (12200) Echocardiography Left<br>Ventricle  |
| 8  | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-32500, SRT, "Right Ventricle")       |
|    |                       |           |  | \$MeasType = DCID (12204) Echocardiography Right<br>Ventricle |
| 9  | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-32300, SRT, "Left Atrium")           |
|    |                       |           |  | \$MeasType = DCID (12205) Echocardiography Left<br>Atrium     |
| 10 | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-32200, SRT, "Right Atrium")          |
|    |                       |           |  | \$MeasType = DCID (12206) Echocardiography Right<br>Atrium    |
| 11 | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-35400, SRT, "Aortic Valve")          |
|    |                       |           |  | \$MeasType = DCID (12211) Echocardiography Aortic<br>Valve    |
| 12 | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-35300, SRT, "Mitral Valve")          |
|    |                       |           |  | \$MeasType = DCID (12207) Echocardiography Mitral<br>Valve    |
| 13 | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-35200, SRT, "Pulmonic Valve")        |
|    |                       |           |  | \$MeasType = DCID (12209) Echocardiography<br>Pulmonic Valve  |
| 14 | CONTAINS              | INCLUDE   | DTID (5202) Echo Section                                   | \$SectionSubject = EV (T-35100, SRT, "Tricuspid Valve")       |
|    |                       |           |  | \$MeasType = DCID (12208) Echocardiography<br>Tricuspid Valve |

|    |          |         |   |   |
|----|----------|---------|---|---|
| 15 | CONTAINS | INCLUDE | DTID (5202) Echo Section                              | \$SectionSubject = EV (T-42000, SRT, "Aorta")                     |
|    |          |         |   | \$MeasType= DCID (12212) Echocardiography Aorta                   |
| 16 | CONTAINS | INCLUDE | DTID (5202) Echo Section                              | \$SectionSubject = EV (T-44000, SRT, "Pulmonary artery")          |
|    |          |         |   | \$MeasType DCID (12210) = Echocardiography Pulmonary Artery       |
| 17 | CONTAINS | INCLUDE | DTID (5202) Echo Section                              | \$SectionSubject = EV (T-48600, SRT, "Vena Cava"                  |
|    |          |         |   | \$MeasType = DCID (12215) Echocardiography Vena Cavae             |
| 18 | CONTAINS | INCLUDE | DTID (5202) Echo Section                              | \$SectionSubject = EV (T-48581, SRT, "Pulmonary Venous Structure" |
|    |          |         |   | \$MeasType = DCID (12214) Echocardiography Pulmonary Veins        |
| 19 | CONTAINS | INCLUDE | DTID (5202) Echo Section                              | \$SectionSubject = EV (P5-30031, SRT, "Cardiac Shunt Study")      |
|    |          |         |   | \$MeasType = DCID (12217) Echocardiography Cardiac Shunt          |
| 20 | CONTAINS | INCLUDE | DTID (SM99210) Adult Echo User Creation Group Section | Ref. Section 9.3.1.1  |

### 9.3.1.1 Adult Echo User Creation Group Section (TID SM99210)

Table 9.3-2

#### User Creation Group Section in Adult Echo SR

|        | REL      | VT        | Concept Name  | Unit / CODE Value | Label | Comments   |
|--------|----------|-----------|---|-------------------|-------|--|
| C-20   | CONTAINS | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |                   |       | *Creation ID: Randomly generated 7-digit unique ID |
| C-20-1 | CONTAINS | NUM       | (99900-Creation ID, MDSN, "User Creation Item Name")  |                   |       |  |

|          |                    |      |                                |                           |  |  |
|----------|--------------------|------|--------------------------------|---------------------------|--|--|
| C-20-1-1 | HAS CONCEPT<br>MOD | CODE | (121401, DCM,<br>"Derivation") | Common CID-<br>Derivation |  |  |
|----------|--------------------|------|--------------------------------|---------------------------|--|--|

### 9.3.2 Cardiac Measurement and Calculation used in Adult Echocardiography SR

- Label – Label of measurement or calculation used in Cardiac Calc. package for the Ultrasound System
- FSite – Finding Site
- Concept – (CV, CSD, "Concept Name")
- Modifier – Additional codes and Modifiers used

**Table 9.3-3  
Cardiac Measurement and Calculation**

| Group   | Label            | Concept  |                                  | Modifiers                    |
|---------|------------------|--|----------------------------------|------------------------------|
| LV (2D) | IVSd             | (18154-5, LN, "Interventricular Septum Diastolic Thickness")               | (T-32600, SRT, "Left Ventricle") | ImageMode=SRT\G-03A2\2D mode |
|         | LVIDd            | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension")           |                                  |                              |
|         | LVPWd            | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness")         |                                  |                              |
|         | IVSs             | (18158-6, LN, "Interventricular Septum Systolic Thickness")                |                                  |                              |
|         | LVIDs            | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")                |                                  |                              |
|         | LVPWs            | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness")          |                                  |                              |
|         | IVS% Thickening  | (18054-7, LN, "Interventricular Septum % Thickening")                      |                                  |                              |
|         | LVPW% Thickening | (18053-9, LN, "Left Ventricle Posterior Wall % Thickening")                |                                  |                              |
|         | IVSd/LVPWd       | (18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio") |                                  |                              |
|         | IVSs/LVPWs       | (18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio") |                                  |                              |

|        |            |  |                                  |                              |
|--------|------------|--|----------------------------------|------------------------------|
|        | RWT        | (99104-09, MDSN, "Relative Wall Thickness")                        |                                  |                              |
|        | LV % FS    | (18051-3, LN, "Left Ventricular Fractional Shortening")            |                                  |                              |
|        | LV EDV     | (18026-5, LN, "Left Ventricular End Diastolic Volume")             |                                  |                              |
|        | LV ESV     | (18148-7, LN, "Left Ventricular End Systolic Volume")              |                                  |                              |
|        | LV SV      | (F-32120, SRT, "Stroke Volume")                                    |                                  |                              |
|        | LV SI      | (F-00078, SRT, "Stroke Index")                                     |                                  |                              |
|        | LV CO      | (F-32100, SRT, "Cardiac Output")                                   |                                  |                              |
|        | LV CI      | (F-32110, SRT, "Cardiac Index")                                    |                                  |                              |
|        | LV EF      | (18043-0, LN, "Left Ventricular Ejection Fraction")                |                                  |                              |
|        | LV Mass(C) | (18087-7, LN, "Left Ventricle Mass")                               |                                  |                              |
|        | LV MI(C)   | (99104-01, MDSN, "Left Ventricular Mass Index")                    |                                  |                              |
| LV (M) | IVSd       | (18154-5, LN, "Interventricular Septum Diastolic Thickness")       | (T-32600, SRT, "Left Ventricle") | Image Mode=SRT\G-0394\M mode |
|        | LVIDd      | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension")   |                                  |                              |
|        | LVPWd      | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") |                                  |                              |
|        | LVIDs      | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")        |                                  |                              |
|        | IVSs       | (18158-6, LN, "Interventricular Septum Systolic Thickness")        |                                  |                              |
|        | LVIDs      | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")        |                                  |                              |
|        | LVPWs      | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness")  |                                  |                              |
|        | MAPSE      | (99104-10, MDSN, "Mitral Annular Plane Systolic Excursion")        |                                  |                              |

|                           |  |
|---------------------------|--|
| IVS%<br>Thickening<br>(M) | (18054-7, LN, "Interventricular Septum % Thickening")                      |
| LVPW%<br>Thickening (M)   | (18053-9, LN, "Left Ventricle Posterior Wall % Thickening")                |
| IVSd/LVPWd<br>(M)         | (18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio") |
| IVSs/LVPWs<br>(M)         | (18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio") |
| RWT                       | (99104-09, MDSN, "Relative Wall Thickness")                                |
| LV % FS (M)               | (18051-3, LN, "Left Ventricular Fractional Shortening")                    |
| LVEDV (M)                 | (18026-5, LN, "Left Ventricular End Diastolic Volume")                     |
| LV ESV (M)                | (18148-7, LN, "Left Ventricular End Systolic Volume")                      |
| LV SV (M)                 | (F-32120, SRT, "Stroke Volume")  |
| LV SI (M)                 | (F-00078, SRT, "Stroke Index")   |
| LV CO (M)                 | (F-32100, SRT, "Cardiac Output")   |
| LV CI (M)                 | (F-32110, SRT, "Cardiac Index")  |
| LV EF (M)                 | (18043-0, LN, "Left Ventricular Ejection Fraction")                        |
| LV Mass(C)                | (18087-7, LN, "Left Ventricle Mass")                                       |
| LV MI(C)                  | (99104-01, MDSN, "Left Ventricular Mass Index")                            |
| LVET                      | (122211, DCM, "Left Ventricular ejection time")                            |
| LVPEP                     | (79989-0, LN, "Left ventricular pre-ejection period")                      |
| LVPEP/ET (M)              | (59088-5, LN, "Pre-Ejection Period/Ejection Time Ratio")                   |
| Vcf                       | (99104-02, MDSN, "Velocity Circumferential Fiber Shortening")              |

|                      |                 |  |                                  |  |
|----------------------|-----------------|--|----------------------------------|--|
| LV Vol.<br>(Simpson) | LVLd A4C        | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | Image Mode=SRT\G-03A2\2D mode<br><br>(A2C)Image View=SRT\G-A19B\Apical two chamber<br><br>(A4C)Image View=SRT\G-A19C\Apical fourchamber<br><br>(Simpson) MeasurementMethod=DC M\125208\Method of Disks, SinglePlane<br><br>(Simpson BP) MeasurementMethod=DC M\125207\Method of Disks, Biplane |
|                      | LVAd A4C        | (G-0375, SRT, "Left Ventricular Diastolic Area")                 |                                  |  |
|                      | LVEDV A4C       | (18026-5, LN, "Left Ventricular End Diastolic Volume")           |                                  |  |
|                      | LVEDV Index A4C | (99104-15, MDSN, "Left Ventricular End Diastolic Volume Index")  |                                  |  |
|                      | LVLs A4C        | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")      |                                  |  |
|                      | LVAs A4C        | (G-0374, SRT, "Left Ventricular Systolic Area")                  |                                  |  |
|                      | LVESV A4C       | (18148-7, LN, "Left Ventricular End Systolic Volume")            |                                  |  |
|                      | LVESV Index A4C | (99104-16, MDSN, "Left Ventricular End Systolic Volume Index")   |                                  |  |
|                      | LV SV A4C       | (F-32120, SRT, "Stroke Volume")                                  |                                  |  |
|                      | LV CO A4C       | (F-32100, SRT, "Cardiac Output")                                 |                                  |  |
|                      | LV EF A4C       | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                                  |  |
|                      | LV SI A4C       | (F-00078, SRT, "Stroke Index")                                   |                                  |  |
|                      | LV CI A4C       | (F-32110, SRT, "Cardiac Index")                                  |                                  |  |
|                      | LVLd A2C        | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") |                                  |  |
|                      | LVAd A2C        | (G-0375, SRT, "Left Ventricular Diastolic Area")                 |                                  |  |
|                      | LVEDV A2C       | (18026-5, LN, "Left Ventricular End Diastolic Volume")           |                                  |  |
|                      | LVEDV Index A2C | (99104-15, MDSN, "Left Ventricular End Diastolic Volume Index")  |                                  |  |
|                      | LVLs A2C        | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")      |                                  |  |
|                      | LVAs A2C        | (G-0374, SRT, "Left Ventricular Systolic Area")                  |                                  |  |

|                 |                    |  |                                  |  |
|-----------------|--------------------|--|----------------------------------|--|
|                 | LVESV A2C          | (18148-7, LN, "Left Ventricular End Systolic Volume")            |                                  |  |
|                 | LVESV Index A2C    | (99104-16, MDSN, "Left Ventricular End Systolic Volume Index")   |                                  |  |
|                 | LV SV A2C          | (F-32120, SRT, "Stroke Volume")                                  |                                  |  |
|                 | LV CO A2C          | (F-32100, SRT, "Cardiac Output")                                 |                                  |  |
|                 | LV EF A2C          | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                                  |  |
|                 | LV SI A2C          | (F-00078, SRT, "Stroke Index")                                   |                                  |  |
|                 | LV CI A2C          | (F-32110, SRT, "Cardiac Index")                                  |                                  |  |
|                 | LVEDV BP           | (18026-5, LN, "Left Ventricular End Diastolic Volume")           |                                  |  |
|                 | LVEDV Index BP     | (99104-15, MDSN, "Left Ventricular End Diastolic Volume Index")  |                                  |  |
|                 | LVESV BP           | (18148-7, LN, "Left Ventricular End Systolic Volume")            |                                  |  |
|                 | LVESV Index BP     | (99104-16, MDSN, "Left Ventricular End Systolic Volume Index")   |                                  |  |
|                 | LV SV BP           | (F-32120, SRT, "Stroke Volume")                                  |                                  |  |
|                 | LV CO BP           | (F-32100, SRT, "Cardiac Output")                                 |                                  |  |
|                 | LV EF BP           | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                                  |  |
|                 | LV SI BP           | (F-00078, SRT, "Stroke Index")                                   |                                  |  |
|                 | LV CI BP           | (F-32110, SRT, "Cardiac Index")                                  |                                  |  |
| LV<br>Vol.(A/L) | LVLd A4C           | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | Image Mode=SRT\G-03A2\2D mode<br><br>(A2C)Image<br>View=SRT\G-A19B\Apical two chamber<br><br>(A4C)<br>Image View=SRT\G-A19C\Apical fourchamber |
|                 | LVAd A4C           | (G-0375, SRT, "Left Ventricular Diastolic Area")                 |                                  |  |
|                 | LVEDV A4C AL       | (18026-5, LN, "Left Ventricular End Diastolic Volume")           |                                  |  |
|                 | LVEDV Index A4C AL | (99104-15, MDSN, "Left Ventricular End Diastolic Volume Index")  |                                  |  |
|                 | LVLs A4C           | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")      |                                  |  |
|                 | LVAs A4C           | (G-0374, SRT, "Left Ventricular Systolic                         |                                  |  |

|                  |                    |  |                                  |  |
|------------------|--------------------|--|----------------------------------|--|
|                  |                    | Area")   |                                  | (A/L)Measurement<br>Method=DCM\125205\Area-Length Single Plane |
|                  | LVESV A4C AL       | (18148-7, LN, "Left Ventricular End Systolic Volume")            |                                  |  |
|                  | LVESV Index A4C AL | (99104-16, MDSN, "Left Ventricular End Systolic Volume Index")   |                                  |  |
|                  | LV SV A4C AL       | (F-32120, SRT, "Stroke Volume")                                  |                                  |  |
|                  | LV CO A4C AL       | (F-32100, SRT, "Cardiac Output")                                 |                                  |  |
|                  | LV EF A4C AL       | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                                  |  |
|                  | LV SI A4C AL       | (F-00078, SRT, "Stroke Index")                                   |                                  |  |
|                  | LV CI A4C AL       | (F-32110, SRT, "Cardiac Index")                                  |                                  |  |
|                  | LVLd A2C           | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") |                                  |  |
|                  | LVAAd A2C          | (G-0375, SRT, "Left Ventricular Diastolic Area")                 |                                  |  |
|                  | LVEDV A2C AL       | (18026-5, LN, "Left Ventricular End Diastolic Volume")           |                                  |  |
|                  | LVEDV Index A2C AL | (99104-15, MDSN, "Left Ventricular End Diastolic Volume Index")  |                                  |  |
|                  | LVLs A2C           | (29438-9, LN, "Left Ventricle Internal Systolic Dimension")      |                                  |  |
|                  | LVAAs A2C          | (G-0374, SRT, "Left Ventricular Systolic Area")                  |                                  |  |
|                  | LVESV A2C AL       | (18148-7, LN, "Left Ventricular End Systolic Volume")            |                                  |  |
|                  | LVESV Index A2C AL | (99104-16, MDSN, "Left Ventricular End Systolic Volume Index")   |                                  |  |
|                  | LV SV A2C AL       | (F-32120, SRT, "Stroke Volume")                                  |                                  |  |
|                  | LV CO A2C AL       | (F-32100, SRT, "Cardiac Output")                                 |                                  |  |
|                  | LV EF A2C AL       | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                                  |  |
|                  | LV SI A2C AL       | (F-00078, SRT, "Stroke Index")                                   |                                  |  |
|                  | LV CI A2C AL       | (F-32110, SRT, "Cardiac Index")                                  |                                  |  |
| LV Vol. (Bullet) | LVAAd SAX MV       | (G-0375, SRT, "Left Ventricular Diastolic Area")                 | (T-32600, SRT, "Left Ventricle") | Image Mode=SRT\G-03A2\2D mode                                  |

|         |                  |   |                                   |   |
|---------|------------------|---|-----------------------------------|---|
|         | LVAAs SAX MV     | (G-0374, SRT, "Left Ventricular Systolic Area")                               |                                   | (Bullet)  |
|         | LVLd Apical      | (18077-8, LN, "Left Ventricle diastolic major axis")                          |                                   | Measurement   |
|         | LVLs Apical      | (18076-0, LN, "Left Ventricle systolic major axis")                           |                                   | Method=DCM\125228\Bullet Method                         |
|         | LVEDV Bullet     | (18026-5, LN, "Left Ventricular End Diastolic Volume")                        |                                   |   |
|         | LVESV Bullet     | (18148-7, LN, "Left Ventricular End Systolic Volume")                         |                                   |   |
|         | LV SV Bullet     | (F-32120, SRT, "Stroke Volume")   |                                   |   |
|         | LV CO Bullet     | (F-32100, SRT, "Cardiac Output")  |                                   |   |
|         | LV EF Bullet     | (18043-0, LN, "Left Ventricular Ejection Fraction")                           |                                   |   |
|         | LV SI Bullet     | (F-00078, SRT, "Stroke Index")  |                                   |   |
|         | LV CI Bullet     | (F-32110, SRT, "Cardiac Index")   |                                   |   |
| LV Mass | LVAAd SAX PM Epi | (G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view")      | (T-32600, SRT, "Left Ventricle")  |   |
|         | LVAAd SAX PM     | (99102-01, MDSN, "Left Ventricle Epicardial Diastolic Area")                  |                                   | Image Mode=SRT\G-03A2\2D mode                           |
|         | LVLd apical      | (18077-8, LN, "Left Ventricle diastolic major axis")                          |                                   | (AL)  |
|         | LV Myo Thick     | (99104-03, MDSN, "Myocardial Thickness")                                      |                                   | Measurement Method =DCM\125205\Area-Length Single Plane |
|         | LV MI AL         | (18087-7, LN, "Left Ventricle Mass")  |                                   |   |
|         | LV TE a          | (99104-01, MDSN, "Left Ventricular Mass Index")                               |                                   | (TE)  |
|         | LV TE d          | (G-0377, SRT, "Left Ventricle Semi-major Axis Diastolic Dimension")           |                                   | Measurement Method =DCM\125222\Left                     |
|         | LV Mass TE       | (G-0378, SRT, "Left Ventricle Truncated Semi-major Axis Diastolic Dimension") |                                   | Ventricle Mass by Truncated Ellipse                     |
|         | LV MI TE         | (99104-01, MDSN, "Left Ventricular Mass Index")                               |                                   |   |
| RV (2D) | RVAWd            | (18153-7, LN, "Right Ventricle Anterior Wall Diastolic Thickness")            | (T-32500, SRT, "Right Ventricle") | Image Mode=SRT\G-03A2\2D mode                           |

|                   |             |  |                                   |                               |
|-------------------|-------------|--|-----------------------------------|-------------------------------|
|                   | RVIDd       | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension")                  |                                   |                               |
|                   | RVAd        | (99105-01, MDSN, "Right Ventricular Diastolic Area")                             |                                   |                               |
|                   | RVAWs       | (18157-8, LN, "Right Ventricular Anterior Wall Systolic Thickness")              |                                   |                               |
|                   | RVIDs       | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")                   |                                   |                               |
|                   | MPA Diam    | (18020-8, LN, "Main Pulmonary Artery Diameter")                                  |                                   |                               |
|                   | RVAs        | (99105-02, MDSN, "Right Ventricular Systolic Area")                              |                                   |                               |
|                   | RV Major    | (99105-03, MDSN, "Right Ventricular Major Diameter")                             |                                   |                               |
|                   | RV Minor    | (99105-04, MDSN, "Right Ventricular Minor Diameter")                             |                                   |                               |
|                   | RV FAC      | (79936-1, LN, Right ventricular fractional area change")                         |                                   |                               |
| RV (M)            | RVAWd       | (18153-7, LN, "Right Ventricle Anterior Wall Diastolic Thickness")               | (T-32500, SRT, "Right Ventricle") | Image Mode=SRT\G-0394\M mode  |
|                   | RVIDd       | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension")                  |                                   |                               |
|                   | RVAWs       | (18157-8, LN, "Right Ventricular Anterior Wall Systolic Thickness")              |                                   |                               |
|                   | RVIDs       | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")                   |                                   |                               |
|                   | TAPSE       | (99105-16, MDSN, "Tricuspid Annular Plane Systolic Excursion")                   |                                   |                               |
|                   | RVPEP       | (20280-4, LN, "Pressure Half-Time")  |                                   |                               |
|                   | RVET        | (99105-05, MDSN, "Ejection Time")  |                                   |                               |
|                   | RVPEP/ET(M) | (99105-06, MDSN, "Ratio of Right Ventricle Pressure Half-Time to Ejection Time") |                                   |                               |
| RV Vol. (Simpson) | RVLd A4C    | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension")                  | (T-32500, SRT, "Right Ventricle") | Image Mode=SRT\G-03A2\2D mode |
|                   | RVAd A4C    | (99105-01, MDSN, "Right Ventricular  |                                   |                               |

|                 |   |
|-----------------|---|
|                 | Diastolic Area")  |
| RVEDV A4C       | (8822-9, LN, "Right Ventricular ED Volume")                     |
| RVEDV Index A4C | (99105-22, MDSN, "Right Ventricular ED Volume Index")           |
| RVLs A4C        | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")  |
| RVAs A4C        | (99105-02, MDSN, "Right Ventricular Systolic Area")             |
| RVESV A4C       | (99105-02, MDSN, "Right Ventricular Systolic Area")             |
| RVESV Index A4C | (99105-23, MDSN, "Right Ventricular ES Volume Index")           |
| RV SV A4C       | (F-32120, SRT, "Stroke Volume")                                 |
| RV CO A4C       | (F-32100, SRT, "Cardiac Output")                                |
| RV EF A4C       | (18043-0, LN, "Left Ventricular Ejection Fraction")             |
| RV SI A4C       | (F-00078, SRT, "Stroke Index")                                  |
| RV CI A4C       | (F-32110, SRT, "Cardiac Index")                                 |
| RVLd A2C        | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") |
| RVAd A2C        | (99105-01, MDSN, "Right Ventricular Diastolic Area")            |
| RVEDV A2C       | (8822-9, LN, "Right Ventricular ED Volume")                     |
| RVEDV Index A2C | (99105-22, MDSN, "Right Ventricular ED Volume Index")           |
| RVLs A2C        | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")  |
| RVAs A2C        | (99105-02, MDSN, "Right Ventricular Systolic Area")             |
| RVESV A2C       | (8824-5, LN, "Right Ventricular ES Volume")                     |
| RVESV Index A2C | (99105-23, MDSN, "Right Ventricular ES Volume Index")           |

(A2C)Image  
View=SRT\G-  
A19B\Apical two chamber

(A4C)  
Image View=SRT\G-  
A19C\Apical fourchamber

(Simpson)  
MeasurementMethod=DC  
M\125208\Method of  
Disks, SinglePlane

(Simpson BP)  
MeasurementMethod=DC  
M\125207\Method of  
Disks, Biplane

|                  |                    |   |   |
|------------------|--------------------|---|---|
|                  |                    |   |   |
|                  | RV SV A2C          | (F-32120, SRT, "Stroke Volume")                                 |   |
|                  | RV CO A2C          | (F-32100, SRT, "Cardiac Output")                                |   |
|                  | RV EF A2C          | (18043-0, LN, "Left Ventricular Ejection Fraction")             |   |
|                  | RV SI A2C          | (F-00078, SRT, "Stroke Index")                                  |   |
|                  | RV CI A2C          | (F-32110, SRT, "Cardiac Index")                                 |   |
|                  | RVEDV BP           | (8824-5, LN, "Right Ventricular ES Volume")                     |   |
|                  | RVEDV Index BP     | (99105-23, MDSN, "Right Ventricular ES Volume Index")           |   |
|                  | RVESV BP           | (8824-5, LN, "Right Ventricular ES Volume")                     |   |
|                  | RVESV Index BP     | (99105-23, MDSN, "Right Ventricular ES Volume Index")           |   |
|                  | RV SV BP           | (F-32120, SRT, "Stroke Volume")                                 |   |
|                  | RV CO BP           | (F-32100, SRT, "Cardiac Output")                                |   |
|                  | RV EF BP           | (18043-0, LN, "Left Ventricular Ejection Fraction")             |   |
|                  | RV SI BP           | (F-00078, SRT, "Stroke Index")                                  |   |
|                  | RV CI BP           | (F-32110, SRT, "Cardiac Index")                                 |   |
| RV Vol.<br>(A/L) | RVLd A4C           | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle")             |
|                  | RVAd A4C           | (99105-01, MDSN, "Right Ventricular Diastolic Area")            | Image Mode=SRT\G-03A2\2D mode                 |
|                  | RVEDV A4C AL       | (8822-9, LN, "Right Ventricular ED Volume")                     | (A2C)Image View=SRT\G-A19B\Apical two chamber |
|                  | RVEDV Index A4C AL | (99105-22, MDSN, "Right Ventricular ED Volume Index")           | (A4C)   |
|                  | RVLs A4C           | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")  | Image View=SRT\G-A19C\Apical fourchamber      |
|                  | RVAs A4C           | (99105-02, MDSN, "Right Ventricular Systolic Area")             | (A/L)Measurement                              |
|                  | RVESV A4C AL       | (99105-02, MDSN, "Right Ventricular Systolic Area")             | Method=DCM\125205\Area-Length Single Plane    |

|    |                       |  |                               |                                       |
|----|-----------------------|--|-------------------------------|---------------------------------------|
|    | RVESV Index<br>A4C AL | (99105-23, MDSN, "Right Ventricular ES Volume Index")            |                               |                                       |
|    | RV SV A4C             | (F-32120, SRT, "Stroke Volume")                                  |                               |                                       |
|    | RV CO A4C             | (F-32100, SRT, "Cardiac Output")                                 |                               |                                       |
|    | RV EF A4C             | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                               |                                       |
|    | RV SI A4C             | (F-00078, SRT, "Stroke Index")                                   |                               |                                       |
|    | RV CI A4C             | (F-32110, SRT, "Cardiac Index")                                  |                               |                                       |
|    | RVLd A2C              | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension")  |                               |                                       |
|    | RVAd A2C              | (99105-01, MDSN, "Right Ventricular Diastolic Area")             |                               |                                       |
|    | RVEDV A2C<br>AL       | (8822-9, LN, "Right Ventricular ED Volume")                      |                               |                                       |
|    | RVEDV Index<br>A2C AL | (99105-22, MDSN, "Right Ventricular ED Volume Index")            |                               |                                       |
|    | RVLs A2C              | (20305-9, LN, "Right Ventricular Internal Systolic Dimension")   |                               |                                       |
|    | RVAs A2C              | (99105-02, MDSN, "Right Ventricular Systolic Area")              |                               |                                       |
|    | RVESV A2C             | (8824-5, LN, "Right Ventricular ES Volume")                      |                               |                                       |
|    | RVESV Index<br>A2C    | (99105-23, MDSN, "Right Ventricular ES Volume Index")            |                               |                                       |
|    | RV SV A2C             | (F-32120, SRT, "Stroke Volume")                                  |                               |                                       |
|    | RV CO A2C             | (F-32100, SRT, "Cardiac Output")                                 |                               |                                       |
|    | RV EF A2C             | (18043-0, LN, "Left Ventricular Ejection Fraction")              |                               |                                       |
|    | RV SI A2C             | (F-00078, SRT, "Stroke Index")                                   |                               |                                       |
|    | RV CI A2C             | (F-32110, SRT, "Cardiac Index")                                  |                               |                                       |
| LA | LA Diam(2D)           | (29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension") | (T-32300, SRT, "Left Atrium") | (2D)<br>Image Mode=SRT\G-03A2\2D mode |
|    | LA/Ao(2D)             | (17985-3, LN, "Left Atrium to Aortic Root Ratio")                |                               |                                       |

|                   |                     |  |                               |  |
|-------------------|---------------------|--|-------------------------------|--|
|                   | LA Diam(M)          | (29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension") |                               | (M)                                      |
|                   | LA/Ao(M)            | (17985-3, LN, "Left Atrium to Aortic Root Ratio")                |                               | Image Mode=SRT\G-0394\M mode             |
|                   | LA Major            | (99106-01, MDSN, "Left Atrium Major Diameter")                   |                               |  |
|                   | LA Minor            | (99106-02, MDSN, "Left Atrium Minor Diameter")                   |                               |  |
|                   | LA Area             | (17977-0, LN, "Left Atrium Systolic Area")                       |                               |  |
|                   | LA Volume(2D)       | (G-0383, SRT, "Left Atrium Systolic Volume")                     |                               |  |
|                   | LA Volume Index(2D) | (99106-07, MDSN, "Left Atrium Volume Index")                     |                               |  |
|                   | LA Volume(M)        | (G-0383, SRT, "Left Atrium Systolic Volume")                     |                               |  |
|                   | LA Volume Index(M)  | (99106-07, MDSN, "Left Atrium Volume Index")                     |                               |  |
|                   | LAAAd A2C           | (99106-03, MDSN, "Left Atrium Diastolic Area")                   |                               |  |
|                   | LAAAs A2C           | (17977-0, LN, "Left Atrium Systolic Area")                       |                               |  |
|                   | LAAAd A4C           | (99106-03, MDSN, "Left Atrium Diastolic Area")                   |                               |  |
|                   | LAAAs A4C           | (17977-0, LN, "Left Atrium Systolic Area")                       |                               |  |
|                   | LAEDV BP AL         | (122407, DCM, "Left Atrial ED Volume")                           |                               |  |
|                   | LAESV BP AL         | (122408, DCM, "Left Atrial ES Volume")                           |                               |  |
| LA Vol. (Simpson) | LALd A4C            | (99106-08, MDSN, "Left Atrium Diastolic Length")                 | (T-32300, SRT, "Left Atrium") | Image Mode =SRT\G-03A2\2D mode           |
|                   | LAAAd A4C           | (99106-03, MDSN, "Left Atrium Diastolic Area")                   |                               | (A2C)                                    |
|                   | LAEDV A4C           | (122407, DCM, "Left Atrial ED Volume")                           |                               | Image View=SRT\G-A19B\Apical two chamber |
|                   | LAEDV Index A4C     | (99106-10, MDSN, "Left Atrial ED Volume Index")                  |                               |  |
|                   | LALs A4C            | (99106-09, MDSN, "Left Atrium Systolic                           |                               | (A4C)                                    |

|              |                    |  |                               |   |
|--------------|--------------------|--|-------------------------------|---|
|              |                    | Length")   |                               | Image View=SRT\G-A19C\Apical four chamber   |
|              | LAA A4C            | (17977-0, LN, "Left Atrium Systolic Area")       |                               |   |
|              | LAESV A4C          | (122408, DCM, "Left Atrial ES Volume")           |                               |   |
|              | LAESV Index A4C    | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                               | (Simpson)<br>Measurement Method =DCM\125208\Method of Disks, Single Plane               |
|              | LALd A2C           | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                               |   |
|              | LAA d A2C          | (99106-08, MDSN, "Left Atrium Diastolic Length") |                               | (Simpson BP)<br>Measurement Method =DCM\125207\Method of Disks, Biplane                 |
|              | LAEDV A2C          | (122407, DCM, "Left Atrial ED Volume")           |                               |   |
|              | LAEDV Index A2C    | (99106-10, MDSN, "Left Atrial ED Volume Index")  |                               |   |
|              | LALs A2C           | (99106-09, MDSN, "Left Atrium Systolic Length")  |                               |   |
|              | LAA A2C            | (17977-0, LN, "Left Atrium Systolic Area")       |                               |   |
|              | LAESV A2C          | (122408, DCM, "Left Atrial ES Volume")           |                               |   |
|              | LAESV Index A2C    | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                               |   |
|              | LAEDV BP           | (122407, DCM, "Left Atrial ED Volume")           |                               |   |
|              | LAEDV Index BP     | (99106-10, MDSN, "Left Atrial ED Volume Index")  |                               |   |
|              | LAESV BP           | (122408, DCM, "Left Atrial ES Volume")           |                               |   |
|              | LAESV Index BP     | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                               |   |
| LA Vol.(A/L) | LALd A4C           | (99106-08, MDSN, "Left Atrium Diastolic Length") | (T-32300, SRT, "Left Atrium") | Image Mode =SRT\G-03A2\2D mode<br><br>(A2C)<br>Image View=SRT\G-A19B\Apical two chamber |
|              | LAA d A4C          | (99106-03, MDSN, "Left Atrium Diastolic Area")   |                               |   |
|              | LAEDV A4C AL       | (122407, DCM, "Left Atrial ED Volume")           |                               |   |
|              | LAEDV Index A4C AL | (99106-10, MDSN, "Left Atrial ED Volume Index")  |                               | (A4C)   |
|              | LALs A4C           | (99106-09, MDSN, "Left Atrium Systolic           |                               | Image View=SRT\G-   |

|       |                    |  |                         |  |
|-------|--------------------|--|-------------------------|--|
|       |                    | Length")   |                         | A19C\Apical four chamber                                     |
|       | LAA A4C            | (17977-0, LN, "Left Atrium Systolic Area")       |                         |  |
|       | LAESV A4C AL       | (122408, DCM, "Left Atrial ES Volume")           |                         | (Simpson)  |
|       | LAESV Index A4C AL | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                         | Measurement Method =DCM\125208\Method of Disks, Single Plane |
|       | LALd A2C           | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                         |  |
|       | LAA d A2C          | (99106-08, MDSN, "Left Atrium Diastolic Length") |                         | (Simpson BP)   |
|       | LAEDV A2C AL       | (122407, DCM, "Left Atrial ED Volume")           |                         | Measurement Method =DCM\125207\Method of Disks, Biplane      |
|       | LAEDV Index A2C AL | (99106-10, MDSN, "Left Atrial ED Volume Index")  |                         |  |
|       | LALs A2C           | (99106-09, MDSN, "Left Atrium Systolic Length")  |                         |  |
|       | LAA A2C            | (17977-0, LN, "Left Atrium Systolic Area")       |                         |  |
|       | LAESV A2C AL       | (122408, DCM, "Left Atrial ES Volume")           |                         |  |
|       | LAESV Index A2C AL | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                         |  |
|       | LAEDV BP AL        | (122407, DCM, "Left Atrial ED Volume")           |                         |  |
|       | LAEDV Index BP AL  | (99106-10, MDSN, "Left Atrial ED Volume Index")  |                         |  |
|       | LAESV BP AL        | (122408, DCM, "Left Atrial ES Volume")           |                         |  |
|       | LAESV Index BP AL  | (99106-11, MDSN, "Left Atrial ES Volume Index")  |                         |  |
| Aorta | Aorta Annulus      | (79940-3, LN, "Aortic annulus diameter")         | (T-42000, SRT, "Aorta") | (2D)   |
|       | Ao Diam(2D)        | (18015-8, LN, "Aortic Root Diameter")            |                         | Image Mode=SRT\G-03A2\2D mode                                |
|       | Ao Diam(M)         | (18015-8, LN, "Aortic Root Diameter")            |                         | (M)  |
|       | Ao Arch Diam       | (18011-7, LN, "Aortic Arch Diameter")            |                         | Image Mode=SRT\G-0394\M mode                                 |
|       | Asc Ao Diam        | (18012-5, LN, "Ascending Aortic Diameter")       |                         |  |
|       | Desc Ao Diam       | (18013-3, LN, "Descending Aortic Diameter")      |                         |  |
|       | Ao Isth Diam       | (18014-1, LN, "Aortic Isthmus Diameter")         |                         |  |

|                   |                  |   |                                |  |
|-------------------|------------------|---|--------------------------------|--|
|                   | Ao ST Junct Diam | (99109-01, MDSN, "Aortic ST Junct Diameter")                                  |                                |  |
|                   | Ao Sinus Diam    | (99109-02, MDSN, "Aortic Sinus Diameter")                                     |                                |  |
|                   | Ao ST/Ao         | (99109-03, MDSN, "Ratio of Aortic ST Junct Diameter to Aortic Root Diameter") |                                |  |
| RA                | RA Major         | (99107-02, MDSN, "Right Atrium Major Diameter")                               | (T-32200, SRT, "Right Atrium") | (2D)<br>Image Mode=SRT\G-03A2\2D mode  |
|                   | RA Minor         | (99107-03, MDSN, "Right Atrium Minor Diameter")                               |                                |  |
|                   | RA Area          | (99107-12, MDSN, "Right Atrium Area")   |                                |  |
|                   | RAAd             | (99107-01, MDSN, "Right Atrium Diastolic Area")                               |                                |  |
|                   | RAAs             | (17988-7, LN, "Right Atrium Systolic Area")                                   |                                |  |
|                   | RAEDV            | (99107-04, MDSN, "Right Atrium Diastolic Volume")                             |                                |  |
|                   | RAESV            | (99107-05, MDSN, "Right Atrium Systolic Volume")                              |                                |  |
| RA Vol. (Simpson) | RALd A4C         | (99107-07, MDSN, "Right Atrium Diastolic Length")                             | (T-32200, SRT, "Right Atrium") | Image Mode =SRT\G-03A2\2D mode<br><br>(A2C)<br>Image View=SRT\G-A19B\Apical two chamber<br><br>(A4C)<br>Image View=SRT\G-A19C\Apical four chamber<br><br>(Simpson)<br>Measurement Method =DCM\125208\Method of Disks, Single Plane |
|                   | RAAd A4C         | (99107-01, MDSN, "Right Atrium Diastolic Area")                               |                                |  |
|                   | RAEDV A4C        | (99107-04, MDSN, "Right Atrium Diastolic Volume")                             |                                |  |
|                   | RAEDV Index A4C  | (99107-10, MDSN, "Right Atrial ED Volume Index")                              |                                |  |
|                   | RALs A4C         | (99107-08, MDSN, "Right Atrium Systolic Length")                              |                                |  |
|                   | RAAs A4C         | (17988-7, LN, "Right Atrium Systolic Area")                                   |                                |  |
|                   | RAESV A4C        | (99107-05, MDSN, "Right Atrium Systolic Volume")                              |                                |  |
|                   | RAESV Index A4C  | (99107-11, MDSN, "Right Atrial ES Volume Index")                              |                                |  |
|                   | RALd A2C         | (99107-07, MDSN, "Right Atrium  |                                |  |

|               |                    |   |                                |   |  |
|---------------|--------------------|---|--------------------------------|---|--|
|               |                    | Diastolic Length")                                |                                | (Simpson BP)  |  |
|               | RAAd A2C           | (99107-01, MDSN, "Right Atrium Diastolic Area")   |                                | Measurement Method<br>=DCM\125207\Method of<br>Disks, Biplane |  |
|               | RAEDV A2C          | (99107-04, MDSN, "Right Atrium Diastolic Volume") |                                |   |  |
|               | RAEDV Index A2C    | (99107-10, MDSN, "Right Atrial ED Volume Index")  |                                |   |  |
|               | RALs A2C           | (99107-08, MDSN, "Right Atrium Systolic Length")  |                                |   |  |
|               | RAAs A2C           | (17988-7, LN, "Right Atrium Systolic Area")       |                                |   |  |
|               | RAESV A2C          | (99107-05, MDSN, "Right Atrium Systolic Volume")  |                                |   |  |
|               | RAESV Index A2C    | (99107-11, MDSN, "Right Atrial ES Volume Index")  |                                |   |  |
|               | RAEDV BP           | (99107-04, MDSN, "Right Atrium Diastolic Volume") |                                |   |  |
|               | RAEDV Index BP     | (99107-10, MDSN, "Right Atrial ED Volume Index")  |                                |   |  |
|               | RAESV BP           | (99107-05, MDSN, "Right Atrium Systolic Volume")  |                                |   |  |
|               | RAESV Index BP     | (99107-11, MDSN, "Right Atrial ES Volume Index")  |                                |   |  |
| RA Vol. (A/L) | RALd A4C           | (99107-07, MDSN, "Right Atrium Diastolic Length") | (T-32200, SRT, "Right Atrium") |   |  |
|               | RAAd A4C           | (99107-01, MDSN, "Right Atrium Diastolic Area")   |                                |   |  |
|               | RAEDV A4C AL       | (99107-04, MDSN, "Right Atrium Diastolic Volume") |                                |   |  |
|               | RAEDV Index A4C AL | (99107-10, MDSN, "Right Atrial ED Volume Index")  |                                |   |  |
|               | RALs A4C           | (99107-08, MDSN, "Right Atrium Systolic Length")  |                                |   |  |
|               | RAAs A4C           | (17988-7, LN, "Right Atrium Systolic Area")       |                                |   |  |
|               | RAESV A4C          | (99107-05, MDSN, "Right Atrium Systolic           |                                |   |  |

|      |                       |  |                                      |  |
|------|-----------------------|--|--------------------------------------|--|
|      | AL                    | Volume")   |                                      |  |
|      | RAESV Index<br>A4C AL | (99107-11, MDSN, "Right Atrial ES<br>Volume Index")  |                                      |  |
|      | RALd A2C              | (99107-07, MDSN, "Right Atrium<br>Diastolic Length") |                                      |  |
|      | RAAd A2C              | (99107-01, MDSN, "Right Atrium<br>Diastolic Area")   |                                      |  |
|      | RAEDV A2C<br>AL       | (99107-04, MDSN, "Right Atrium<br>Diastolic Volume") |                                      |  |
|      | RAEDV Index<br>A2C AL | (99107-10, MDSN, "Right Atrial ED<br>Volume Index")  |                                      |  |
|      | RALs A2C              | (99107-08, MDSN, "Right Atrium Systolic<br>Length")  |                                      |  |
|      | RAAs A2C              | (17988-7, LN, "Right Atrium Systolic<br>Area")       |                                      |  |
|      | RAESV A2C<br>AL       | (99107-05, MDSN, "Right Atrium Systolic<br>Volume")  |                                      |  |
|      | RAESV Index<br>A2C AL | (99107-11, MDSN, "Right Atrial ES<br>Volume Index")  |                                      |  |
| LVOT | LVOT Diam             | (G-038F, SRT, "Cardiovascular Orifice<br>Diameter")  | (T-32600, SRT, "Left<br>Ventricle")  | FSite=SRT\T-32650\Left<br>Ventricle Outflow Tract      |
|      | LVOT Area             | (G-038E, SRT, "Cardiovascular Orifice<br>Area")      |                                      |  |
|      | LVOT Vmax             | (11726-7, LN, "Peak Velocity")                       |                                      |  |
|      | LVOT Pgmax            | (20247-3, LN, "Peak Gradient")                       |                                      |  |
|      | LVOT Vmean            | (20352-1, LN, "Mean Velocity")                       |                                      |  |
|      | LVOT Pgmean           | (20256-4, LN, "Mean Gradient")                       |                                      |  |
|      | LVOT VTI              | (20354-7, LN, "Velocity Time Integral")              |                                      |  |
|      | LVOT SV               | (F-32120, SRT, "Stroke Volume")                      |                                      |  |
|      | LVOT CO               | (F-32100, SRT, "Cardiac Output")                     |                                      |  |
|      | LVOT CI               | (F-32110, SRT, "Cardiac Index")                      |                                      |  |
|      | LVOT SI               | (F-00078, SRT, "Stroke Index")                       |                                      |  |
| RVOT | RVOT Diam             | (G-038F, SRT, "Cardiovascular Orifice<br>Diameter")  | (T-32500, SRT, "Right<br>Ventricle") | FSite=SRT\T-<br>32550\Right Ventricle<br>Outflow Tract |
|      | RVOT Area             | (G-038E, SRT, "Cardiovascular Orifice<br>Area")      |                                      |  |

|    |                |  |                                |  |
|----|----------------|--|--------------------------------|--|
|    |                | Area")   |                                |  |
|    | RVOT Vmax      | (11726-7, LN, "Peak Velocity")                     |                                |  |
|    | RVOT Vmean     | (20352-1, LN, "Mean Velocity")                     |                                |  |
|    | RVOT Pgmax     | (20247-3, LN, "Peak Gradient")                     |                                |  |
|    | RVOT Pgmean    | (20256-4, LN, "Mean Gradient")                     |                                |  |
|    | RVOT VTI       | (20354-7, LN, "Velocity Time Integral")            |                                |  |
|    | RVOT SV        | (F-32120, SRT, "Stroke Volume")                    |                                |  |
|    | RVOT CO        | (F-32100, SRT, "Cardiac Output")                   |                                |  |
|    | RVOT CI        | (F-32110, SRT, "Cardiac Index")                    |                                |  |
|    | RVOT SI        | (F-00078, SRT, "Stroke Index")                     |                                |  |
| AV | AV Cusp        | (17996-0, LN, "Aortic Valve Cusp Separation")      | (T-35400, SRT, "Aortic Valve") | (2D)<br>Image Mode=SRT\G-03A2\2D mode                            |
|    | AV Cusp(M)     | (17996-0, LN, "Aortic Valve Cusp Separation")      |                                | (M)<br>Image Mode=SRT\G-0394\M mode                              |
|    | AV Diam        | (G-038F, SRT, "Cardiovascular Orifice Diameter")   |                                | (Color Doppler)<br>Image Mode=SRT\R-409E2\Doppler Color Flow     |
|    | AV Area        | (G-038E, SRT, "Cardiovascular Orifice Area")       |                                | (AV)<br>FDirect=SRT\R-42047\Antegrade Flow                       |
|    | AVA Planimetry | (G-038E, SRT, "Cardiovascular Orifice Area")       |                                | (AR)<br>FDirect=SRT\R-42E61\Regurgitant Flow                     |
|    | AVA Index      | (99108-11, MDSN, "AVA Index")                      |                                | Measurement Method =DCM\125216\Proximal Isovelocity Surface Area |
|    | AV Vmax        | (11726-7, LN, "Peak Velocity")                     |                                |  |
|    | AVA(Vmax)      | (G-038E, SRT, "Cardiovascular Orifice Area")       |                                |  |
|    | AV PGmax       | (20247-3, LN, "Peak Gradient")                     |                                |  |
|    | AV Vmean       | (20352-1, LN, "Mean Velocity")                     |                                |  |
|    | AV PGmean      | (20256-4, LN, "Mean Gradient")                     |                                |  |
|    | AV PHT         | (20280-4, LN, "Pressure Half-Time")                |                                |  |
|    | AV VTI         | (20354-7, LN, "Velocity Time Integral")            |                                |  |
|    | AV AccT        | (20168-1, LN, "Acceleration Time")                 |                                |  |
|    | AV Acc         | (99108-01, MDSN, "Aortic Valve Flow Acceleration") |                                |  |
|    | AV DecT        | (20217-6, LN, "Deceleration Time")                 |                                |  |
|    | AV Dec         | (20216-8, LN, "Deceleration Slope")                |                                |  |
|    | AV ET          | (18041-4, LN, Aortic Valve Ejection Time)          |                                |  |

|                   |   |  |
|-------------------|---|--|
| AV AccT/ET        | (G-0382, SRT, "Ratio of Aortic Valve Acceleration Time to Ejection Time") | Measurement Method<br>=DCM\125220\Planimetry |
| R-R Int           | (99108-02, MDSN, "Aortic R to R Interval")                                | Measurement Method<br>=DCM\125214\Continuity |
| AVA (VTI)         | (G-038E, SRT, "Cardiovascular Orifice Area")                              | Equation by Peak Velocity                    |
| LVOT VTI / AV VTI | (99108-12, MDSN, "LVOT to AV Velocity Time Integral Ratio")               | Velocity                                     |
| AR VCW            | (99108-08, MDSN, "Aortic Regurgitant Vena Contracta Width")               | Measurement Method<br>=DCM\125215\Continuity |
| AR Vmax           | (11726-7, LN, "Peak Velocity")  | Equation by Velocity                         |
| AR Vmean          | (20352-1, LN, "Mean Velocity")  | Time Integral                                |
| AR PGmax          | (20247-3, LN, "Peak Gradient")  |  |
| AR PGmean         | (20256-4, LN, "Mean Gradient")  | Measurement Method<br>=DCM\125210\Area by    |
| AR ed Vmax        | (99108-03, MDSN, "Aortic Regurgitant End-Diastolic Velocity")             | Pressure Half-Time                           |
| AR ed PGmax       | (99108-04, MDSN, "Aortic Regurgitant End-Diastolic Pressure Gradient")    |  |
| AR PHT            | (20280-4, LN, "Pressure Half-Time")                                       |  |
| AR VTI            | (20354-7, LN, "Velocity Time Integral")                                   |  |
| AR AccT           | (20168-1, LN, "Acceleration Time")  |  |
| AR Acc            | (99108-06, MDSN, "Aortic Regurgitant Flow Acceleration")                  |  |
| AR DecT           | (20217-6, LN, "Deceleration Time")  |  |
| AR Dec            | (20216-8, LN, "Deceleration Slope")                                       |  |
| AR PISA Rad       | (G-038F, SRT, "Cardiovascular Orifice Diameter")                          |  |
| AR Alias Vel.     | (99108-05, MDSN, "Aortic Regurgitant Aliasing Velocity")                  |  |
| AR Flow Rate      | (34141-2, LN, "Peak Instantaneous Flow Rate")                             |  |
| AR ERO            | (G-038E, SRT, "Cardiovascular Orifice Area")                              |  |
| AR Volume         | (33878-0, LN, "Volume Flow")  |  |
| AR Fraction       | (G-0390-4, SRT, "Regurgitant Fraction")                                   |  |

|    |                      |  |                                |  |
|----|----------------------|--|--------------------------------|--|
|    | AVO                  | (99108-09, MDSN, "AVO")                                    |                                |  |
|    | AVC                  | (99108-10, MDSN, "AVC")                                    |                                |  |
|    | AR Jet Diam          | (125333, DCM, "Regurgitation Jet Width")                   |                                |  |
|    | AR jet ratio by dia. | (79944-5, LN, "Aortic regurgitant jet width/LVOT width %") |                                |  |
|    | AR Jet Area          | (125332, DCM, "Regurgitation Jet Area")                    |                                |  |
|    | AR jet ratio by area | (79943-7, LN, "Aortic regurgitant jet area/LVOT area %")   |                                |  |
| MV | MV D-E Excursion     | (99114-02, MDSN, "Mitral Valve D-E Excursion")             | (T-35300, SRT, "Mitral Valve") | (2D)<br>Image Mode=SRT\G-03A2\2D mode                            |
|    | MV D-E Slope         | (99114-03, MDSN, "Mitral Valve D-E Slope")                 |                                |  |
|    | MV E-F Slope         | (18040-6, LN, "Mitral Valve E-F Slope by M-Mode")          |                                | (M)<br>Image Mode=SRT\G-0394\M mode                              |
|    | MV A-C Interval      | (99114-01, MDSN, "Mitral Valve A-C Interval")              |                                |  |
|    | MV EPSS              | (18036-4, LN, "Mitral Valve EPSS, E wave")                 |                                | (Color Doppler)<br>Image Mode=SRT\R-409E2\Doppler Color Flow     |
|    | MV Ann Diam          | (G-038F, SRT, "Cardiovascular Orifice Diameter")           |                                |  |
|    | MV Diam1             | (G-038F, SRT, "Cardiovascular Orifice Diameter")           |                                | (MV)   |
|    | MV Diam2             | (G-038F, SRT, "Cardiovascular Orifice Diameter")           |                                | FDirect=SRT\R-42047\Antegrade Flow                               |
|    | MV Area              | (G-038E, SRT, "Cardiovascular Orifice Area")               |                                | (MR)   |
|    | MV Area (Ann Diam)   | (G-038E, SRT, "Cardiovascular Orifice Area")               |                                | FDirect=SRT\R-42E61\Regurgitant Flow                             |
|    | MVA Planimetry       | (G-038E, SRT, "Cardiovascular Orifice Area")               |                                | Measurement Method =DCM\125216\Proximal Isovelocity Surface Area |
|    | MV Vp                | (99114-04, MDSN, "Mitral Valve Propagation Velocity")      |                                |  |
|    | MV Peak E            | (18037-2, LN, "Mitral Valve E-Wave Peak Velocity")         |                                | Measurement Method   |

|                     |   |
|---------------------|---|
| MV Peak A           | (17978-8, LN, "Mitral Valve A-Wave Peak")                           |
| MV E/A              | (18038-0, LN, "Mitral Valve E to A Ratio")                          |
| MV Vmax             | (11726-7, LN, "Peak Velocity")                                      |
| MVA(Vmax)           | (G-038E, SRT, "Cardiovascular Orifice Area")                        |
| MV PGmax            | (18057-0, LN, "Mitral Valve Diastolic Peak Instantaneous Gradient") |
| MV Vmean            | (20352-1, LN, "Mean Velocity")                                      |
| MV PGmean           | (20256-4, LN, "Mean Gradient")                                      |
| MV PHT              | (20280-4, LN, "Pressure Half-Time")                                 |
| MVA(PHT)            | (G-038E, SRT, "Cardiovascular Orifice Area")                        |
| MV VTI              | (20354-7, LN, "Velocity Time Integral")                             |
| MVA(VTI)            | (G-038E, SRT, "Cardiovascular Orifice Area")                        |
| MV AccT             | (20168-1, LN, "Acceleration Time")                                  |
| MV Acc              | (99114-06, MDSN, "Mitral Valve Flow Acceleration")                  |
| MV DecT             | (20217-6, LN, "Deceleration Time")                                  |
| MV Dec              | (20216-8, LN, "Deceleration Slope")                                 |
| MV AccT/DecT        | (G-0386, SRT, "Mitral Valve AT/DT Ratio")                           |
| MV A Dur            | (G-0385, SRT, "Mitral Valve A-Wave Duration")                       |
| P Vein A -MV A Dur. | (99110-01, MDSN, "P Vein A Dur minus MV A Dur")                     |
| MV SV               | (F-32120, SRT, "Stroke Volume")                                     |
| MV SV(Ann Diam)     | (F-32120, SRT, "Stroke Volume")                                     |
| MV CO               | (F-32100, SRT, "Cardiac Output")                                    |
| MV ET               | (99114-07, MDSN, "Mitral Valve Ejection Time")                      |
| R-R Int             | (99114-09, MDSN, "Mitral Valve R to R Interval")                    |

=DCM\125220\Planimetry  
Measurement Method  
=DCM\125214\Continuity  
Equation by Peak  
Velocity  
Measurement Method  
=DCM\125215\Continuity  
Equation by Velocity  
Time Integral  
Measurement Method  
=DCM\125210\Area by  
Pressure Half-Time

|    |                           |  |                                   |                                       |
|----|---------------------------|--|-----------------------------------|---------------------------------------|
|    | MR VCW                    | (99114-10, MDSN, "Mitral Regurgitant Vena Contracta Width")                            |                                   |                                       |
|    | MR Vmax                   | (11726-7, LN, "Peak Velocity")   |                                   |                                       |
|    | MR PGmax                  | (20247-3, LN, "Peak Gradient")   |                                   |                                       |
|    | MR Vmean                  | (20352-1, LN, "Mean Velocity")   |                                   |                                       |
|    | MR PGmean                 | (20256-4, LN, Mean Gradient")  |                                   |                                       |
|    | MR VTI                    | (20354-7, LN, "Velocity Time Integral")  |                                   |                                       |
|    | MR dp/dt Interval         | (99114-11, MDSN, "Mitral Regurgitant dp/dt Interval")                                  |                                   |                                       |
|    | MR dp/dt                  | (18035-6, LN, "Mitral Regurgitation dP/dt derived from Mitral Regurgitation velocity") |                                   |                                       |
|    | MR PISA Rad               | (G-038F, SRT, "Cardiovascular Orifice Diameter")                                       |                                   |                                       |
|    | MR Alias Vel.             | (99114-12, MDSN, "Mitral Regurgitant Aliasing Velocity")                               |                                   |                                       |
|    | MR Flow Rate              | (34141-2, LN, "Peak Instantaneous Flow Rate")  |                                   |                                       |
|    | MR ERO                    | (G-038E, SRT, "Cardiovascular Orifice Area")   |                                   |                                       |
|    | MR Volume                 | (33878-0, LN, "Volume Flow")   |                                   |                                       |
|    | MR Fraction               | (G-0390, SRT, "Regurgitant Fraction")  |                                   |                                       |
|    | MR Volume(Continuity VTI) | (80063-1, LN, "Mitral regurgitation volume (Continuity VTI)")                          |                                   |                                       |
|    | MR Jet Area               | (125332, DCM, "Regurgitation Jet Area")  |                                   |                                       |
|    | MR jet ratio by area      | (99114-13, MDSN, "Mitral regurgitant jet area/LA area %")                              |                                   |                                       |
|    | MV E/Vp                   | (99114-15, MSDN, "Mitral Valve E to Propagation Velocity Ratio")                       |                                   |                                       |
|    | MV DecT (PHT)             | (20217-6, LN, "Deceleration Time")   |                                   |                                       |
| TV | TV D-E Excursion          | (99115-09, MDSN, "Tricuspid Valve D-E Excursion")                                      | (T-35100, SRT, "Tricuspid Valve") | (2D)<br>Image Mode=SRT\G-03A2\2D mode |
|    | TV D-E Slope              | (99115-10, MDSN, "Tricuspid Valve D-E  |                                   |                                       |

|                 |   |
|-----------------|---|
|                 | Slope")   |
| TV E-F Slope    | (99115-11, MDSN, "Tricuspid Valve E-F Slope")         |
| TV A-C Interval | (99115-08, MDSN, "Tricuspid Valve A-C Interval")      |
| TV Ann Diam     | (G-038F, SRT, "Cardiovascular Orifice Diameter")      |
| TV Diam1        | (G-038F, SRT, "Cardiovascular Orifice Diameter")      |
| TV Diam2        | (G-038F, SRT, "Cardiovascular Orifice Diameter")      |
| TVA Planimetry  | (G-038E, SRT, "Cardiovascular Orifice Area")          |
| TV Area         | (G-038E, SRT, "Cardiovascular Orifice Area")          |
| TV Peak E       | (18031-5, LN, "Tricuspid Valve E Wave Peak Velocity") |
| TV Peak A       | (18030-7, LN, "Tricuspid Valve A Wave Peak Velocity") |
| TV E/A          | (18039-8, LN, "Tricuspid Valve E to A Ratio")         |
| TV Vmax         | (11726-7, LN, "Peak Velocity")                        |
| TVA (Vmax)      | (G-038E, SRT, "Cardiovascular Orifice Area")          |
| TV Vmean        | (20352-1, LN, "Mean Velocity")                        |
| TV PGmax        | (20247-3, LN, "Peak Gradient")                        |
| TV PGmean       | (20256-4, LN, "Mean Gradient")                        |
| TV PHT          | (20280-4, LN, "Pressure Half-Time")                   |
| TV VTI          | (20354-7, LN, "Velocity Time Integral")               |
| TVA(VTI)        | (G-038E, SRT, "Cardiovascular Orifice Area")          |
| TV AccT         | (20168-1, LN, "Acceleration Time")                    |
| TV Acc          | (99115-01, MDSN, "Tricuspid Valve Flow Acceleration") |
| TV DecT         | (20217-6, LN, "Deceleration Time")                    |

(M)  
Image Mode=SRT\G-0394\M mode

(Color Doppler)  
Image Mode=SRT\R-409E2\Doppler Color Flow

(TV)  
FDirect=SRT\R-42047\Antegrade Flow

(TR)  
FDirect=SRT\R-42E61\Regurgitant Flow

Measurement Method =DCM\125216\Proximal Isovelocity Surface Area

Measurement Method =DCM\125220\Planimetry

Measurement Method =DCM\125214\Continuity Equation by Peak Velocity

Measurement Method =DCM\125215\Continuity Equation by Velocity Time Integral

|               |  |
|---------------|--|
| TV Dec        | (20216-8, LN, "Deceleration Slope")                            |
| TV SV         | (F-32120, SRT, "Stroke Volume")                                |
| TV CO         | (F-32100, SRT, "Cardiac Output")                               |
| TV A Dur      | (99115-02, MDSN, "Tricuspid Valve A-Wave Duration")            |
| Q to TV Open  | (20296-0, LN, "Time from Q wave to Tricuspid Valve Opens")     |
| R-R Int       | (99115-12, MDSN, "Tricuspid Valve R to R Interval")            |
| TR VCW        | (99115-03, MDSN, "Tricuspid Regurgitant Vena Contracta Width") |
| TR Vmax       | (11726-7, LN, "Peak Velocity")                                 |
| TR PGmax      | (20247-3, LN, "Peak Gradient")                                 |
| TR Vmean      | (20352-1, LN, "Mean Velocity")                                 |
| TR PGmean     | (20256-4, LN, "Mean Gradient")                                 |
| TR VTI        | (20354-7, LN, "Velocity Time Integral")                        |
| RAP           | (18070-3, LN, "Right Atrium Systolic Pressure")                |
| TR RVSP       | (G-0380, SRT, "Right Ventricular Peak Systolic Pressure")      |
| TR dp/dt Int  | (99115-05, MDSN, "Tricuspid Regurgitant dp/dp Interval")       |
| TR dp/dt      | (18034-9, LN, "Tricuspid Regurgitation dP/dt")                 |
| TR PISA Rad   | (G-038F, SRT, "Cardiovascular Orifice Diameter")               |
| TR Alias Vel. | (99115-06, MDSN, "Tricuspid Regurgitant Aliasing Velocity")    |
| TR Flow Rate  | (34141-2, LN, "Peak Instantaneous Flow Rate")                  |
| TR ERO        | (G-038E, SRT, "Cardiovascular Orifice Area")                   |
| TR Volume     | (33878-0, LN, "Volume Flow")                                   |
| TR Fraction   | (G-0390, SRT, "Regurgitant Fraction")                          |
| TR Jet Area   | (125332, DCM, "Regurgitation Jet Area")                        |

Measurement Method  
=DCM\125210\Area by  
Pressure Half-Time

|               |   |   |                                  |   |
|---------------|---|---|----------------------------------|---|
|               | TR jet ratio by area  | (99115-12, MDSN, "Tricuspid regurgitant jet area/RA area %")                |                                  |   |
| PV            | PV Ann Diam   | (G-038F, SRT, "Cardiovascular Orifice Diameter")                            | (T-35200, SRT, "Pulmonic Valve") | (2D)<br>Image Mode=SRT\G-03A2\2D mode                               |
|               | PV Area   | (G-038E, SRT, "Cardiovascular Orifice Area")                                |                                  | (M)<br>Image Mode=SRT\G-0394\M mode                                 |
|               | PVA Planimetry  | (G-038E, SRT, "Cardiovascular Orifice Area")                                |                                  | (Color Doppler)<br>Image Mode=SRT\R-409E2\Doppler Color Flow        |
|               | PV Vmax   | (11726-7, LN, "Peak Velocity")  |                                  | (PV)<br>FDirect=SRT\R-42047\Antegrade Flow                          |
|               | PVA (Vmax)  | (G-038E, SRT, "Cardiovascular Orifice Area")                                |                                  | (PR)FDirect=SRT\R-42E61\Regurgitant Flow                            |
|               | PV Vmean  | (20352-1, LN, "Mean Velocity")  |                                  | Measurement Method =DCM\125216\Proximal Isovelocity Surface Area    |
|               | PV Pgmax  | (20247-3, LN, "Peak Gradient")  |                                  | Measurement Method =DCM\125220\Planimetry                           |
|               | PV Pgmean   | (20256-4, LN, "Mean Gradient")  |                                  | Measurement Method =DCM\125214\Continuity Equation by Peak Velocity |
|               | PV PHT  | (20280-4, LN, "Pressure Half-Time")   |                                  | Measurement Method  |
|               | PV VTI  | (20354-7, LN, "Velocity Time Integral")                                     |                                  |   |
|               | PVA(VTI)  | (G-038E, SRT, "Cardiovascular Orifice Area")                                |                                  |   |
|               | PV AccT   | (20168-1, LN, "Acceleration Time")  |                                  |   |
|               | PV Acc  | (99116-01, MDSN, "Pulmonic Valve Flow Acceleration")                        |                                  |   |
|               | PV DecT   | (20217-6, LN, "Deceleration Time")  |                                  |   |
|               | PV Dec  | (20216-8, LN, "Deceleration Slope")   |                                  |   |
|               | PV ET   | (18042-2, LN, "Pulmonic Valve Ejection Time")                               |                                  |   |
|               | PV AccT/ET  | (G-0388, SRT, "Ratio of Pulmonic Valve Acceleration Time to Ejection Time") |                                  |   |
| R-R Int       | (99116-02, MDSN, "Pulmonic Valve R to R Interval")            |   |                                  |   |
| Q to PV Close | (20295-2, LN, "Time from Q wave to Pulmonic Valve Closes")    |   |                                  |   |
| MPA Vmax      | (G-038A, SRT, "Main Pulmonary Artery Velocity")               |   |                                  |   |
| PR VCW        | (99116-06, MDSN, "Pulmonic Regurgitant Vena Contracta Width") |   |                                  |   |

|           |               |   |                                  |  |
|-----------|---------------|---|----------------------------------|--|
|           | PR Vmax       | (11726-7, LN, "Peak Velocity")                                      |                                  | =DCM\125215\Continuity Equation by Velocity Time Integral<br><br>Measurement Method =DCM\125210\Area by Pressure Half-Time |
|           | PR PGmax      | (20247-3, LN, "Peak Gradient")                                      |                                  |  |
|           | PRend Vmax    | (79918-9, LN, "Pulmonic regurgitation end diastolic velocity")      |                                  |  |
|           | PRend PGmax   | (79934-6, LN, "Pulmonic regurgitation end diastolic peak gradient") |                                  |  |
|           | PR Vmean      | (20352-1, LN, "Mean Velocity")                                      |                                  |  |
|           | PR PGmean     | (20256-4, LN, "Mean Gradient")                                      |                                  |  |
|           | PV PHT        | (20280-4, LN, "Pressure Half-Time")                                 |                                  |  |
|           | PR AccT       | (20168-1, LN, "Acceleration Time")                                  |                                  |  |
|           | PR Acc        | (99116-03, MDSN, "Pulmonic Regurgitant Flow Acceleration")          |                                  |  |
|           | PR DecT       | (20217-6, LN, "Deceleration Time")                                  |                                  |  |
|           | PR Dec        | (20216-8, LN, "Deceleration Slope")                                 |                                  |  |
|           | PR VTI        | (20354-7, LN, "Velocity Time Integral")                             |                                  |  |
|           | PR PISA Rad   | (99116-04, MDSN, "Pulmonic Regurgitant PISA Radius")                |                                  |  |
|           | PR Alias Vel. | (99116-05, MDSN, "Pulmonic Regurgitant Aliasing Velocity")          |                                  |  |
|           | PR Flow Rate  | (34141-2, LM, "Peak Instantaneous Flow Rate")                       |                                  |  |
|           | PR ERO        | (G-038E, SRT, "Cardiovascular Orifice Area")                        |                                  |  |
|           | PR Volume     | (33878-0, LM, "Volume Flow")  |                                  |  |
|           | PR Fraction   | (G-0390, SRT, "Regurgitant Fraction")                               |                                  |  |
|           | mPAP          | (99116-09, MDSN, "Mean Pulmonary Artery Pressure")                  |                                  |  |
|           | dPAP          | (99116-10, MDSN, "Diastolic Pulmonary Artery Pressure")             |                                  |  |
|           | PR RVSP       | (F-03DFE, SRT, "Right Ventricular Systolic Pressure")               |                                  |  |
| Tei Index | LV IVCT       | (G-037E, SRT, "Left Ventricular Isovolumic Contraction Time")       | (T-32600, SRT, "Left Ventricle") |  |
|           | LV ET         | (122211, DCM, "Left Ventricular ejection time")                     |                                  |  |

|               |               |   |  |  |
|---------------|---------------|---|--|--|
|               | LV IVRT       | (18071-1, LN, "Left Ventricular IsoVolumic Relaxation Time")              |  |  |
|               | LV MCO        | (99104-04, MDSN, "Left Ventricle Total Systolic Time")                    |  |  |
|               | LV Tei Index  | (G-037F, SRT, "Left Ventricular Index of Myocardial Performance")         |  |  |
|               | RV IVCT       | (99105-09, MDSN, "Right Ventricle IsoVolumic Contraction Time")           | (T-32500, SRT, "Right Ventricle")            |  |
|               | RV ET         | (122213, DCM, "Right Ventricular ejection time")                          |  |  |
|               | RV IVRT       | (99105-10, MDSN, "Right Ventricle IsoVolumic Relaxation Time")            |  |  |
|               | RV MCO        | (99105-07, MDSN, "Right Ventricle Total Systolic Time")                   |  |  |
|               | RV Tei Index  | (99105-08, MDSN, "Right Ventricle Myocardial Performance Index")          |  |  |
| Pulm. Veins   | PVein S Vmax  | (29450-4, LN, "Pulmonary Vein Systolic Peak Velocity")                    | (T-48581, SRT, "Pulmonary Venous Structure") |  |
|               | PVein D Vmax  | (29451-2, LN, "Pulmonary Vein Diastolic Peak Velocity")                   |  |  |
|               | P Vein S/D    | (29452-0, LN, "Pulmonary Vein Systolic to Diastolic Ratio")               |  |  |
|               | P Vein A Vmax | (29453-8, LN, "Pulmonary Vein Atrial Contraction Reversal Peak Velocity") |  |  |
|               | P Vein A Dur  | (G-038B, SRT, "Pulmonary Vein A-Wave Duration")                           |  |  |
| Hepatic Veins | H Vein S Vmax | (29471-0, LN, "Hepatic Vein Systolic Peak Velocity")                      | (T-48720, SRT, "Hepatic Vein")               |  |
|               | H Vein D Vmax | (29472-8, LN, "Hepatic Vein Diastolic Peak Velocity")                     |  |  |
|               | H Vein S/D    | (29473-6, LN, "Hepatic Vein Systolic to Diastolic Ratio")                 |  |  |
|               | H Vein A Vmax | (29474-4, LN, "Hepatic Vein Atrial Contraction Reversal Peak Velocity")   |  |  |
|               | H Vein A Dur  | (99112-01, MDSN, "Hepatic Vein A-Wave Duration")                          |  |  |

|                 |   |                                  |
|-----------------|---|----------------------------------|
| Peak E'         | (G-037A, SRT, "Left Ventricular Peak Early Diastolic Tissue Velocity")              | (T-32600, SRT, "Left Ventricle") |
| Septal E'       | (99104-17, MDSN, "Left Ventricular Septal Peak Early Diastolic Tissue Velocity")    |                                  |
| Lateral E'      | (99104-18, MDSN, "Left Ventricular Lateral Peak Early Diastolic Tissue Velocity")   |                                  |
| Avg E'          | (99104-19, MDSN, "Left Ventricular Average Peak Early Diastolic Tissue Velocity")   |                                  |
| Peak A'         | (G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")            |                                  |
| Septal A'       | (99104-20, MDSN, "LV Septal Peak Diastolic Tissue Velocity During Atrial Systole")  |                                  |
| Lateral A'      | (99104-21, MDSN, "LV Lateral Peak Diastolic Tissue Velocity During Atrial Systole") |                                  |
| Peak S          | (G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity")                     |                                  |
| Septal S'       | (99104-22, MDSN, "Left Ventricular Septal Peak Systolic Tissue Velocity")           |                                  |
| Lateral S'      | (99104-23, MDSN, "Left Ventricular Lateral Peak Systolic Tissue Velocity")          |                                  |
| MV E/E'         | (G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")        |                                  |
| MV E/Septal E'  | (99104-26, MDSN, "MV E to LV Septal Peak A' Ratio")                                 |                                  |
| MV E/Lateral E' | (99104-27, MDSN, "MV E to LV Lateral Peak A' Ratio")                                |                                  |
| MV E/Avg E'     | (99104-28, MDSN, "MV E to LV Average Peak A' Ratio")                                |                                  |

|       |               |  |                                   |                                  |  |
|-------|---------------|--|-----------------------------------|----------------------------------|--|
|       | E'/A'         | (99104-06, MDSN, "LV Peak E' to LV Peak A' Ratio")                                       |                                   |                                  |  |
|       | Septal E'/A'  | (99104-24, MDSN, "LV Septal Peak E' to LV Peak A' Ratio")                                |                                   |                                  |  |
|       | Lateral E'/A' | (99104-25, MDSN, "LV Lateral Peak E' to LV Peak A' Ratio")                               |                                   |                                  |  |
|       | PCWP          | (118433006, SCT, "Pulmonary Capillary Wedge Pressure")                                   |                                   |                                  |  |
|       | LV AccT       | (99104-07, MDSN, "Left Ventricle Acceleration Time")                                     |                                   |                                  |  |
|       | LV DecT       | (99104-08, MDSN, "Left Ventricle Deceleration Time")                                     |                                   |                                  |  |
|       | RV Peak E'    | (99105-11, MDSN, "Right Ventricle Peak Early Diastolic Tissue Velocity")                 | (T-32500, SRT, "Right Ventricle") |                                  |  |
|       | RV Peak A'    | (99105-12, MDSN, "Right Ventricle Peak Diastolic Tissue Velocity During Atrial Systole") |                                   |                                  |  |
|       | RV Peak S'    | (99105-13, MDSN, "Right Ventricle Peak Systolic Tissue Velocity")                        |                                   |                                  |  |
|       | TV E/E'       | (99105-14, MDSN, "TV Peak E to RV Peak E' Ratio")  |                                   |                                  |  |
|       | RV E'/A'      | (99105-15, MDSN, "RV Peak E' to RV Peak A' Ratio")                                       |                                   |                                  |  |
|       | RV AccT       | (20168-1, LN, "Acceleration Time")   |                                   |                                  |  |
|       | RV DecT       | (20217-6, LN, "Deceleration Time")   |                                   |                                  |  |
| Qp/Qs | Systemic VTl  | (20354-7, LN, "Velocity Time Integral")  |                                   | (T-32600, SRT, "Left Ventricle") |  |
|       | Pulm. VTl(P)  | (20354-7, LN, "Velocity Time Integral")  |                                   |                                  |  |
|       | Systemic SV   | (F-32120, SRT, "Stroke Volume")  |                                   |                                  |  |
|       | Systemic SI   | (F-00078, SRT, "Stroke Index")   |                                   |                                  |  |
|       | Systemic CO   | (F-32100, SRT, "Cardiac Output")   |                                   |                                  |  |
|       | Systemic CI   | (F-32110, SRT, "Cardiac Index")  |                                   |                                  |  |
|       | Pulmonic VTl  | (20354-7, LN, "Velocity Time Integral")  | (T-32500, SRT, "Right Ventricle") |                                  |  |
|       | Pulmonic SV   | (F-32120, SRT, "Stroke Volume")  |                                   |                                  |  |
|       | Pulmonic SI   | (F-00078, SRT, "Stroke Index")   |                                   |                                  |  |

|             |  |  |
|-------------|--|--|
| Pulmonic CO | (F-32100, SRT, "Cardiac Output")             |  |
| Pulmonic CI | (F-32110, SRT, "Cardiac Index")              |  |
| Qp/Qs       | (29462-9, LN, "Pulmonary-to-Systemic Shunt") |  |

## 9.4 UROLOGY STRUCTURED REPORT TEMPLATE

### 9.4.1 Urology Ultrasound Report Templates (TID SM99400)

**Table 9.4-1  
Urology Ultrasound Report Procedure Templates**

| No | Rel With Parent | VT        | Concept Name                                       | Comments  | V7 Label         |
|----|-----------------|-----------|--|---|------------------|
| 1  |                 | CONTAINER | EV (99400, MDSN, "Urology Ultrasound Report")      |   |                  |
| 2  | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observation Context                    |   |                  |
| 3  | CONTAINS        | INCLUDE   | DTID (5001) Patient Characteristics                |   |                  |
| 4  | CONTAINS        | INCLUDE   | DTID (5016) LWH Volume Group                       | \$GroupName = EV (99017-0, MDSN, " WG Prostate Volume")           | WG Prostate Vol. |
| 5  | CONTAINS        | INCLUDE   | DTID (5016) LWH Volume Group                       | \$GroupName = EV (99017-1, MDSN, " T-Zone Volume")                | T-Zone Vol.      |
| 6  | CONTAINS        | INCLUDE   | DTID (5016) LWH Volume Group                       | \$GroupName = EV (99017-2, MDSN, " Bladder Volume")               | Bladder Vol.     |
| 7  | CONTAINS        | INCLUDE   | DTID (SM99401) Residual Volume Group               | \$GroupName = EV (99017-3, MDSN, " Residual Volume")              | Residual Vol.    |
| 8  | CONTAINS        | INCLUDE   | DTID (5016) LWH Volume Group                       | \$GroupName = EV (99017-6, MDSN, " Right Renal Volume")           | Rt. Renal Vol.   |
| 9  | CONTAINS        | INCLUDE   | DTID (5016) LWH Volume Group                       | \$GroupName = EV (99017-7, MDSN, " Left Renal Volume")            | Lt. Renal Vol.   |
| 10 | CONTAINS        | INCLUDE   | DTID (SM99402) PSA Measurement Group               | \$GroupName = EV (99017-9, MDSN, "Prostate Specific Antigen")     | PSA              |
| 11 | CONTAINS        | INCLUDE   | DTID (5104) Vascular Measurement Group             | \$AnatomyGroup = EV (99017-8, MDSN, "General Urology Measurement" | General          |
| 12 | CONTAINS        | INCLUDE   | DTID (SM99410) Urology User Creation Group Section | Ref. Section 9.4.1.7  |                  |

### 9.4.1.1 Observation Context (TID 1001)

**Table 9.4-2  
OBSERVATION CONTEXT IN UROLOGY SR**

|     | REL             | VT    | Concept Name                          | Unit / CODE Value   | V7 Label                 |
|-----|-----------------|-------|---------------------------------------|---|--------------------------|
| D-1 | HAS OBS CONTEXT | CODE  | (121005, DCM, "Observer Type")        | (121006, DCM, "Person")   |                          |
| D-2 | HAS OBS CONTEXT | PNAME | (121008, DCM, "Person Observer Name") |   | Ref. Physician           |
| D-3 | HAS OBS CONTEXT | CODE  | (121024, DCM, "Subject Class")        | (121025 ,DCM,"Patient")   |                          |
| D-4 | HAS OBS CONTEXT | PNAME | (121029,DCM, "Subject Name")          |   | Last Name,<br>First Name |
| D-5 | HAS OBS CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")    |   | BirthDate                |
| D-6 | HAS OBS CONTEXT | CODE  | (121032,DCM, "Subject Sex")           | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown sex") | Gender                   |

### 9.4.1.2 Patient Characteristics (TID 5001)

**Table 9.4-3  
PATIENT CHARACTERISTICS IN UROLOGY SR**

|       | REL      | VT        | Concept Name                           | Unit / CODE Value                                    | V7 Label |
|-------|----------|-----------|--|--|----------|
| D-7   | CONTAINS | CONTAINER | (121118,DCM "Patient Characteristics") |  |          |
| D-7-1 | CONTAINS | NUM       | (8302-2, LN, "Patient Height")         | (cm, UCUM, "centimeter")<br>(mm, UCUM, "millimeter") | Height   |
| D-7-2 | CONTAINS | NUM       | (29463-7, LN, "Patient Weight")        | (kg, UCUM, "kilograms")                              | Weight   |

### 9.4.1.3 LHW Volume Group (TID 5016)

**Table 9.4-4  
LHW VOLUME GROUP IN UROLOGY SR**

|       | REL      | VT        | Concept Name           | Unit / CODE Value                              | V7 Label |
|-------|----------|-----------|------------------------|--|----------|
| D-8   |          | CONTAINER | \$GroupName            |  |          |
| D-8-1 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Volume                       | Vol.     |
| D-8-2 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Length                       | Length   |
|       |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |
| D-8-3 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Width                        | Width    |
|       |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |
| D-8-4 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Height                       | Height   |
|       |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |

#### 9.4.1.4 Residual Volume Group (TID SM99401)

**Table 9.4-5  
Residual Volume Group in UROLOGY SR**

|       | REL      | VT        | Concept Name                 | Unit / CODE Value                                       | V7 Label |
|-------|----------|-----------|------------------------------|---|----------|
| D-9   |          | CONTAINER | \$GroupName                  |   |          |
| D-9-1 | CONTAINS | INCLUDE   | DTID (5016) LWH Volume Group | \$GroupName = EV (99017-4,<br>MDSN, "Pre Void Volume")  | Pre      |
| D-9-2 | CONTAINS | INCLUDE   | DTID (5016) LWH Volume Group | \$GroupName = EV (99017-5,<br>MDSN, "Post Void Volume") | Post     |
| D-9-3 | CONTAINS | INCLUDE   | DTID (5016) LWH Volume Group | \$GroupName = EV (99017-12,<br>MDSN, "Void Volume")     | Void     |

#### 9.4.1.5 PSA Measurement Group (TID SM99402)

**Table 9.4-6  
PSA Measurement Group in UROLOGY SR**

|        | REL      | VT        | Concept Name                            | Unit / CODE Value | V7 Label |
|--------|----------|-----------|---|-------------------|----------|
| D-10   |          | CONTAINER | \$GroupName                             |                   |          |
| D-10-1 | CONTAINS | NUM       | (99017-10, MDSN, "Predicted PSA by WG") |                   | PREDPSA  |

|        |          |     |   |  |       |
|--------|----------|-----|---|--|-------|
| D-10-2 | CONTAINS | NUM | (99017-11, MDSN, "Predicted PSA by T-Zone") |  | TZPSA |
|--------|----------|-----|---|--|-------|

#### 9.4.1.6 Vascular Measurement Group (TID 5104)

**Table 9.4-7**  
**Vascular Measurement Group in UROLOGY SR**

|        | REL      | VT        | Concept Name           | Unit / CODE Value  | V7 Label |
|--------|----------|-----------|------------------------|--|----------|
| D-11   |          | CONTAINER | \$AnatomyGroup         | \$AnatomyGroup=MDSN\99017-8\General Urology Measurement      | General  |
| D-11-1 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = DCID (12119)<br>Vascular Ultrasound Property |          |
|        |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type               |          |

#### 9.4.1.7 Urology User Creation Group Section (TID SM99410)

**Table 9.4-8**  
**User Creation Group in UROLOGY SR**

|          | REL             | VT        | Concept Name  | Unit / CODE Value     | V7 Label | Comments   |
|----------|-----------------|-----------|---|-----------------------|----------|--|
| D-12     | CONTAINS        | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |                       |          | *Creation ID: Randomly generated 7-digit unique ID |
| D-12-1   | CONTAINS        | NUM       | (99900-Creation ID, MDSN, "User Creation Item Name")  |                       |          |  |
| D-12-1-1 | HAS CONCEPT MOD | CODE      | (121401, DCM, "Derivation")                           | Common CID-Derivation |          |  |

#### 9.4.2 Urology Measurement and Calculation used in Urology SR

**Table 9.4-9**  
**Urology Measurement and Calculation Items**

| <b>V7 Label</b>    | <b>DICOM SR Concept Name (CSD\CV\CM)</b> |
|--------------------|--|
| WG Prostate Length | SRT\G-A22A\Length                        |
| WG Prostate Height | DCM\121207\Height                        |
| WG Prostate Width  | SRT\G-A220\Width                         |
| Vol.               | SRT\G-D705\Volume                        |
| T-Zone Length      | SRT\G-A22A\Length                        |
| T-Zone Height      | DCM\121207\Height                        |
| T-Zone Width       | SRT\G-A220\Width                         |
| Vol.               | SRT\G-D705\Volume                        |
| Bladder Length     | SRT\G-A22A\Length                        |
| Bladder Height     | DCM\121207\Height                        |
| Bladder Width      | SRT\G-A220\Width                         |
| Vol.               | SRT\G-D705\Volume                        |
| Pre L              | SRT\G-A22A\Length                        |
| Pre H              | DCM\121207\Height                        |
| Pre W              | SRT\G-A220\Width                         |
| Pre Vol.           | SRT\G-D705\Volume                        |
| Post L             | SRT\G-A22A\Length                        |
| Post H             | DCM\121207\Height                        |
| Post W             | SRT\G-A220\Width                         |
| Post Vol.          | SRT\G-D705\Volume                        |
| Void Vol.          | SRT\G-D705\Volume                        |
| Rt. Renal Length   | SRT\G-A22A\Length                        |
| Rt. Renal Height   | DCM\121207\Height                        |
| Rt. Renal Width    | SRT\G-A220\Width                         |
| Vol.               | SRT\G-D705\Volume                        |
| Rt. Renal Pelvis   | MDSN\99005-13\Right Pelvis               |
| Lt. Renal Length   | SRT\G-A22A\Length                        |
| Lt. Renal Height   | DCM\121207\Height                        |
| Lt. Renal Width    | SRT\G-A220\Width                         |
| Vol.               | SRT\G-D705\Volume                        |
| Lt. Renal Pelvis   | MDSN\99005-14\Left Pelvis                |
| PSA                | MDSN\99017-9\Prostate Specific Antigen   |
| PREDPSA            | MDSN\99017-10\Predicted PSA by WG        |
| TZPSA              | MDSN\99017-11\Predicted PSA by T-Zone    |

## 9.5 SMALL PARTS STRUCTURED REPORT TEMPLATE

### 9.5.1 SmallParts Ultrasound Report Templates (TID SM99500)

**Table 9.5-1  
SmallParts Ultrasound Report Procedure Templates**

| No | Rel With Parent | VT        | Concept Name                                       | Comments | V7 Label |
|----|-----------------|-----------|--|----------|----------|
| 1  |                 | CONTAINER | EV (SM99500, MDSN, "SmallParts Ultrasound Report") |          |          |
| 2  | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observation Context                    |          |          |
| 3  | CONTAINS        | INCLUDE   | DTID (5001) Patient Characteristics                |          |          |
| 4  | CONTAINS        | INCLUDE   | DTID (SM99501) Small Parts Measurement Group       |          |          |
| 8  | CONTAINS        | INCLUDE   | DTID (SM99501) Small Parts Measurement Group       |          |          |
| 9  | CONTAINS        | INCLUDE   | DTID (SM99501) Small Parts Measurement Group       |          |          |

#### 9.5.1.1 Observation Context (TID 1001)

**Table 9.5-2  
OBSERVATION CONTEXT IN UROLOGY SR**

|     | REL             | VT    | Concept Name                          | Unit / CODE Value       | V7 Label                 |
|-----|-----------------|-------|---------------------------------------|-------------------------|--------------------------|
| E-1 | HAS OBS CONTEXT | CODE  | (121005, DCM, "Observer Type")        | (121006, DCM, "Person") |                          |
| E-2 | HAS OBS CONTEXT | PNAME | (121008, DCM, "Person Observer Name") |                         | Ref. Physician           |
| E-3 | HAS OBS CONTEXT | CODE  | (121024, DCM, "Subject Class")        | (121025 ,DCM,"Patient") |                          |
| E-4 | HAS OBS CONTEXT | PNAME | (121029,DCM, "Subject Name")          |                         | Last Name,<br>First Name |
| E-5 | HAS OBS CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")    |                         | BirthDate                |

|     |                 |      |                             |   |        |
|-----|-----------------|------|-----------------------------|---|--------|
| E-6 | HAS OBS CONTEXT | CODE | (121032,DCM, "Subject Sex") | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown sex") | Gender |
|-----|-----------------|------|-----------------------------|---|--------|

### 9.5.1.2 Patient Characteristics (TID 5001)

**Table 9.5-3  
PATIENT CHARACTERISTICS IN SMALL PARTS SR**

|       | REL      | VT        | Concept Name                           | Unit / CODE Value                                    | V7 Label |
|-------|----------|-----------|--|--|----------|
| E-7   | CONTAINS | CONTAINER | (121118,DCM "Patient Characteristics") |  |          |
| E-7-1 | CONTAINS | NUM       | (8302-2, LN, "Patient Height")         | (cm, UCUM, "centimeter")<br>(mm, UCUM, "millimeter") | Height   |
| E-7-2 | CONTAINS | NUM       | (29463-7, LN, "Patient Weight")        | (kg, UCUM, "kilograms")                              | Weight   |

### 9.5.1.3 Small Parts Measurement Group (TID SM99501)

**Table 9.5-4  
Small Parts Measurement Group in SMALL PARTS SR**

|       | REL                   | VT        | Concept Name                          | Unit / CODE Value  | V7 Label |
|-------|-----------------------|-----------|---------------------------------------|--|----------|
| E-8   |                       | CONTAINER | (121070, DCM, "Findings")             |  |          |
| E-8-1 | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding Site")         | (T-B6000, SRT, "Thyroid")  |          |
|       |                       |           |                                       | (T-94000, SRT, "Testis")   |          |
|       |                       |           |                                       | (G-A139, SRT, "Superficial")   |          |
| E-8-2 | HAS<br>CONCEPT<br>MOD | CODE      | (G-C171, SRT, "Laterality")           | (G-A100, SRT, "Right")<br>(G-A101, SRT, "Left")<br>(G-A103, SRT, "Unilateral") |          |
| E-8-3 | CONTAINS              | INCLUDE   | DTID (SM99102) Mass Measurement Group |  | Mass     |
| E-8-4 | CONTAINS              | INCLUDE   | DTID (5016) LWH Volume Group          | \$GroupName = EV (99100-6, MDSN, "Thyroid Volume")                             | Vol.     |
|       |                       |           |                                       | \$GroupName = EV (99100-9, MDSN, "Testis Volume")                              |          |

|         |          |           |  |   |         |
|---------|----------|-----------|--|---|---------|
|         |          |           |  | \$GroupName = EV (99100-11, MDSN, "Superficial Volume")   |         |
| E-8-5   | CONTAINS | INCLUDE   | DTID (5104) Vascular Measurement Group                 | \$AnatomyGroup = EV (99100-7, MDSN, "Thyroid Flow")<br>\$AnatomyGroup = EV(99100-10, MDSN, "Testis Flow")<br>\$AnatomyGroup = EV (99100-12, MDSN, "Superficial Flow") | Flow    |
| E-8-6   | CONTAINS | CONTAINER | (99100-16, MDSN, "Isthmus")                            |   | Isthmus |
| E-8-6-1 | CONTAINS | NUM       | (99100-16, MDSN, "Isthmus")                            |   |         |
| E-8-7   | CONTAINS | INCLUDE   | DTID (SM99510) Small Parts User Creation Group Section | Ref. Section 9.5.1.7  |         |

#### 9.5.1.4 Mass Measurement Group (TID SM99102)

**Table 9.5-5  
Mass Measurement Group in SMALL PARTS SR**

|         | REL             | VT        | Concept Name               | Unit / CODE Value  | V7 Label            |
|---------|-----------------|-----------|----------------------------|--|---------------------|
| E-8-2   | CONTAINS        | CONTAINER | (M-03000, SRT, "Mass")     |  |                     |
| E-8-2-1 | HAS OBS CONTEXT | TEXT      | (12510, DCM, "Identifier") |  | "1", "2" ...<br>"5" |
| E-8-2-2 | CONTAINS        | INCLUDE   | DTID (300) Measurement     | \$Measurement = \$Volume   | Vol.                |
| E-8-2-3 | CONTAINS        | INCLUDE   | DTID (300) Measurement     | \$Measurement = \$Length<br>\$Derivation = DCID (3627)<br>Measurement Type | L                   |
| E-8-2-4 | CONTAINS        | INCLUDE   | DTID (300) Measurement     | \$Measurement = \$Depth<br>\$Derivation = DCID (3627)<br>Measurement Type  | D                   |
| E-8-2-5 | CONTAINS        | INCLUDE   | DTID (300) Measurement     | \$Measurement = \$Width<br>\$Derivation = DCID (3627)<br>Measurement Type  | W                   |

### 9.5.1.5 LHW Volume Group (TID 5016)

**Table 9.5-6  
LHW VOLUME GROUP IN SMALL PARTS SR**

|         | REL      | VT        | Concept Name           | Unit / CODE Value                              | V7 Label |
|---------|----------|-----------|------------------------|--|----------|
| E-8-3   |          | CONTAINER | \$GroupName            |  |          |
| E-8-3-1 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Volume                       | Vol.     |
| E-8-3-2 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Length                       | L        |
|         |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |
| E-8-3-3 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Width                        | W        |
|         |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |
| E-8-3-4 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = \$Height                       | H        |
|         |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type |          |

### 9.5.1.6 Vascular Measurement Group (TID 5104)

**Table 9.5-7  
Vascular Measurement Group in SMALL PARTS SR**

|         | REL      | VT        | Concept Name           | Unit / CODE Value  | V7 Label |
|---------|----------|-----------|------------------------|--|----------|
| E-8-4   |          | CONTAINER | \$AnatomyGroup         |  | Flow     |
| E-8-4-1 | CONTAINS | INCLUDE   | DTID (300) Measurement | \$Measurement = DCID (12119)<br>Vascular Ultrasound Property |          |
|         |          |           |                        | \$Derivation = DCID (3627)<br>Measurement Type               |          |

### 9.5.1.7 Small Parts User Creation Group Section (TID SM99510)

**Table 9.5-8  
User Creation Group in SMALL PARTS SR**

|  | REL | VT | Concept Name | Unit / CODE | V7 | Comments |
|--|-----|----|--------------|-------------|----|----------|
|--|-----|----|--------------|-------------|----|----------|

|           |                    |           |   | Value                     | Label |   |
|-----------|--------------------|-----------|---|---------------------------|-------|---|
| E-8-7     | CONTAINS           | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |                           |       | *Creation ID:<br>Randomly<br>generated 7-digit<br>unique ID |
| E-8-7-1   | CONTAINS           | NUM       | (99900-Creation ID, MDSN, "User Creation Item Name")  |                           |       |   |
| E-8-7-1-1 | HAS CONCEPT<br>MOD | CODE      | (121401, DCM, "Derivation")                           | Common CID-<br>Derivation |       |   |

## 9.5.2 SmallParts Measurement and Calculation used in Small Parts SR

### 9.5.2.1 Thyroid Measurement and Calculation used in Small Parts SR

**Table 9.5-9  
Thyroid Measurement and Calculation Items**

| V7 Label   | DICOM SR Concept Name (CDS CV CM) |
|------------|-----------------------------------|
| Thyroid    | SRT\T-B6000\Thyroid               |
| Mass1      | SRT\M-03000\Mass                  |
| Mass1 L    | SRT\G-A22A\Length                 |
| Mass1 D    | DCM\G-D785\Depth                  |
| Mass1 W    | SRT\G-A220\Width                  |
| Mass1 Vol. | SRT\G-D705\Volume                 |
| Mass2      | SRT\M-03000\Mass                  |
| Mass2 L    | SRT\G-A22A\Length                 |
| Mass2 D    | DCM\G-D785\Depth                  |
| Mass2 W    | SRT\G-A220\Width                  |
| Mass2 Vol. | SRT\G-D705\Volume                 |
| Mass3      | SRT\M-03000\Mass                  |
| Mass3 L    | SRT\G-A22A\Length                 |
| Mass3 D    | DCM\G-D785\Depth                  |
| Mass3 W    | SRT\G-A220\Width                  |
| Mass3 Vol. | SRT\G-D705\Volume                 |
| Mass4      | SRT\M-03000\Mass                  |

|                  |  |
|------------------|--|
| Mass4 L          | SRT\G-A22A\Length                                  |
| Mass4 D          | DCM\G-D785\Depth                                   |
| Mass4 W          | SRT\G-A220\Width                                   |
| Mass4 Vol.       | SRT\G-D705\Volume                                  |
| Mass5            | SRT\M-03000\Mass                                   |
| Mass5 L          | SRT\G-A22A\Length                                  |
| Mass5 D          | DCM\G-D785\Depth                                   |
| Mass5 W          | SRT\G-A220\Width                                   |
| Mass5 Vol.       | SRT\G-D705\Volume                                  |
| Thyroid Vol.     | MDSN\99100-6\Thyroid Volume                        |
| Thyroid L        | SRT\G-A22A\Length                                  |
| Thyroid H        | DCM\121207\Height                                  |
| Thyroid W        | SRT\G-A220\Width                                   |
| Thyroid Flow     | MDSN\99100-7\Thyroid Flow                          |
| PS               | LN\11726-7\Peak Systolic Velocity                  |
| ED               | LN\11653-3\End Diastolic Velocity                  |
| MD               | LN\11665-7\Minimum Diastolic Velocity              |
| TAmx             | LN\11692-1\Time averaged peak velocity             |
| TAmx             | LN\20352-1\Time averaged mean velocity             |
| PGmax            | LN\20247-3\Peak Gradient                           |
| PGmean           | LN\20256-4\Mean Gradient                           |
| S/D              | LN\12144-2\Systolic to Diastolic Velocity Ratio    |
| D/S              | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |
| RI               | LN\12023-8\Resistivity Index                       |
| PI               | LN\12008-9\Pulsatility Index                       |
| %StA             | SRT\R-101BA\Lumen Area Stenosis                    |
| %StA Outer Area  | SRT\G-0366\Vessel lumen cross-sectional area       |
| %StA Inner Area  | SRT\R-1025D\Vessel Intimal Cross-Sectional Area    |
| %StD             | SRT\R-101BB\Lumen Diameter Stenosis                |
| %StD Outer Dist. | SRT\G-0364\Vessel lumen diameter                   |
| %StD Inner Dist. | SRT\R-1025C\Vessel Intimal Diameter                |
| Vesl. Area       | MDSN\99200-02\Vessel Area                          |
| Vol. Flow(A)     | SRT\G-0365\Vessel outside diameter                 |
| Vesl. Dist       | LN\33878-0\Volume flow                             |
| Vol. Flow(D)     | LN\33878-0\Volume flow                             |

## 9.5.2.2 Testis Measurement and Calculation used in Small Parts SR

**Table 9.5-9  
Testis Measurement and Calculation Items**

| V7 Label    | DICOM SR Concept Name (CDS CV CM) |
|-------------|-----------------------------------|
| Testis      | SRT\T-94000\Testis                |
| Mass1       | SRT\M-03000\Mass                  |
| Mass1 L     | SRT\G-A22A\Length                 |
| Mass1 D     | DCM\G-D785\Depth                  |
| Mass1 W     | SRT\G-A220\Width                  |
| Mass1 Vol.  | SRT\G-D705\Volume                 |
| Mass2       | SRT\M-03000\Mass                  |
| Mass2 L     | SRT\G-A22A\Length                 |
| Mass2 D     | DCM\G-D785\Depth                  |
| Mass2 W     | SRT\G-A220\Width                  |
| Mass2 Vol.  | SRT\G-D705\Volume                 |
| Mass3       | SRT\M-03000\Mass                  |
| Mass3 L     | SRT\G-A22A\Length                 |
| Mass3 D     | DCM\G-D785\Depth                  |
| Mass3 W     | SRT\G-A220\Width                  |
| Mass3 Vol.  | SRT\G-D705\Volume                 |
| Mass4       | SRT\M-03000\Mass                  |
| Mass4 L     | SRT\G-A22A\Length                 |
| Mass4 D     | DCM\G-D785\Depth                  |
| Mass4 W     | SRT\G-A220\Width                  |
| Mass4 Vol.  | SRT\G-D705\Volume                 |
| Mass5       | SRT\M-03000\Mass                  |
| Mass5 L     | SRT\G-A22A\Length                 |
| Mass5 D     | DCM\G-D785\Depth                  |
| Mass5 W     | SRT\G-A220\Width                  |
| Mass5 Vol.  | SRT\G-D705\Volume                 |
| Testis Vol. | MDSN\99100-9\Testis Volume        |
| Testis L    | SRT\G-A22A\Length                 |

|                  |  |
|------------------|--|
| Testis H         | DCM\121207\Height                                  |
| Testis W         | SRT\G-A220\Width                                   |
| Testis Flow      | MDSN\99100-10\Testis Flow                          |
| PS               | LN\11726-7\Peak Systolic Velocity                  |
| ED               | LN\11653-3\End Diastolic Velocity                  |
| MD               | LN\11665-7\Minimum Diastolic Velocity              |
| TAmx             | LN\11692-1\Time averaged peak velocity             |
| TAmx             | LN\20352-1\Time averaged mean velocity             |
| PGmax            | LN\20247-3\Peak Gradient                           |
| PGmean           | LN\20256-4\Mean Gradient                           |
| S/D              | LN\12144-2\Systolic to Diastolic Velocity Ratio    |
| D/S              | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |
| RI               | LN\12023-8\Resistivity Index                       |
| PI               | LN\12008-9\Pulsatility Index                       |
| %StA             | SRT\R-101BA\Lumen Area Stenosis                    |
| %StA Outer Area  | SRT\G-0366\Vessel lumen cross-sectional area       |
| %StA Inner Area  | SRT\R-1025D\Vessel Intimal Cross-Sectional Area    |
| %StD             | SRT\R-101BB\Lumen Diameter Stenosis                |
| %StD Outer Dist. | SRT\G-0364\Vessel lumen diameter                   |
| %StD Inner Dist. | SRT\R-1025C\Vessel Intimal Diameter                |
| Vesl. Area       | MDSN\99200-02\Vessel Area                          |
| Vol. Flow(A)     | SRT\G-0365\Vessel outside diameter                 |
| Vesl. Dist       | LN\33878-0\Volume flow                             |
| Vol. Flow(D)     | LN\33878-0\Volume flow                             |

### 9.5.2.3 Superficial Measurement and Calculation used in Small Parts SR

**Table 9.5-10  
Superficial Measurement and Calculation Items**

| V7 Label    | DICOM SR Concept Name (CDS CV CM) |
|-------------|-----------------------------------|
| Superficial | SRT\G-A139\Superficial            |
| Mass1       | SRT\M-03000\Mass                  |
| Mass1 L     | SRT\G-A22A\Length                 |
| Mass1 D     | DCM\G-D785\Depth                  |

|                  |   |
|------------------|---|
| Mass1 W          | SRT\G-A220\Width                                |
| Mass1 Vol.       | SRT\G-D705\Volume                               |
| Mass2            | SRT\M-03000\Mass                                |
| Mass2 L          | SRT\G-A22A\Length                               |
| Mass2 D          | DCM\G-D785\Depth                                |
| Mass2 W          | SRT\G-A220\Width                                |
| Mass2 Vol.       | SRT\G-D705\Volume                               |
| Mass3            | SRT\M-03000\Mass                                |
| Mass3 L          | SRT\G-A22A\Length                               |
| Mass3 D          | DCM\G-D785\Depth                                |
| Mass3 W          | SRT\G-A220\Width                                |
| Mass3 Vol.       | SRT\G-D705\Volume                               |
| Mass4            | SRT\M-03000\Mass                                |
| Mass4 L          | SRT\G-A22A\Length                               |
| Mass4 D          | DCM\G-D785\Depth                                |
| Mass4 W          | SRT\G-A220\Width                                |
| Mass4 Vol.       | SRT\G-D705\Volume                               |
| Mass5            | SRT\M-03000\Mass                                |
| Mass5 L          | SRT\G-A22A\Length                               |
| Mass5 D          | DCM\G-D785\Depth                                |
| Mass5 W          | SRT\G-A220\Width                                |
| Mass5 Vol.       | SRT\G-D705\Volume                               |
| Superficial Vol. | MDSN\99100-11\Superficial Volume                |
| Superficial L    | SRT\G-A22A\Length                               |
| Superficial H    | DCM\121207\Height                               |
| Superficial W    | SRT\G-A220\Width                                |
| Superficial Flow | MDSN\99100-12\Superficial Flow                  |
| PS               | LN\11726-7\Peak Systolic Velocity               |
| ED               | LN\11653-3\End Diastolic Velocity               |
| MD               | LN\11665-7\Minimum Diastolic Velocity           |
| TAmx             | LN\11692-1\Time averaged peak velocity          |
| TAmx             | LN\20352-1\Time averaged mean velocity          |
| PGmax            | LN\20247-3\Peak Gradient                        |
| PGmean           | LN\20256-4\Mean Gradient                        |
| S/D              | LN\12144-2\Systolic to Diastolic Velocity Ratio |

|                  |  |
|------------------|--|
| D/S              | MDSN\99200-01\Diastolic to Systolic Velocity Ratio |
| RI               | LN\12023-8\Resistivity Index                       |
| PI               | LN\12008-9\Pulsatility Index                       |
| %StA             | SRT\R-101BA\Lumen Area Stenosis                    |
| %StA Outer Area  | SRT\G-0366\Vessel lumen cross-sectional area       |
| %StA Inner Area  | SRT\R-1025D\Vessel Intimal Cross-Sectional Area    |
| %StD             | SRT\R-101BB\Lumen Diameter Stenosis                |
| %StD Outer Dist. | SRT\G-0364\Vessel lumen diameter                   |
| %StD Inner Dist. | SRT\R-1025C\Vessel Intimal Diameter                |
| Vesl. Area       | MDSN\99200-02\Vessel Area                          |
| Vol. Flow(A)     | SRT\G-0365\Vessel outside diameter                 |
| Vesl. Dist       | LN\33878-0\Volume flow                             |
| Vol. Flow(D)     | LN\33878-0\Volume flow                             |

## 9.6 BREAST IMAGING STRUCTURED REPORT TEMPLATE

### 9.6.1 Breast Imaging Report (TID 4200)

**Table 9.6-1  
Breast Imaging Report**

| No | Rel With Parent | VT        | Concept Name   | Comments | V7 Label |
|----|-----------------|-----------|--|----------|----------|
| 1  |                 | CONTAINER | (111400, DCM, "Breast Imaging Report")               |          |          |
| 2  | HAS CONCEPT MOD | INCLUDE   | DTID (1204) Language of Content Item and Descendants |          |          |
| 3  | CONTAINS        | INCLUDE   | DTID (4202) Breast Imaging Report Narrative          |          |          |
| 4  | CONTAINS        | INCLUDE   | DTID (4208) Breast Imaging Report Supplementary Data |          |          |

#### 9.6.1.1 Language of Content Item and Descendants (TID 1204)

**Table 9.6-2  
Breast Imaging Procedure Reported**

| No  | REL             | VT   | Concept Name  | Unit / CODE Value        | V7 Label |
|-----|-----------------|------|---|--------------------------|----------|
| F-1 | HAS CONCEPT MOD | CODE | (121049,DCM,"Language of Content Item and Descendants") | (en, RFC3066, "English") |          |

#### 9.6.1.2 Breast Imaging Report Narrative (TID 4202)

**Table 9.6-3  
Breast Imaging Report Report Narrative**

|     | REL      | VT        | Concept Name                       | Unit / CODE Value | V7 Label |
|-----|----------|-----------|------------------------------------|-------------------|----------|
| F-2 | CONTAINS | CONTAINER | (111412, DCM, "Narrative Summary") |                   |          |

|         |                    |           |                                     |                                |  |
|---------|--------------------|-----------|-------------------------------------|--------------------------------|--|
| F-2-1   | CONTAINS           | CONTAINER | (121058, DCM, "Procedure reported") |                                |  |
| F-2-1-1 | HAS OBS<br>CONTEXT | INCLUDE   | DTID (1002) Observer Context        |                                |  |
| F-2-1-7 | CONTAINS           | TEXT      | (121058, DCM, "Procedure reported") | "Ultrasonography of<br>breast" |  |

### 9.6.1.3 Observer Context (TID 1002)

**Table 9.6-4  
OBSERVER CONTEXT IN BREAST SR**

|         | REL                | VT    | Concept Name                             | Unit / CODE Value   | V7 Label               |
|---------|--------------------|-------|--|---|------------------------|
| F-2-1-1 | HAS OBS<br>CONTEXT | CODE  | (121005, DCM, "Observer Type")           | (121006, DCM, "Person")   |                        |
| F-2-1-2 | HAS OBS<br>CONTEXT | PNAME | (121008, DCM, "Person Observer<br>Name") |   | Ref.<br>Physician      |
| F-2-1-3 | HAS OBS<br>CONTEXT | CODE  | (121024, DCM, "Subject Class")           | (121025 ,DCM,"Patient")   |                        |
| F-2-1-4 | HAS OBS<br>CONTEXT | PNAME | (121029,DCM, "Subject Name")             |   | LastName,<br>FirstName |
| F-2-1-5 | HAS OBS<br>CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")       |   | BirthDate              |
| F-2-1-6 | HAS OBS<br>CONTEXT | CODE  | (121032,DCM, "Subject Sex")              | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown sex") | Gender                 |

### 9.6.1.4 Breast Imaging Report Supplementary Data (TID 4208)

**Table 9.6-5  
Breast Imaging Report Supplementary Data**

|     | REL      | VT        | Concept Name                        | Unit / CODE<br>Value | V7<br>Label |
|-----|----------|-----------|-------------------------------------|----------------------|-------------|
| F-3 | CONTAINS | CONTAINER | (111414, DCM, "Supplementary Data") |                      |             |

|       |          |         |   |  |  |
|-------|----------|---------|---|--|--|
| F-3-1 | CONTAINS | INCLUDE | DTID (4206) Breast Imaging Report Finding Section |  |  |
|-------|----------|---------|---|--|--|

### 9.6.1.5 Breast Imaging Report Finding Section (TID 4206)

**Table 9.6-6  
Breast Imaging Report Finding Section**

|           | REL             | VT        | Concept Name                                  | Unit / CODE Value                | V7 Label                 | Comments          |
|-----------|-----------------|-----------|---|----------------------------------|--------------------------|-------------------|
| F-3-1     | CONTAINS        | CONTAINER | (121070, DCM, "Findings")                     |                                  |                          |                   |
| F-3-1-1   | CONTAINS        | INCLUDE   | DTID (4201) Breast Imaging Procedure Reported |                                  |                          |                   |
| F-3-1-2   | CONTAINS        | CODE      | (121071, DCM, "Finding")                      | (111099, DCM, "Selected region") | Seed Point, Seed Ellipse |                   |
|           |                 |           |   | (M-03000, SRT, "Mass")           | Mass                     |                   |
| F-3-1-2-1 | HAS OBS CONTEXT | TEXT      | (125010, DCM, "Identifier")                   |                                  | '1','2','3',,,           |                   |
| F-3-1-2-2 | HAS PROPERTIES  | INCLUDE   | DTID (4203) Breast Imaging Assessment         |                                  |                          | Ref. Table 9.6-8  |
| F-3-1-2-3 | HAS PROPERTIES  | INCLUDE   | DTID (1400) Linear Measurement                |                                  |                          | Ref. Table 9.6-9  |
| F-3-1-2-4 | HAS PROPERTIES  | INCLUDE   | DTID (1401) Area Measurement                  |                                  |                          | Ref. Table 9.6-10 |
| F-3-1-2-5 | HAS PROPERTIES  | INCLUDE   | DTID (1402) Volume Measurement                |                                  | Vol.                     | Ref. Table 9.6-11 |
| F-3-1-2-6 | HAS PROPERTIES  | NUM       | (99180-1, MDSN, "Angle")                      |                                  | Angle                    |                   |
| F-3-1-2-7 | HAS PROPERTIES  | CODE      | (M-020F9, SRT, "Shape")                       | DCID (6004) Mammography          | Shape                    | Ref. Table 9.6-15 |

|            |                |      |   |  |                        |                   |
|------------|----------------|------|---|--|------------------------|-------------------|
|            |                |      |   | Characteristics of Shape                             |                        |                   |
| F-3-1-2-8  | HAS PROPERTIES | CODE | (111037, DCM, "Margins")                                    | DCID (6006)<br>Mammography Characteristics of Margin | Margin                 | Ref. Table 9.6-16 |
| F-3-1-2-9  | HAS PROPERTIES | CODE | (111009, DCM, "Calcification Type")                         | DCID (6010)<br>Mammography Calcification Types       | Calcifications         | Ref. Table 9.6-17 |
| F-3-1-2-10 | HAS PROPERTIES | CODE | (G-C189, SRT, "Associated Finding")                         | DCID (6056)<br>Associated Findings for Breast        | Special Case           | Ref. Table 9.6-18 |
| F-3-1-2-11 | HAS PROPERTIES | CODE | (111354, DCM, "Orientation")                                | DCID (6152)<br>Orientation                           | Orientation            | Ref. Table 9.6-19 |
| F-3-1-2-12 | HAS PROPERTIES | CODE | (111357, DCM, "Lesion boundary")                            | DCID (6153)<br>Lesion boundary                       | Lesion boundary        | Ref. Table 9.6-20 |
| F-3-1-2-13 | HAS PROPERTIES | CODE | (111360, DCM, "Echo pattern")                               | DCID (6154) Echo pattern                             | Echo pattern           | Ref. Table 9.6-21 |
| F-3-1-2-14 | HAS PROPERTIES | CODE | (111366, DCM, "Posterior acoustic features")                | DCID (6155)<br>Posterior acoustic features           | Posterior Feature      | Ref. Table 9.6-22 |
| F-3-1-2-15 | HAS PROPERTIES | CODE | (111371, DCM, "Identifiable effect on surrounding tissues") | DCID (6015)<br>Single Image Finding from BI-RADS®    | Surrounding Tissue     | Ref. Table 9.6-23 |
| F-3-1-2-16 | HAS PROPERTIES | CODE | (111372, DCM, "Vascularity")                                | DCID (6157)<br>Vascularity                           | Vascularity            | Ref. Table 9.6-24 |
| F-3-1-2-17 | HAS PROPERTIES | CODE | (99600-01, MDSN, "Vascular Abnormalities")                  | DCID (6157)<br>Vascularity                           | Vascular Abnormalities | Ref. Table 9.6-24 |
| F-3-1-2-18 | HAS PROPERTIES | CODE | (99600-02, MDSN, "Elasticity Assessment")                   | DCID (99601)<br>Elasticity Assessment                | Elasticity Assessment  | Ref. Table 9.6-25 |

|         |          |         |  |  |  |                      |
|---------|----------|---------|--|--|--|----------------------|
| F-3-1-3 | CONTAINS | INCLUDE | DTID (SM99610) Breast User<br>Creation Group Section |  |  | Ref. Table<br>9.6-12 |
|---------|----------|---------|--|--|--|----------------------|

### 9.6.1.6 Breast Imaging Procedure Reported (TID 4201)

**Table 9.6-7  
Breast Imaging Procedure Reported**

|           | REL                   | VT   | Concept Name                        | Unit / CODE Value   | V7 Label |
|-----------|-----------------------|------|-------------------------------------|---|----------|
| F-3-1-1   | CONTAINS              | CODE | (121058, DCM, "Procedure reported") | (P5-B8500, SRT,<br>"Ultrasonography of breast")                 |          |
| F-3-1-1-1 | HAS<br>CONCEPT<br>MOD | CODE | (G-C171, SRT, "Laterality")         | (T-04030, SRT, "Left breast")<br>(T-04020, SRT, "Right breast") | Lt, Rt   |

### 9.6.1.7 Breast Imaging Assessment (TID 4203)

**Table 9.6-8  
Breast Imaging Assessment**

|           | REL                | VT   | Concept Name                            | Unit / CODE Value                     | Comments             |
|-----------|--------------------|------|---|---------------------------------------|----------------------|
| F-3-1-2-2 | HAS OBS<br>CONTEXT | CODE | (111005, DCM, "Assessment<br>Category") | DCID (6027) Mammography<br>Assessment | Ref. Table<br>9.6-14 |

### 9.6.1.8 Linear Measurement Template (TID 1400)

**Table 9.6-9  
Linear Measurement**

|           | REL      | VT  | Concept Name                           | Unit / CODE Value                                  | Comments             |
|-----------|----------|-----|--|--|----------------------|
| F-3-1-2-3 | CONTAINS | NUM | Context ID 7470 Linear<br>Measurements | Units = DCID (7460) Units of<br>Linear Measurement | Ref. Table<br>9.6-13 |

### 9.6.1.9 Area Measurement Template (TID 1401)

**Table 9.6-10  
Area Measurement**

|           | REL      | VT  | Concept Name          | Unit / CODE Value                             | Comments |
|-----------|----------|-----|-----------------------|---|----------|
| F-3-1-2-4 | CONTAINS | NUM | (G-A166, SRT, "Area") | Units = DCID (7461) Units of Area Measurement |          |

### 9.6.1.10 Volume Measurement Template (TID 1402)

**Table 9.6-11  
Volume Measurement**

|           | REL      | VT  | Concept Name            | Unit / CODE Value                               | Comments |
|-----------|----------|-----|-------------------------|---|----------|
| F-3-1-2-5 | CONTAINS | NUM | (G-D705, SRT, "Volume") | Units = DCID (7462) Units of Volume Measurement |          |

### 9.6.1.11 Breast User Creation Group Section (TID SM99610)

**Table 9.6-12  
User Creation Group in Breast SR**

|             | REL             | VT        | Concept Name  | Unit / CODE Value     | V7 Label | Comments   |
|-------------|-----------------|-----------|---|-----------------------|----------|--|
| F-3-1-3     | CONTAINS        | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |                       |          | *Creation ID: Randomly generated 7-digit unique ID |
| F-3-1-3-1   | CONTAINS        | NUM       | (99900-Creation ID, MDSN, "User Creation Item Name")  |                       |          |  |
| F-3-1-3-1-1 | HAS CONCEPT MOD | CODE      | (121401, DCM, "Derivation")                           | Common CID-Derivation |          |  |

## 9.6.2 Standard Extended Context Groups in Breast SR

**Table 9.6-13**  
**Context ID 7470 Linear Measurements**

| CSD | CV     | CM       | V7 Label  |
|-----|--------|----------|-----------|
| SRT | G-A22A | Length   | Length, L |
| DCM | 121206 | Distance | Distance  |
| SRT | G-A220 | Width    | Width, W  |
| SRT | G-D785 | Depth    | Depth, D  |
| DCM | 121207 | Height   | Height, H |

**Table 9.6-14**  
**Context ID 6027 Mammography Assessment**

| CSD | CV           | CM   | V7 Label   |
|-----|--------------|--|--|
| BI  | II.AC.a      | 0 - Need additional imaging evaluation                     | Unknown  |
| BI  | II.AC.b.1    | 1-Negative   | Negative   |
| BI  | II.AC.b.2    | 2-Benign Finding   | Benign   |
| BI  | II.AC.b.3    | 3-Probably Benign Finding-short interval follow-up         | Probably benign<br>Possibly Benign               |
| BI  | II.AC.b.4    | 4-Suspicious abnormality, biopsy should be considered      | Possibly Malignant                               |
| BI  | MA.II.A.5.4A | 4A-Low suspicion   | Low suspicion of Malignancy                      |
| BI  | MA.II.A.5.4B | 4B-Intermediate suspicion                                  | Intermediate suspicion of malignancy             |
| BI  | MA.II.A.5.4C | 4C-Moderate suspicion                                      | Moderate concern, but not classic for malignancy |
| BI  | II.AC.b.5    | 5-Highly suggestive of malignancy, take appropriate action | Highly suggestive of malignancy                  |
| BI  | MA.II.A.5.6  | 6-Known biopsy proven malignancy                           | Known biopsy-proven malignancy                   |

**Table 9.6-15**  
**Context ID 6004 Mammography Characteristics of Shape**

| CSD | CV | CM | V7 Label |
|-----|----|----|----------|
|-----|----|----|----------|

|     |         |                    |           |
|-----|---------|--------------------|-----------|
| SRT | M-02100 | Round shape        | Round     |
| SRT | M-02120 | Ovoid shape (Oval) | Oval      |
| SRT | G-A402  | Irregular          | Irregular |

**Table 9.6-16**

**Context ID 6006 Mammography Characteristics of Margin**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>             | <b>V7 Label</b> |
|------------|-----------|-----------------------|-----------------|
| SRT        | F-01741   | Circumscribed lesion  | Circumscribed   |
| SRT        | F-01742   | Microlobulated lesion | Microlobulated  |
| SRT        | F-01744   | Indistinct lesion     | Indistinct      |
| SRT        | F-01745   | Spiculated lesion     | Spiculated      |
| DCM        | 111343    | Angular margins       | Angular         |

**Table 9.6-17**

**Context ID 6010 Mammography Calcification Types**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                        | <b>V7 Label</b>                         |
|------------|-----------|----------------------------------|---|
| SRT        | R-41198   | Unknown                          | Unknown                                 |
| DCM        | 111345    | Macrocalcifications              | Macrocalcifications                     |
| DCM        | 111347    | Calcifications outside of a mass | Micro. Out Of Mass<br>Outside of a mass |
| DCM        | 111346    | Calcifications within a mass     | Macro. In Mass<br>In a mass             |
| MDSN       | 99018-12  | Intraductal                      | Intraductal                             |

**Table 9.6-18**

**Context ID 6056 Associated Findings for Breast**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                        | <b>V7 Label</b>        |
|------------|-----------|----------------------------------|------------------------|
| SRT        | R-41198   | Unknown                          | Unknown                |
| DCM        | 111129    | Clustered microcysts             | Clustered Microcysts   |
| DCM        | 111130    | Complicated cyst                 | Complicated Cysts      |
| SRT        | D7-90382  | Sebaceous cyst of skin of breast | Mass in or on skin     |
| SRT        | M-30400   | Foreign body                     | Foreign Body           |
| SRT        | T-C4351   | Intra-mammary lymph node         | Intra-mam. Lymph Nodes |
| SRT        | T-C4710   | Axillary lymph node              | Axillary Lymph Node    |
| MDSN       | 99018-24  | Foreign Body+Implants            | Foreign Body+Implants  |

|      |          |                        |                        |
|------|----------|------------------------|------------------------|
| MDSN | 99018-21 | Simple Cyst            | Simple Cyst            |
| MDSN | 99018-22 | Post-surg. Fluid Coll. | Post-surg. Fluid Coll. |
| MDSN | 99018-23 | Fat Necrosis           | Fat Necrosis           |

**Table 9.6-19**

**Context ID DCID 6152 Orientation**

| CSD | CV     | CM           | V7 Label     |
|-----|--------|--------------|--------------|
| DCM | 111355 | Parallel     | Parallel     |
| DCM | 111356 | Not parallel | Not parallel |

**Table 9.6-20**

**Context ID 6153 Lesion boundary**

| CSD | CV     | CM               | V7 Label         |
|-----|--------|------------------|------------------|
| DCM | 111358 | Abrupt interface | Abrupt Interface |
| DCM | 111359 | Echogenic halo   | Echogenic Halo   |

**Table 9.6-21**

**Context ID 6154 Echo Pattern**

| CSD  | CV       | CM            | V7 Label   |
|------|----------|---------------|--|
| DCM  | 111361   | Anechoic      | Anechoic   |
| DCM  | 111362   | Hyperechoic   | Hyperechoic                                      |
| DCM  | 111363   | Complex       | Complex Echogenicity<br>Complex Cystic And Solid |
| DCM  | 111364   | Hypoechoic    | Hypoechoic                                       |
| DCM  | 111365   | Isoechoic     | Isoechoic  |
| MDSN | 99018-13 | Heterogeneous | Heterogeneous                                    |

**Table 9.6-22**

**Context ID 6155 Posterior acoustic features**

| CSD | CV     | CM   | V7 Label                                       |
|-----|--------|--|--|
| DCM | 111367 | No posterior acoustic features               | No Posterior Findings<br>No Posterior Features |
| DCM | 111368 | Posterior enhancement                        | Enhancement                                    |
| DCM | 111369 | Posterior shadowing                          | Shadowing                                      |
| DCM | 111370 | Combined posterior enhancement and shadowing | Combined Pattern                               |

**Table 9.6-23**  
**Context ID 6015 Single Image Finding from BI-RADS®**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                          | <b>V7 Label</b>       |
|------------|-----------|------------------------------------|-----------------------|
| SRT        | R-41198   | Unknown                            | Unknown               |
| MDSN       | 99018-11  | Duct Changes                       | Duct Changes          |
| DCM        | 111111    | Cooper's ligament changes          | Cooper Ligament Chg.  |
| SRT        | M-36300   | Edema                              | Edema                 |
| SRT        | F-01795   | Architectural distortion of breast | Architectural Changes |
| SRT        | F-0179A   | Skin thickening of breast          | Skin Thickening       |
| SRT        | F-01799   | Skin retraction of breast          | Skin Retraction       |

**Table 9.6-24**  
**Context ID 6157 Vascularity**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>   | <b>V7 Label</b>      |
|------------|-----------|---|----------------------|
| DCM        | 111373    | Vascularity not present                               | None<br>Absent       |
| DCM        | 111374    | Vascularity not assessed                              | Unknown              |
| DCM        | 111375    | Vascularity present in lesion                         | In Lesion            |
| DCM        | 111376    | Vascularity present immediately adjacent to lesion    | Adjacent To Lesion   |
| DCM        | 111377    | Diffusely increased vascularity in surrounding tissue | Diffuse Vascularity  |
| MDSN       | 99018-14  | Internal Vascularity                                  | Internal Vascularity |
| MDSN       | 99018-15  | Vessels In Rim  | Vessels In Rim       |
| MDSN       | 99018-16  | AVMs  | AVMs                 |
| MDSN       | 99018-17  | Mondor Disease  | Mondor Disease       |

**Table 9.6-25**  
**Context ID 99601 Elasticity Assessment**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>    | <b>V7 Label</b> |
|------------|-----------|--------------|-----------------|
| MDSN       | 99018-18  | Soft         | Soft            |
| MDSN       | 99018-19  | Intermediate | Intermediate    |
| MDSN       | 99018-20  | Hard         | Hard            |

## 9.7 THYROID IMAGING STRUCTURED REPORT TEMPLATE

### 9.7.1 Thyroid Imaging Report (TID SM99700)

**Table 9.7-1**  
**Thyroid Imaging Report**

| No | Rel With Parent       | VT        | Concept Name   | Comments                       | V7 Label |
|----|-----------------------|-----------|--|--------------------------------|----------|
| 1  |                       | CONTAINER | (99700, MDSN, " Soft Tissue Neck and Head Imaging Report") |                                |          |
| 2  | HAS<br>CONCEPT<br>MOD | CODE      | (121049,DCM,"Language of Content Item and Descendants")    | (en,<br>RFC3066,<br>"English") |          |
| 3  | HAS OBS<br>CONTEXT    | INCLUDE   | DTID (1002) Observer Context                               |                                |          |
| 4  | CONTAINS              | INCLUDE   | DTID (SM99701) Thyroid Imaging Finding Section             |                                |          |

#### 9.7.1.1 Observer Context (TID 1002)

**Table 9.7-2**  
**OBSERVER CONTEXT IN THYROID SR**

|     | REL                | VT    | Concept Name                          | Unit / CODE Value   | V7 Label               |
|-----|--------------------|-------|---------------------------------------|---|------------------------|
| G-3 | HAS OBS<br>CONTEXT | CODE  | (121005, DCM, "Observer Type")        | (121006, DCM, "Person")   |                        |
| G-4 | HAS OBS<br>CONTEXT | PNAME | (121008, DCM, "Person Observer Name") |   | Ref.<br>Physician      |
| G-5 | HAS OBS<br>CONTEXT | CODE  | (121024, DCM, "Subject Class")        | (121025 ,DCM,"Patient")   |                        |
| G-6 | HAS OBS<br>CONTEXT | PNAME | (121029,DCM, "Subject Name")          |   | LastName,<br>FirstName |
| G-7 | HAS OBS<br>CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")    |   | BirthDate              |
| G-8 | HAS OBS<br>CONTEXT | CODE  | (121032,DCM, "Subject Sex")           | (M, DCM, "Male")<br>(F, DCM, "Female")<br>(U, DCM, "Unknown sex") | Gender                 |

### 9.7.1.2 Thyroid Imaging Finding Section (TID SM99701)

Table 9.7-3

#### Thyroid Imaging Report Finding Section

|         | REL                   | VT        | Concept Name                           | Unit / CODE Value                             | V7 Label       | Comments  |
|---------|-----------------------|-----------|--|---|----------------|---|
| G-9     | CONTAINS              | CONTAINER | (121070, DCM, "Findings")              |   |                |   |
| G-9-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding Site")          | EV (T-D1600,<br>SRT, "Neck")                  |                |   |
| G-9-2   | HAS<br>CONCEPT<br>MOD | CODE      | (G-0373, SRT, "Image Mode")            | EV (G-03A2, SRT,<br>"2D mode")                |                |   |
| G-9-3   | CONTAINS              | TEXT      | (121106, DCM, "Comment")               |   |                |   |
| G-9-4   | CONTAINS              | CONTAINER | (T-B6000, SRT, "Thyroid")              |   |                |   |
| G-9-4-1 | HAS OBS<br>CONTEXT    | TEXT      | (125010, DCM, "Identifier")            |   | '1','2','3',,, |   |
| G-9-4-2 | HAS<br>PROPERTIES     | INCLUDE   | DTID (1400) Linear<br>Measurement      |   |                | Ref. Table<br>9.6-9   |
| G-9-4-3 | HAS<br>PROPERTIES     | INCLUDE   | DTID (1401) Area<br>Measurement        |   |                | Ref. Table<br>9.6-10  |
| G-9-4-4 | HAS<br>PROPERTIES     | CODE      | (M-020F9, SRT, "Shape")                | DCID (99701)<br>Thyroid Shape                 | Shape          | Ref.<br>Standard<br>Extended<br>Context<br>Groups in<br>Thyroid SR<br>Table 9.7-4 |
| G-9-4-5 | HAS<br>PROPERTIES     | CODE      | (111037, DCM, "Margins")               | DCID (99702)<br>Thyroid Margin                | Margin         | Ref. Table<br>9.7-5   |
| G-9-4-6 | HAS<br>PROPERTIES     | CODE      | (111009, DCM, "Calcification<br>Type") | DCID (99703)<br>Thyroid<br>Calcification Type | Calcifications | Ref. Table<br>9.7-6   |

|          |                   |      |  |  |                        |                      |
|----------|-------------------|------|--|--|------------------------|----------------------|
| G-9-4-7  | HAS<br>PROPERTIES | CODE | (111354, DCM, "Orientation")               | DCID (6152)<br>Orientation   | Orientation            | Ref. Table<br>9.6-19 |
| G-9-4-8  | HAS<br>PROPERTIES | CODE | (99700-05, MDSN, "Taller than<br>Wide")    | DCID (99706)<br>Taller than Wide   | Taller than<br>Wide    | Ref. Table<br>9.7-9  |
| G-9-4-9  | HAS<br>PROPERTIES | CODE | (99700-01, MDSN,<br>"Composition")         | DCID (99704)<br>Thyroid<br>Composition   | Composition            | Ref. Table<br>9.7-7  |
| G-9-4-10 | HAS<br>PROPERTIES | CODE | (99700-06, MDSN, "Nodule<br>Composition")  | DCID (99707)<br>Nodule<br>Composition  | Nodule<br>Composition  | Ref. Table<br>9.7-10 |
| G-9-4-11 | HAS<br>PROPERTIES | CODE | (99700-02, MDSN,<br>"Echogenicity")        | DCID (99705)<br>Thyroid<br>Echogenicity  | Echogenicity           | Ref. Table<br>9.7-8  |
| G-9-4-12 | HAS<br>PROPERTIES | CODE | (99700-03, MDSN,<br>"Spongiform")          | (99019-14, MDSN,<br>"Appearance")<br>(99019-15, MDSN,<br>"Nonappearance")                                      | Spongiform             |                      |
| G-9-4-13 | HAS<br>PROPERTIES | CODE | (99700-04, MDSN, "Central<br>Vascularity") | (99018-19, MDSN,<br>"Unselected")<br>(99019-14, MDSN,<br>"Appearance")<br>(99019-15, MDSN,<br>"Nonappearance") | Central<br>Vascularity |                      |
| G-9-4-14 | HAS<br>PROPERTIES | CODE | (99700-07, MDSN, "Colour<br>flow")         | DCID (99708)<br>Colour flow  | Colour flow            | Ref. Table<br>9.7-11 |
| G-9-4-15 | HAS<br>PROPERTIES | CODE | (110830, DCM, "Elasticity")                | DCID (99601)<br>Elasticity<br>Assessment   | Elasticity             | Ref. Table<br>9.6-25 |
| G-9-4-16 | HAS<br>PROPERTIES | CODE | (99700-08, MDSN, "Cystic<br>Component")    | DCID (99709)<br>Cystic Component   | Cystic<br>Component    | Ref. Table<br>9.7-12 |
| G-9-4-17 | HAS<br>PROPERTIES | CODE | (99700-09, MDSN, "Halo")                   | DCID (99710)<br>Halo   | Halo                   | Ref. Table<br>9.7-13 |
| G-9-4-18 | HAS<br>PROPERTIES | CODE | (99700-10, MDSN, "Extent")                 | DCID (99711)<br>Extent   | Extent                 | Ref. Table<br>9.7-14 |
| G-9-4-19 | HAS<br>PROPERTIES | CODE | (99700-11, MDSN,<br>"Lymphadenopathy")     | DCID (99712)<br>Lymphadenopathy  | Lymphadeno<br>pathy    | Ref. Table<br>9.7-15 |

## 9.7.2 Standard Extended Context Groups in Thyroid SR

**Table 9.7-4**  
**Context ID 99701 Thyroid Shape**

| CSD  | CV       | CM             | V7 Label       |
|------|----------|----------------|----------------|
| SRT  | G-A402   | Irregular      | Irregular      |
| MDSN | 99019-01 | Ovoid to round | Ovoid to round |

**Table 9.7-5**  
**Context ID 99702 Thyroid Margin**

| CSD  | CV       | CM                             | V7 Label                       |
|------|----------|--------------------------------|--------------------------------|
| MDSN | 99019-02 | Well-defined smooth            | Well-defined smooth            |
| MDSN | 99019-03 | Microlobulated/Spiculated      | Microlobulated/Spiculated      |
| MDSN | 99019-04 | Ill-defined                    | Ill-defined                    |
| MDSN | 99019-20 | Well defined                   | Well defined                   |
| MDSN | 99019-21 | Irregular/Lobulated/Spiculated | Irregular/Lobulated/Spiculated |
| MDSN | 99019-22 | Irregular/Lobulated            | Irregular/Lobulated            |
| SRT  | F-01745  | Spiculated lesion              | Spiculated lesion              |

**Table 9.7-6**  
**Context ID 99703 Thyroid Calcification Type**

| CSD  | CV       | CM                  | V7 Label                                   |
|------|----------|---------------------|--|
| MDSN | 99019-05 | Calcifications      | Calcifications                             |
| DCM  | 111345   | Macrocalcifications | Macrocalcifications<br>Macro-calcification |
| MDSN | 99019-07 | Microcalcifications | Microcalcifications<br>Micro-calcification |
| MDSN | 99019-06 | No calcifications   | No calcifications                          |
| MDSN | 99019-23 | Rim/Egg shell       | Rim/Egg shell                              |

**Table 9.7-7**  
**Context ID 99704 Thyroid Composition**

| CSD  | CV       | CM    | V7 Label |
|------|----------|-------|----------|
| MDSN | 99019-19 | Solid | Solid    |

|      |          |                      |                      |
|------|----------|----------------------|----------------------|
| MDSN | 99019-08 | Partially cystic     | Partially cystic     |
| MDSN | 99019-09 | Predominantly Solid  | Predominantly Solid  |
| MDSN | 99019-10 | Predominantly cystic | Predominantly cystic |
| MDSN | 99019-11 | Cystic               | Cystic               |

**Table 9.7-8**  
**Context ID 99705 Thyroid Echogenicity**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>         | <b>V7 Label</b>                           |
|------------|-----------|-------------------|---|
| MDSN       | 99019-12  | Hyper/Isoechoic   | Hyper/Isoechoic                           |
| DCM        | 111364    | Hypoechoic        | Hypoechoic<br>Hypo-echoic                 |
| MDSN       | 99019-13  | Marked Hypoechoic | Marked Hypoechoic<br>Markedly hypo-echoic |
| DCM        | 111362    | Hyperechoic       | Hyper-echoic                              |
| DCM        | 111365    | Isoechoic         | Iso-echoic                                |

**Table 9.7-9**  
**Context ID 99706 Taller than Wide**

| <b>CSD</b> | <b>CV</b> | <b>CM</b> | <b>V7 Label</b> |
|------------|-----------|-----------|-----------------|
| MDSN       | 99019-27  | AP > TR   | AP > TR         |
| MDSN       | 99019-28  | AP <= TR  | AP <= TR        |

**Table 9.7-10**  
**Context ID 99707 Nodule Composition**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>               | <b>V7 Label</b>         |
|------------|-----------|-------------------------|-------------------------|
| MDSN       | 99019-19  | Solid                   | Solid                   |
| MDSN       | 99019-08  | Partially cystic        | Partially cystic        |
| MDSN       | 99019-24  | Mixed solid/cystic      | Mixed solid/cystic      |
| MDSN       | 99019-25  | Micro-cystic/spongiform | Micro-cystic/spongiform |
| MDSN       | 99019-11  | Cystic                  | Cystic                  |

**Table 9.7-11**  
**Context ID 99708 Colour flow**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>  | <b>V7 Label</b> |
|------------|-----------|------------|-----------------|
| MDSN       | 99018-19  | Unselected | Unselected      |

|      |          |            |            |
|------|----------|------------|------------|
| MDSN | 99019-29 | Central    | Central    |
| MDSN | 99019-30 | Peripheral | Peripheral |
| MDSN | 99019-31 | Mixed      | Mixed      |
| MDSN | 99019-32 | None       | None       |

**Table 9.7-12**  
**Context ID 99709 Cystic Component**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                | <b>V7 Label</b>          |
|------------|-----------|--------------------------|--------------------------|
| MDSN       | 99018-19  | Unselected               | Unselected               |
| MDSN       | 99019-33  | Ring down sign - colloid | Ring down sign - colloid |
| MDSN       | 99019-34  | No Ring down sign        | No Ring down sign        |

**Table 9.7-13**  
**Context ID 99710 Halo**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>          | <b>V7 Label</b>    |
|------------|-----------|--------------------|--------------------|
| MDSN       | 99018-19  | Unselected         | Unselected         |
| SRT        | R-4089B   | Absent             | Absent             |
| MDSN       | 99019-35  | Interrupted        | Interrupted        |
| MDSN       | 99019-36  | Regular/continuous | Regular/continuous |

**Table 9.7-14**  
**Context ID 99711 Extent**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                                 | <b>V7 Label</b>                                 |
|------------|-----------|---|---|
| MDSN       | 99018-19  | Unselected                                | Unselected                                      |
| MDSN       | 99019-37  | Retrosternal extension/tracheal deviation | Retrosternal<br>extension/tracheal<br>deviation |
| MDSN       | 99019-32  | None                                      | None  |

**Table 9.7-15**  
**Context ID 99712 Lymphadenopathy**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>            | <b>V7 Label</b>      |
|------------|-----------|----------------------|----------------------|
| MDSN       | 99018-19  | Unselected           | Unselected           |
| MDSN       | 99019-39  | Suspected malignancy | Suspected malignancy |
| MDSN       | 99019-32  | None                 | None                 |

## 9.8 FETAL ECHO STRUCTURED REPORT TEMPLATE

### 9.8.1 FetalEcho Ultrasound Report (TID 5220)

**Table 9.8-1**  
**FetalEcho Ultrasound Report**

| No | Rel With Parent       | VT        | Concept Name   | Comments                       | V7 Label                    |
|----|-----------------------|-----------|--|--------------------------------|-----------------------------|
| 1  |                       | CONTAINER | (125196, DCM, "Fetal Cardiac Ultrasound Report")         |                                |                             |
| 2  | HAS<br>CONCEPT<br>MOD | CODE      | (121049,DCM,"Language of Content Item and Descendants")  | (en,<br>RFC3066,<br>"English") |                             |
| 3  | HAS OBS<br>CONTEXT    | INCLUDE   | DTID (1001) Observer Context                             |                                |                             |
| 4  | CONTAINS              | INCLUDE   | DTID (3602) Cardiovascular Patient Characteristics       |                                |                             |
| 5  | CONTAINS              | INCLUDE   | DTID (5228) Cardiac Ultrasound Fetal Measurement Section |                                |                             |
| 6  | CONTAINS              | INCLUDE   | DTID (SM99010) OB-GYN User Creation Group Section        |                                | Ref.<br>Section<br>9.1.1.18 |

#### 9.8.1.1 Observation Context (TID 1001)

**Table 9.8-2**  
**OBSERVATION CONTEXT IN FETAL ECHO SR**

|     | REL             | VT    | Concept Name                          | Unit / CODE Value       | V7 Label                 |
|-----|-----------------|-------|---------------------------------------|-------------------------|--------------------------|
| H-1 | HAS OBS CONTEXT | CODE  | (121005, DCM, "Observer Type")        | (121006, DCM, "Person") |                          |
| H-2 | HAS OBS CONTEXT | PNAME | (121008, DCM, "Person Observer Name") |                         | Ref. Physician           |
| H-3 | HAS OBS CONTEXT | CODE  | (121024, DCM, "Subject Class")        | (121025 ,DCM,"Patient") |                          |
| H-4 | HAS OBS CONTEXT | PNAME | (121029,DCM, "Subject Name")          |                         | Last Name,<br>First Name |
| H-5 | HAS OBS CONTEXT | DATE  | (121031,DCM, "Subject Birth Date")    |                         | BirthDate                |

### 9.8.1.2 Cardiovascular Patient Characteristics (TID 3602)

**Table 9.8-3  
Patient Characteristics in FetalEcho SR**

|       | REL      | VT        | Concept Name                             | Unit / CODE Value                            | V7 Label    |
|-------|----------|-----------|--|--|-------------|
| H-6   | CONTAINS | CONTAINER | (121118, DCM, "Patient Characteristics") |  |             |
| H-6-1 | CONTAINS | NUM       | (121033, DCM, "Subject Age")             | Units = DCID (7456) Units of Measure for Age | Description |
| H-6-2 | CONTAINS | CODE      | (121032, DCM, "Subject Sex")             | DCID (7455) Sex                              |             |
| H-6-3 | CONTAINS | NUM       | (8302-2, LN, "Patient Height")           | (cm, UCUM, "centimeter")                     | Height      |
| H-6-4 | CONTAINS | NUM       | (29463-7, LN, "Patient Weight")          | (kg, UCUM, "kilograms")                      | Weight      |
| H-6-5 | CONTAINS | NUM       | (8867-4, LN, "Heart Rate")               |  | HR          |

### 9.8.1.3 Cardiac Ultrasound Fetal Measurement Section (TID 5228)

**Table 9.8-4  
Cardiac Ultrasound Fetal Measurement Section**

|       | REL             | VT        | Concept Name  | Unit / CODE Value  | V7 Label       | Comment          |
|-------|-----------------|-----------|---|--|----------------|------------------|
| H-7   | CONTAINS        | CONTAINER | (125016, DCM, "Fetal Measurements")                                     |  |                |                  |
| H-7-1 | HAS OBS CONTEXT | TEXT      | (11951-1, LN, "Fetus ID")   |  |                |                  |
| H-7-2 | CONTAINS        | INCLUDE   | DTID (300) Measurement  | \$Measurement = DCID (12279) Cardiac Ultrasound Fetal General Measurements   |                | Ref. Table 9.8-6 |
| H-7-3 | CONTAINS        | INCLUDE   | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-F6845, SRT, "Ductus arteriosus")<br>\$Measurement = DCID (12218) Echocardiography Congenital | Duct Art       | Ref. Table 9.8-7 |
| H-7-4 | CONTAINS        | INCLUDE   | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-F680F, SRT, "Ductus venosus")<br>\$Measurement = DCID (12218) Echocardiography               | Ductus Venosus |                  |

|        |          |         |   |  |     |                   |
|--------|----------|---------|---|--|-----|-------------------|
|        |          |         |   | Congenital   |     |                   |
| H-7-5  | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-48710, SRT, "Inferior Vena cava")<br>\$Measurement = DCID (12264) Cardiac Ultrasound Venous Return Systemic Measurements | IVC | Ref. Table 9.8-8  |
| H-7-6  | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-32600, SRT, "Left Ventricle")<br>\$Measurement = DCID (12259) Cardiac Ultrasound Ventricles Measurements                 | LV  | Ref. Table 9.8-9  |
| H-7-7  | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-32500, SRT, "Right Ventricle")<br>\$Measurement = DCID (12259) Cardiac Ultrasound Ventricles Measurements                | RV  |                   |
| H-7-8  | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-32300, SRT, "Left Atrium")<br>\$Measurement = DCID (12265) Cardiac Ultrasound Atria and Atrial Septum Measurements       | LA  | Ref. Table 9.8-10 |
| H-7-9  | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-32200, SRT, "Right Atrium")<br>\$Measurement = DCID (12265) Cardiac Ultrasound Atria and Atrial Septum Measurements      | RA  |                   |
| H-7-10 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-35300, SRT, "Mitral Valve")<br>\$Measurement = DCID (12268) Cardiac Ultrasound Atrioventricular Valves Measurements      | MV  | Ref. Table 9.8-11 |

|        |          |         |   |   |            |                   |
|--------|----------|---------|---|---|------------|-------------------|
| H-7-11 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-35100, SRT, "Tricuspid Valve")<br>\$Measurement = DCID (12268 ) Cardiac Ultrasound Atrioventricular Valves Measurements | TV         |                   |
| H-7-12 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-35400, SRT, "Aortic Valve")<br>\$Measurement = DCID (12271) Cardiac Ultrasound Outflow Tracts Measurements              | AV         | Ref. Table 9.8-12 |
| H-7-13 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-35200, SRT, "Pulmonic Valve")<br>\$Measurement = DCID (12271) Cardiac Ultrasound Outflow Tracts Measurements            | PV         |                   |
| H-7-14 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-44000, SRT, "Pulmonary Artery")<br>\$Measurement = DCID (12260) Cardiac Ultrasound Pulmonary Artery                     | PA         | Ref. Table 9.8-13 |
| H-7-15 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (42000, SRT, "Aorta")<br>\$Measurement = DCID (12274) Cardiac Ultrasound Aorta Measurements                                | Aorta      | Ref. Table 9.8-14 |
| H-7-16 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-42100, SRT, "Ascending aorta")<br>\$Measurement = DCID (12274) Cardiac Ultrasound Aorta Measurements                    | Asc. Aorta |                   |
| H-7-18 | CONTAINS | INCLUDE | DTID (5222) Pediatric, Fetal and Congenital Cardiac Ultrasound Section" | \$SectionSubject = (T-D0765, SRT, "Descending Aorta")<br>\$Measurement = DCID   | Dsc. Aorta |                   |

|  |  |  |  |  |  |
|--|--|--|--|--|--|
|  |  |  |  | (12274) Cardiac Ultrasound<br>Aorta Measurements |  |
|--|--|--|--|--|--|

#### 9.8.1.4 Pediatric, Fetal and Congenital Cardiac Ultrasound Section (TID 5222)

**Table 9.8-5  
Pediatric, Fetal and Congenital Cardiac Ultrasound Section**

|       | REL                   | VT        | Concept Name                       | Unit / CODE Value  | V7 Label |
|-------|-----------------------|-----------|------------------------------------|--|----------|
| 1     | CONTAINS              | CONTAINER | (121070, DCM, "Findings")          |  |          |
| 1-1   | HAS<br>CONCEPT<br>MOD | CODE      | (G-C0E3, SRT, "Finding Site")      | \$SectionSubject   |          |
| 1-2   | CONTAINS              | CONTAINER | (125007, DCM, "Measurement Group") |  |          |
| 1-3   | CONTAINS              | INCLUDE   | DTID (300) Measurement             | \$Measurement = \$Measurement<br>\$Derivation = DCID (3627)<br>Measurement Type<br>\$Method = DCID (12227)<br>Echocardiography Measurement<br>Method |          |
| 1-3-1 | HAS ACQ<br>CONTEXT    | CODE      | (G-0373, SRT, "Image Mode")        | (G-03A2, SRT, "2D mode")<br>(G-0394, SRT, "M mode")  |          |
| 1-3-2 | HAS ACQ<br>CONTEXT    | CODE      | (111031, DCM, "Image View")        | (G-A19B, SRT, "Apical two chamber")<br>(G-A19C, SRT, "Apical four chamber")  |          |

#### 9.8.2 Standard Extended Context Groups in FetalEcho SR

**Table 9.8-6  
Context ID 12279 Cardiac Ultrasound Fetal General Measurements**

| CSD  | CV       | CM                                   | V7 Label    |
|------|----------|--------------------------------------|-------------|
| LN   | 11948-7  | Fetal Heart Rate                     | Fetal HR    |
| MDSN | 99116-08 | PR Interval                          | PR Interval |
| MDSN | 99001-06 | Anterior-Posterior Thoracic Diameter | ThD ap      |
| LN   | 11864-6  | Transverse Thoracic Diameter         | ThD trans   |

|      |          |  |              |
|------|----------|--|--------------|
| MDSN | 99001-08 | Anterior-Posterior Cardiac Diameter            | HrtD ap      |
| MDSN | 99001-09 | Transverse Cardiac Diameter                    | HrtD trans   |
| MDSN | 99001-11 | Cardiothoracic Diameter Ratio                  | CTAR(D)      |
| LN   | 33068-8  | Thoracic Area                                  | ThA          |
| LN   | 59075-2  | Cardiac Cross-sectional Area, transverse by US | HrtA         |
| LN   | 59076-0  | Cardiothoracic Area Ratio                      | CTAR(A)      |
| LN   | 11988-3  | Thoracic Circumference                         | ThC          |
| LN   | 59073-7  | Cardiac Circumference, transverse by US        | HrtC         |
| LN   | 59074-5  | Cardiothoracic Circumference Ratio             | CTAR(C)      |
| LN   | 11988-3  | Thoracic Circumference                         | ThC          |
| LN   | 59073-7  | Cardiac Circumference, transverse by US        | HrtC         |
| LN   | 59074-5  | Cardiothoracic Circumference Ratio             | HrtC/ThC     |
| MDSN | 99001-12 | Cardiac Axis                                   | Cardiac Axis |
| LN   | 11726-7  | Peak Systolic Velocity                         | PLI S Vmax   |
| MDSN | 99008-05 | Peak Diastolic Velocity                        | PLI D Vmax   |
| LN   | 11665-7  | Minimum Diastolic Velocity                     | PLI A Vmax   |
| MDSN | 99200-09 | Preload Index                                  | PLI          |

**Table 9.8-7**  
**Context ID 12218 Echocardiography Congenital**

| <b>CSD</b> | <b>CV</b> | <b>CM</b>                            | <b>V7 Label</b> |
|------------|-----------|--------------------------------------|-----------------|
| LN         | 11726-7   | Peak Systolic Velocity               | PS              |
| LN         | 11653-3   | End Diastolic Velocity               | ED              |
| LN         | 11665-7   | Minimum Diastolic Velocity           | MD              |
| LN         | 20247-3   | Peak Gradient                        | PGmax           |
| LN         | 11692-1   | Time averaged peak velocity          | Vmean           |
| LN         | 12144-2   | Systolic to Diastolic Velocity Ratio | S/D             |
| MDSN       | 99200-01  | Diastolic to Systolic Velocity Ratio | D/S             |
| LN         | 12023-8   | Resistivity Index                    | RI              |
| LN         | 12008-9   | Pulsatility Index                    | PI              |
| MDSN       | 99001-13  | Ductus arteriosus diameter           | Duct Art        |
| LN         | 11726-7   | Peak Systolic Velocity               | S               |
| MDSN       | 99008-05  | Peak Diastolic Velocity              | D               |
| LN         | 11665-7   | Minimum Diastolic Velocity           | a               |
| LN         | 11692-1   | Time averaged peak velocity          | TAmx            |

|      |          |                                      |                    |
|------|----------|--------------------------------------|--------------------|
| LN   | 20352-1  | Time Averaged Mean Velocity          | TAm <sub>ean</sub> |
| LN   | 20256-4  | Mean Gradient                        | PG <sub>mean</sub> |
| LN   | 20247-3  | Peak Gradient                        | PG <sub>max</sub>  |
| LN   | 12144-2  | Systolic to Diastolic Velocity Ratio | S/D                |
| MDSN | 99200-01 | Diastolic to Systolic Velocity Ratio | D/S                |
| MDSN | 99200-00 | Systolic to Atrial Velocity Ratio    | S/a                |
| LN   | 12023-8  | Resistivity Index                    | RI                 |
| LN   | 12008-9  | Pulsatility Index                    | PI                 |
| MDSN | 99200-09 | Preload Index                        | PLI                |
| MDSN | 99200-10 | Peak velocity index for veins        | PVIV               |

**Table 9.8-8**

**Context ID 12264 Cardiac Ultrasound Venous Return Systemic Measurements**

| CSD  | CV       | CM                                   | V7 Label           |
|------|----------|--------------------------------------|--------------------|
| LN   | 11726-7  | Peak Systolic Velocity               | S                  |
| MDSN | 99008-05 | Peak Diastolic Velocity              | D                  |
| LN   | 11665-7  | Minimum Diastolic Velocity           | a                  |
| LN   | 11692-1  | Time averaged peak velocity          | V <sub>mean</sub>  |
| LN   | 20256-4  | Mean Gradient                        | PG <sub>mean</sub> |
| LN   | 20247-3  | Peak Gradient                        | PG <sub>max</sub>  |
| LN   | 12144-2  | Systolic to Diastolic Velocity Ratio | S/D                |
| MDSN | 99200-01 | Diastolic to Systolic Velocity Ratio | D/S                |
| MDSN | 99200-00 | Systolic to Atrial Velocity Ratio    | S/a                |
| LN   | 12023-8  | Resistivity Index                    | RI                 |
| LN   | 12008-9  | Pulsatility Index                    | PI                 |
| MDSN | 99200-09 | Preload Index                        | PLI                |
| MDSN | 99200-10 | Peak velocity index for veins        | PVIV               |
| LN   | 18006-7  | Inferior Vena Cava Diameter          | IVC                |

**Table 9.8-9**

**Context ID 12259 Cardiac Ultrasound Ventricles Measurements**

| CSD | CV      | CM  | Modifier                   | V7 Label             |
|-----|---------|---|----------------------------|----------------------|
| LN  | 29436-3 | Left Ventricle Internal End Diastolic Dimension | ImgMode=SRT\G-03A2\2D mode | LVLd A2C             |
| SRT | G-0375  | Left Ventricular Diastolic Area                 | ImgView=SRT\G-A19B\Apical  | LVA <sub>d</sub> A2C |

|      |          |   |                             |           |
|------|----------|---|-----------------------------|-----------|
| LN   | 18026-5  | Left Ventricular End Diastolic Volume           | two chamber                 | LVEDV A2C |
| LN   | 29438-9  | Left Ventricle Internal Systolic Dimension      |                             | LVLs A2C  |
| SRT  | G-0374   | Left Ventricular Systolic Area                  |                             | LVAAs A2C |
| LN   | 18148-7  | Left Ventricular End Systolic Volume            |                             | LVESV A2C |
| SRT  | F-32120  | Stroke Volume                                   |                             | LV SV A2C |
| SRT  | F-32100  | Cardiac Output                                  |                             | LV CO A2C |
| LN   | 18043-0  | Left Ventricular Ejection Fraction              |                             | LV EF A2C |
| LN   | 29436-3  | Left Ventricle Internal End Diastolic Dimension |                             | LVLd A4C  |
| SRT  | G-0375   | Left Ventricular Diastolic Area                 |                             | LVA d A4C |
| LN   | 18026-5  | Left Ventricular End Diastolic Volume           | ImgMode=SRT\G-03A2\2D       | LVEDV A4C |
| LN   | 29438-9  | Left Ventricle Internal Systolic Dimension      | mode                        | LVLs A4C  |
| SRT  | G-0374   | Left Ventricular Systolic Area                  | ImgView=SRT\G-A19C\Apical   | LVAAs A4C |
| LN   | 18148-7  | Left Ventricular End Systolic Volume            | four chamber                | LVESV A4C |
| SRT  | F-32120  | Stroke Volume                                   |                             | LV SV A4C |
| SRT  | F-32100  | Cardiac Output                                  |                             | LV CO A4C |
| LN   | 18043-0  | Left Ventricular Ejection Fraction              |                             | LV EF A4C |
| LN   | 18026-5  | Left Ventricular End Diastolic Volume           | Measurement                 | LVEDV BP  |
| LN   | 18148-7  | Left Ventricular End Systolic Volume            | Method=DCM\125207\Method of | LVESV BP  |
| SRT  | F-32120  | Stroke Volume                                   | Disks, Biplane              | LV SV BP  |
| SRT  | F-32100  | Cardiac Output                                  | ImgMode=SRT\G-03A2\2D       | LV CO BP  |
| LN   | 18043-0  | Left Ventricular Ejection Fraction              | mode                        | LV EF BP  |
| MDSN | 99104-13 | Left Ventricular Diameter                       |                             | LV Diam   |
| MDSN | 99104-14 | Left Ventricular Length                         |                             | LV Length |
| LN   | 18054-7  | Interventricular Septum % Thickening            |                             | IVS       |
| LN   | 29436-3  | Left Ventricle Internal End Diastolic Dimension |                             | LVIDd     |
| LN   | 29438-9  | Left Ventricle Internal Systolic Dimension      |                             | LVIDs     |
| LN   | 18053-9  | Left Ventricle Posterior Wall % Thickening      | ImgMode=SRT\G-03A2\2D       | LVPW      |
| MDSN | 99104-11 | Left Ventricular Inlet                          | mode                        | LV Inlet  |
| MDSN | 99104-12 | Left Ventricular Area                           |                             | LV Area   |
| MDSN | 99104-09 | Relative Wall Thickness                         |                             | RWT       |
| LN   | 18026-5  | Left Ventricular End Diastolic Volume           |                             | LV Vol. d |
| LN   | 18148-7  | Left Ventricular End Systolic Volume            |                             | LV Vol. s |
| SRT  | F-32120  | Stroke Volume                                   |                             | SV        |

|      |          |   |  |             |
|------|----------|---|--|-------------|
| SRT  | F-32100  | Cardiac Output                                    |  | CO          |
| LN   | 18043-0  | Left Ventricular Ejection Fraction                |  | EF          |
| LN   | 18154-5  | Interventricular Septum Diastolic Thickness       | ImgMode=SRT\G-0394\M mode                      | IVSd        |
| LN   | 29436-3  | Left Ventricle Internal End Diastolic Dimension   |  | LVIDd       |
| LN   | 18152-9  | Left Ventricle Posterior Wall Diastolic Thickness |  | LVPWd       |
| MDSN | 99104-09 | Relative Wall Thickness                           |  | RWT         |
| LN   | 18158-6  | Interventricular Septum Systolic Thickness        |  | IVSs        |
| LN   | 29438-9  | Left Ventricle Internal Systolic Dimension        |  | LVIDs       |
| LN   | 18156-0  | Left Ventricle Posterior Wall Systolic Thickness  |  | LVPWs       |
| LN   | 18026-5  | Left Ventricular End Diastolic Volume             |  | LV EDV      |
| LN   | 18148-7  | Left Ventricular End Systolic Volume              |  | LV ESV      |
| LN   | 18051-3  | Left Ventricular Fractional Shortening            |  | FS          |
| SRT  | F-32120  | Stroke Volume                                     |  | SV          |
| SRT  | F-32100  | Cardiac Output                                    |  | CO          |
| LN   | 18043-0  | Left Ventricular Ejection Fraction                |  | EF          |
| LN   | 18087-7  | Left Ventricle Mass                               |  | LV Mass     |
| SRT  | G-037E   | Left Ventricular Isovolumic Contraction Time      |  | IVCT        |
| DCM  | 122211   | Left Ventricular ejection time                    |  | ET          |
| LN   | 18071-1  | Left Ventricular Isovolumic Relaxation Time       |  | IVRT        |
| MDSN | 99104-04 | Left Ventricle Total Systolic Time                |  | MCO         |
| LN   | 59099-2  | Myocardial Performance Index (Tei)                |  | Tei Index   |
| SRT  | G-038F   | Cardiovascular Orifice Diameter                   | FSite=SRT\T-32650\Left Ventricle Outflow Tract | LVOT Diam   |
| SRT  | G-038E   | Cardiovascular Orifice Area                       |  | LVOT Area   |
| LN   | 11726-7  | Peak Velocity                                     |  | LVOT Vmax   |
| LN   | 20247-3  | Peak Gradient                                     |  | LVOT PGmax  |
| LN   | 20352-1  | Mean Velocity                                     |  | LVOT Vmean  |
| LN   | 20256-4  | Mean Gradient                                     |  | LVOT PGmean |

|      |          |  |  |                               |
|------|----------|--|--|-------------------------------|
| LN   | 20354-7  | Velocity Time Integral                         |  | LVOT VTI                      |
| SRT  | F-32120  | Stroke Volume                                  |  | LVOT SV                       |
| SRT  | F-32100  | Cardiac Output                                 |  | LVOT CO                       |
| SRT  | G-038F   | Cardiovascular Orifice Diameter                | FSite=SRT\T-32550\Right<br>Ventricle Outflow Tract | RVOT Diam                     |
| SRT  | G-038E   | Cardiovascular Orifice Area                    |  | RVOT Area                     |
| LN   | 11726-7  | Peak Velocity                                  |  | RVOT Vmax                     |
| LN   | 20247-3  | Peak Gradient                                  |  | RVOT PGmax                    |
| LN   | 20352-1  | Mean Velocity                                  |  | RVOT Vmean                    |
| LN   | 20256-4  | Mean Gradient                                  |  | RVOT<br>PGmean                |
| LN   | 20354-7  | Velocity Time Integral                         |  | RVOT VTI                      |
| SRT  | F-32120  | Stroke Volume                                  |  | RVOT SV                       |
| SRT  | F-32100  | Cardiac Output                                 |  | RVOT CO                       |
| MDSN | 99105-17 | Right Ventricular Diameter                     |  | ImgMode=SRT\G-03A2\2D<br>mode |
| MDSN | 99105-18 | Right Ventricular Length                       |  | RV Length                     |
| MDSN | 99105-19 | Right Ventricular Inlet                        |  | RV Inlet                      |
| MDSN | 99105-20 | Right Ventricular Area                         |  | RV Area                       |
| MDSN | 99105-21 | Right Ventricular Wall                         |  | RV Wall                       |
| LN   | 20304-2  | Right Ventricular Internal Diastolic Dimension | ImgMode=SRT\G-0394\M mode                          | RVDd                          |

**Table 9.8-10**

**Context ID 12265 Cardiac Ultrasound Atria and Atrial Septum Measurements**

| CSD  | CV       | CM                    | V7 Label |
|------|----------|-----------------------|----------|
| MDSN | 99106-06 | Left Atrium Diameter  | LA Diam  |
| MDSN | 99107-06 | Right Atrium Diameter | RA Diam  |

**Table 9.8-11**

**Context ID 12268 Cardiac Ultrasound Atrioventricular Valves Measurements**

| CSD | CV      | CM                                | V7 Label |
|-----|---------|-----------------------------------|----------|
| LN  | 18037-2 | Mitral Valve E-Wave Peak Velocity | E        |
| LN  | 17978-8 | Mitral Valve A-Wave Peak Velocity | A        |
| LN  | 18038-0 | Mitral Valve E to A Ratio         | E/A      |

|      |          |                                      |                 |
|------|----------|--------------------------------------|-----------------|
| LN   | 80062-3  | Mitral regurgitation Vmax            | MR Vmax         |
| MDSN | 99114-14 | Mitral Valve Diameter                | Mitral Valve    |
| LN   | 18031-5  | Tricuspid Valve E Wave Peak Velocity | E               |
| LN   | 18030-7  | Tricuspid Valve A Wave Peak Velocity | A               |
| LN   | 18039-8  | Tricuspid Valve E to A Ratio         | E/A             |
| LN   | 79921-3  | Tricuspid regurgitation Vmax         | TR Vmax         |
| MDSN | 99115-14 | Tricuspid Valve Diameter             | Tricuspid Valve |

**Table 9.8-12**

**Context ID 12271 Cardiac Ultrasound Outflow Tracts Measurements**

| CSD  | CV       | CM                       | V7 Label        |
|------|----------|--------------------------|-----------------|
| MDSN | 99108-11 | Aortic Valve Diameter    | Aortic Valve    |
| MDSN | 99116-07 | Pulmonary Valve Diameter | Pulmonary Valve |

**Table 9.8-13**

**Context ID 12260 Cardiac Ultrasound Pulmonary Artery**

| CSD  | CV       | CM                                   | V7 Label |
|------|----------|--------------------------------------|----------|
| LN   | 11726-7  | Peak Systolic Velocity               | PS       |
| LN   | 11653-3  | End Diastolic Velocity               | ED       |
| LN   | 11665-7  | Minimum Diastolic Velocity           | MD       |
| LN   | 20247-3  | Peak Gradient                        | PGmax    |
| LN   | 11692-1  | Time averaged peak velocity          | TAMean   |
| LN   | 12144-2  | Systolic to Diastolic Velocity Ratio | S/D      |
| MDSN | 99200-01 | Diastolic to Systolic Velocity Ratio | D/S      |
| LN   | 12023-8  | Resistivity Index                    | RI       |
| LN   | 12008-9  | Pulsatility Index                    | PI       |
| LN   | 18020-8  | Main Pulmonary Artery Diameter       | MPA Diam |
| MDSN | 99117-01 | Pulmonary Artery Diameter            | PA Diam  |
| LN   | 18021-6  | Right Pulmonary Artery Diameter      | RPA      |
| LN   | 18019-0  | Left Pulmonary Artery Diameter       | LPA      |

**Table 9.8-14**

**Context ID 12274 Cardiac Ultrasound Aorta Measurements**

| CSD | CV      | CM                     | Finding Site | V7 Label |
|-----|---------|------------------------|--------------|----------|
| LN  | 11726-7 | Peak Systolic Velocity |              | PS       |

|      |          |                                      |                                       |                                    |
|------|----------|--------------------------------------|---------------------------------------|------------------------------------|
| LN   | 11653-3  | End Diastolic Velocity               |                                       | ED                                 |
| LN   | 11665-7  | Minimum Diastolic Velocity           |                                       | MD                                 |
| LN   | 20247-3  | Peak Gradient                        |                                       | PGmax                              |
| LN   | 11692-1  | Time averaged peak velocity          |                                       | TAm <sub>ean</sub>                 |
| LN   | 12144-2  | Systolic to Diastolic Velocity Ratio |                                       | S/D                                |
| MDSN | 99200-01 | Diastolic to Systolic Velocity Ratio |                                       | D/S                                |
| LN   | 12023-8  | Resistivity Index                    |                                       | RI                                 |
| LN   | 12008-9  | Pulsatility Index                    |                                       | PI                                 |
| LN   | 18012-5  | Ascending Aortic Diameter            |                                       | Asc. Aorta                         |
| LN   | 18013-3  | Descending Aortic Diameter           |                                       | Desc. Aorta                        |
| SRT  | T-42310  | Aortic Isthmus                       |                                       | Ao Isthmus                         |
| LN   | 11726-7  | Peak Systolic Velocity               | SRT\T-42304\Transverse<br>Aortic Arch | Transverse<br>Aortic Arch -<br>PSV |

## 9.9 PEDIATRIC STRUCTURED REPORT TEMPLATE

### 9.9.1 Pediatric Ultrasound Report (TID SM99800)

**Table 9.9-1  
Pediatric Ultrasound Report**

| No    | Rel With Parent | VT        | Concept Name  | Unit / CODE Value  | V7 Label                                       |
|-------|-----------------|-----------|---|--|--|
| 1     |                 | CONTAINER | (99800, MDSN, "Pediatric Ultrasound Report")  |  |  |
| 2     | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observer Context  |  |  |
| 3     | CONTAINS        | CONTAINER | (121118, DCM, "Patient Characteristics")  |  |  |
| 3-1   | CONTAINS        | NUM       | (8302-2, LN, "Patient Height")  |  |  |
| 3-2   | CONTAINS        | NUM       | (29463-7, LN, "Patient Weight")   |  |  |
| 3-3   | CONTAINS        | NUM       | (8867-4, LN, "Heart Rate")  |  | HR   |
| 3-4   | CONTAINS        | NUM       | (8277-6, LN, "Body Surface Area")   |  | BSA  |
| 3-4-1 | INFERRED FROM   | CODE      | (8278-4, LN, "Body Surface Area Formula")   | (122241, DCM, "BSA = $0.007184 * WT^{0.425} * HT^{0.725}$ ") |  |
| 4     | CONTAINS        | CONTAINER | (99800-01, MDSN, "Right Hip")<br>(99800-02, MDSN, "Left Hip")<br>(99800-03, MDSN, "Coverage of femoral head") |  | Rt. Hip<br>Lt. Hip<br>Coverage of femoral head |
| 4-1   | CONTAINS        | NUM       | DCID (99800) Pediatric Hip Ultrasound Measurements  | Ref.Table 9.9-3  |  |
| 5     | CONTAINS        | INCLUDE   | DTID (SM99810) Pediatric User Creation Group Section  | Ref.Table 9.9-3  |  |

#### 9.9.1.1 Pediatric User Creation Group Section (TID SM99810)

**Table 9.9-2  
User Creation Group in PEDIATRIC SR**

|  | REL | VT | Concept Name | Unit / CODE | V7 | Comments |
|--|-----|----|--------------|-------------|----|----------|
|--|-----|----|--------------|-------------|----|----------|

|       |                    |           |   | Value                     | Label |   |
|-------|--------------------|-----------|---|---------------------------|-------|---|
| 5     | CONTAINS           | CONTAINER | (99900-Creation ID, MDSN, "User Creation Group Name") |                           |       | *Creation ID:<br>Randomly<br>generated 7-digit<br>unique ID |
| 5-1   | CONTAINS           | NUM       | (99900-Creation ID, MDSN, "User Creation Item Name")  |                           |       |   |
| 5-1-1 | HAS CONCEPT<br>MOD | CODE      | (121401, DCM, "Derivation")                           | Common CID-<br>Derivation |       |   |

## 9.9.2 Standard Extended Context Groups in Pediatric SR

**Table 9.9-3**  
**Context ID 99800 Pediatric Hip Ultrasound Measurements**

| CSD  | CV       | CM          | V7 Label  |
|------|----------|-------------|-----------|
| MDSN | 99800-04 | Alpha Angle | Alpha     |
| MDSN | 99800-05 | Beta Angle  | Beta      |
| MDSN | 99800-13 | Hip Type    | Hip Type  |
| MDSN | 99800-06 | Hip Angle   | Hip Angle |

## 9.10 MUSCULOSKELETAL STRUCTURED REPORT TEMPLATE

### 9.10.1 Musculoskeletal Ultrasound Report (TID SM99900)

**Table 9.10-1**  
**Musculoskeletal Ultrasound Report**

| No  | Rel With Parent    | VT        | Concept Name                                       | Comments | V7 Label |
|-----|--------------------|-----------|--|----------|----------|
| 1   |                    | CONTAINER | (99900, MDSN, "Musculoskeletal Ultrasound Report") |          |          |
| 2   | HAS OBS<br>CONTEXT | INCLUDE   | DTID (1001) Observer Context                       |          |          |
| 3   | CONTAINS           | CONTAINER | (121118, DCM, "Patient Characteristics")           |          |          |
| 3-1 | CONTAINS           | NUM       | (8302-2, LN, "Patient Height")                     |          |          |

|         |                 |           |   |  |     |
|---------|-----------------|-----------|---|--|-----|
| 3-2     | CONTAINS        | NUM       | (29463-7, LN, "Patient Weight")           |  |     |
| 3-3     | CONTAINS        | NUM       | (8867-4, LN, "Heart Rate")                |  | HR  |
| 3-4     | CONTAINS        | NUM       | (8277-6, LN, "Body Surface Area")         |  | BSA |
| 3-4-1   | INFERRED FROM   | CODE      | (8278-4, LN, "Body Surface Area Formula") | (122241, DCM, "BSA = 0.007184*WT^0.4 25*HT^0.725") |     |
| 4       | CONTAINS        | CONTAINER | (121070, DCM, "Findings")                 |  |     |
| 4-1     | HAS CONCEPT MOD | CODE      | (G-C0E3, SRT, "Finding Site")             | (99900-01, MDSN, "Shoulder")                       |     |
|         |                 |           |   | (99900-02, MDSN, "Wrist")                          |     |
|         |                 |           |   | (99900-03, MDSN, "Knee")                           |     |
|         |                 |           |   | (99900-04, MDSN, "Ankle")                          |     |
| 4-2     | HAS CONCEPT MOD | CODE      | (G-C171, SRT, "Laterality")               | (G-A100, SRT, "Right")                             |     |
|         |                 |           |   | (G-A101, SRT, "Left")                              |     |
|         |                 |           |   | (G-A102, SRT, "Unilateral")                        |     |
| 4-3     | CONTAINS        | CONTAINER | (125007, DCM, "Measurement Group")        |  |     |
| 4-3-1   | HAS OBS CONTEXT | TEXT      | (125010, DCM, "Identifier")               | 1~10   |     |
| 4-3-2   | CONTAINS        | NUM       | (121206, DCM, "Distance")                 |  |     |
| 4-3-2-1 | HAS CONCEPT MOD | CODE      | (121401, DCM, "Derivation")               | Common CID-Derivation                              |     |

END OF DOCUMENT

